


OBJECT FOLDING BOX GINEPINK RELIEF VALUEMED PHARMA		ValueMED PHARMA	color 1	C
ARTWORK 6031			color 2	M
VERSION N. 0		DIE CUT	color 3	Y
date 30/05/2019		DD0030	color 4	K
DESIGNER		TECHNICAL CHECK	p. finito 10011	
		SIGNED - OTHER	SIGNED - CUSTOMER	



6031  
10011

GYNEPINK<sup>®</sup> relief  
CE Medical Device

GYNEPINK<sup>®</sup> relief  
CE Medical Device

GYNEPINK<sup>®</sup> relief

CE Medical Device

5 vaginal douches | 100 ml



CE Medical Device

Use by: see the bottom of the box

LOT Lot: see the bottom of the box

Net Content: 500 ml e

ValueMED Pharma Srl  
Via 28 Luglio, 211 - 47893 Borgo Maggiore  
Rep. San Marino - www.valuedmedpharma.com



8 054956 421015

LOT





OBJECT <b>TASSELLO</b> <b>GINEPINK RELIEF VALUEMED PHARMA</b>				
ARTWORK <b>6033</b> date 29/05/2019			VERSION N. 0	DIE CUT <b>56 x 35 mm</b>
DESIGNER	TECHNICAL CHECK	SIGNED - OTHER	SIGNED - CUSTOMER	p. finito    10011

# GYNEPINK<sup>®</sup> relief

CE Medical Device

## CONTENT

5 vaginal douches 100 ml each.

## COMPOSITION

Water, polysorbate 20, Aloe Vera gel 10x, Boric acid, Lactic acid, sodium benzoate, potassium sorbate, Dipotassium glycyrrhizate, bisodic EDTA, Lavender (*Lavandula officinalis*) flowers essential oil, Clove (*Syzygium aromaticum*) flower buds essential oil, alpha-Bisabolol.

## INDICATIONS

GYNEPINK<sup>®</sup> relief is based on Boric acid, Aloe Vera gel, essential oils, Lactic acid and alpha-Bisabolol. GYNEPINK<sup>®</sup> relief effectively cleanse vaginal mucosa and decreases the adhesion properties of undesired microorganisms. Boric acid is both antifungal and antibacterial, while Aloe Vera favors the healing process of intra-vaginal wounds. Clove oil gives relief from itching symptoms while Lactic acid rebalances pH to its physiological acidity.

## WARNINGS

Do not use in case of evident allergy and/or hypersensitivity to any of the ingredients. Keep out from the reach of children. Prolonged use of topical/local skin products can lead to sensitization processes: if irritation occurs, it is recommended to discontinue the treatment and consult a physician. Do not use after the expiry date. The expiry date is referred to product in its intact correctly stored package. Do not use the product if packaging is damaged. No interactions known at the given instructions of use. Do not use during pregnancy or lactation.

## INSTRUCTIONS OF USE

1. before using the product, it is recommended to carefully clean hands
2. shake well the bottle before use
3. keep the bottle from its neck and rotate the cap until the seal is broken
4. apply the cannula onto the bottle neck using a firm pressure to lock it

5. find a comfortable position, anywhere you like (sitting on the toilet, in the bath tub or under the shower: it's important that position ensures the vaginal douche to be flushed away
6. squeeze the bottle; repeat if necessary to empty the bottle
7. dry yourself externally without rinse
8. the product is for single-use only. Discard the bottle after use

Use the vaginal douche preferably before bedtime.

## SHELF-LIFE

3 years. Do not use the product after the Expiry date.

## HANDLING AND STORAGE CONDITIONS

Store in cool dry place away from direct light and heat sources. Place out of sight and reach of children.

## MANUFACTURER

Company name: ValueMed Pharma Srl  
Address: Via 28 Luglio, 211 - 47893 Borgo Maggiore  
Country: Republic of San Marino

Last review  
June 2019

6032

Fronte

Retro

OBJECT LEAFLET GINEPINK RELIEF VALUEMED PHARMA					
ARTWORK 6032 date 30/05/2019			VERSION N. 0	DIE CUT 100x37,5 mm FRONT	
DESIGNER	TECHNICAL CHECK	SIGNED - OTHER	SIGNED - CUSTOMER	color 1	K
				p. finito	10011