



RIQASnet - Login

RIQASnet QA programme is designed to monitor the performance of 17 routine coagulation tests and provides factor specific. Two flexible yet cost effective programme options are available.

Accredited to ISO/IEC 17043:2010

- Lyophilised for enhanced stability
- 100% human plasma
- Monthly reporting
- Choose between 5 routine parameters or all 17 parameters
- Submit results and view reports online via RIQASNet
- Register up to five instruments at no extra cost

Kit No	Kit Size	Frequency	Cycle Start	Parameters
RIQ91352a	12 x 1ml	Monthly (1 x 12 month cycle)	January	5 selected parameters (aPTT, PT, TT, Fibrinogen and A1III)
RIQ91352b	12 x 1ml	Monthly (1 x 12 month cycle)	January	All 17 parameters

Routine Parameters

- Activated Partial Thromboplastin Time (aPTT)
- Antithrombin III (ATIII)
- Fibrinogen
- Prothrombin Time (PT)
- Thrombin Time (TT)

Parameters

- Activated Partial Thromboplastin Time (aPTT)
- Antithrombin III (ATIII)
- D-Dimer (Pilot)
- Factor II
- Factor IX
- Factor V
- Factor VII
- Factor VIII
- Factor X
- Factor XI
- Factor XII
- Fibrinogen
- Plasminogen
- Protein C
- Protein S
- Prothrombin Time (PT)
- Thrombin Time (TT)



RQ9135
 COAGULATION PROGRAMME
 PROGRAMME COAGULATION
 PROGRAMMA COAGULAZIONE
 PROGRAMA DE COAGULACION
 PROGRAMA DE COAGULAÇÃO
 KREŠEJIMO PROGRAMA
 โปรแกรม COAGULATION
 PROGRAM KOAGULOLOGICZNY:
 血凝程序
 CHƯƠNG TRÌNH NGOẠI KIỂM ĐỒNG MÁU
 KOAGULASYON PROGRAMI



ENGLISH
 COAGULATION PROGRAMME: RQ9135

CONFIRMATION OF KIT CHARACTERISTICS AND RECEIPT DATE

Please confirm that the correct number of samples are present and that your samples have the appearance as indicated in the CHARACTERISTICS section below. Please confirm that none of the vials are broken and notify your local Rando representative immediately if there are any discrepancies. Finally, please log on to www.rando.com to confirm the exact date on which you received this kit.

CHARACTERISTICS

The pack contains 6 vials of lyophilised material (6 x 1 ml). The vials are labelled with the sample number.

PREPARATION/STORAGE/STABILITY OF SAMPLES

The samples are sealed under vacuum. Open the vial very carefully, avoiding any loss of material and using a calibrated pipette reconstitute with an accurately measured 1 ml volume of freshly double distilled water at +20°C to +25°C. Replace the rubber stopper and ensure that samples are dissolved completely by swirling gently (ideally place on a roller for half an hour prior to analysis). Do not shake the vials. Analyse samples within 6 hours of reconstitution. The samples should be treated in the same way as patient samples. The samples should be stored at 2 - 8°C when not in use.

SAFETY

Warning: Potentially Biohazardous Material

Human source material from which this product has been derived has been tested at donor level for the Human Immunodeficiency Virus (HIV 1, HIV 2) antibody, Hepatitis B surface Antigen (HbsAg), and Hepatitis C Virus (HCV) antibody and found to be NON-REACTIVE. FDA approved methods have been used to conduct these tests. However, since no method can offer complete assurance as to the absence of infectious agents, this material and all patient samples should be handled as though capable of transmitting infectious diseases and disposed of accordingly.

PLEASE NOTE Some users of ACL TOP instruments when used in combination with Hemosil Synthasil reagents may be unable to attain a clot with extended read time samples.

For IN VITRO use only. Do not pipette by mouth. Exercise the normal precautions required for handling laboratory reagents.

*** IMPORTANT NOTE:** Results must arrive at RIQAS by 17:00 hrs GMT on the FINAL DATE. If the RECOMMENDED ANALYSIS DATE gives insufficient time, we suggest that the sample is analysed earlier to ensure you meet the deadline. If you are faxing results, please transmit 3 working days before the FINAL DATE. Late results will not be accepted after the final date for the next sample.

COAGULATION PROGRAMME / PROGRAMME COAGULATION / PROGRAMMA COAGULAZIONE /
 PROGRAMA DE COAGULACION / PROGRAMA DE COAGULAÇÃO / KREŠEJIMO PROGRAMA / โปรแกรม Coagulation
 PROGRAM KOAGULOLOGICZNY / 血凝程序 / CHƯƠNG TRÌNH NGOẠI KIỂM ĐỒNG MÁU / KOAGULASYON PROGRAMI / PROGRAM KOAGULÁCIE

RETURN OF RESULTS / RETOUR DES RESULTATS / INVIO DEI RISULTATI / ENVIO DE RESULTADOS / REZULTATŲ GRAŽINIMAS /
 การส่งผลการ RIQAS กลับ / PRZESYŁANIE WYNIKÓW / 回复结果 / GỬI TRẢ KẾT QUẢ / SONUÇ GÖNDERİM TARİHLERİ / ODOSLANIE VÝSLEDKOV

CYCLE 12B / CICLO 12B / CIKLAS 12B / CYCLE 12B / CYKL 12B / 循环 12B / CHU KỶ 12B / DÖNEM 12B / CYKLUS 12B

SAMPLE NO/ N°ECH CAMPIONE N° / MUESTRA NO / AMOSTRA N° / MÉGINIO NR. / ตัวอย่างตรวจทั้ง / NUMER PRÓBK I / 样品号 MẪU SỐ/ ÖRNEK NO / Č.VZORKY	RECOMMENDED ANALYSIS DATE/ DATE RECOMMANDEE POUR L'ANALYSE / DATA DI ANALISI RACCOMANDATA/ FECHA RECOMENDADA DE ANALISIS/ DATA RECOMENDADA PARA ANÁLISE/ ZALECANA DATA OZNACZENIA / / REKOMENDUOJAMA TYRIMO ATLIKIMO DATA /วันที่แนะนำให้ทำการวิเคราะห์ / 推荐的分析日期 / NGÀY KHUYẾN CÁO PHẦN TÍCH/ ONERILEN ANALİZ TARİHİ / ODPORÁČANÝ DÁTUM TESTOVANIA	*FINAL DATE / *DATE FINALE/ * ULTIMA DATA / *FECHA FINAL / *DATA FINAL / * GALUTINÉ DATA *วันสุดท้ายของการส่งผลกลับ / DATA FINALNA / * *最终日期 / HẠN CUỐI GỬI TRẢ KẾT QUẢ/ *SON SONUÇ GÖNDERİM TARİHİ / * FINÁLNY DÁTUM
--	--	---

7	13.07.20	20.07.20
8	10.08.20	17.08.20
9	14.09.20	21.09.20
10	12.10.20	19.10.20
11	09.11.20	16.11.20
12	14.12.20	21.12.20

