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I am Paul Cornelis Verhulst, holder of Passport no. BL6439JJ2 issued on 5 October 2016 in the Netherlands. I am the General Director cum Legal Representative of Medtronic Vietnam Company Limited, a company duly incorporated under the Enterprise Registration Certificate No. 0313651200 issued for the first time by Department of Planning and Investment of Ho Chi Minh City on 4 February 2016, as amended from time to time (the "**Company**"). I hereby appoint the following person:

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Unless otherwise instructed by me, the Authorized Representative is authorized to, on my behalf and for the benefit of the Company, conduct any and all necessary actions and handle (including sign) all the documents described in Appendix A attached to this Power of Attorney ("**Authorized Documents**").

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I hereby undertake to confirm and ratify each and everything that the Authorized Representative shall lawfully do or purport to do by virtue of this Power of Attorney, including the confirmation and ratification of anything done honestly and bona fide on my behalf, to perform and complete the aforementioned tasks.

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The Authorized Representative must comply with all relevant policies, processes and procedures of the Company and applicable laws when exercising her authority granted and duties appointed under this Power of Attorney.

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The Authorized Representative is not allowed to sub-authorize any scope of tasks authorized in this Power of Attorney to any other person.

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This Power of Attorney shall have full force and effect from the date of its signing until otherwise revoked in writing.

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Paul Cornelis Verhulst

Paul Cornelis Verhulst

General Director cum Legal Representative of the Company

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For the authorized party

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A handwritten signature in blue ink, appearing to read "Tmt", with a horizontal line underneath.

Tran Thi Thanh Tra

Trần Thị Thanh Trà

Quality Assurance and Regulatory Affairs Lead of the Company

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Rebar™ Reinforced Micro Catheter

CAUTION

Federal (USA) law restricts this device to sale, distribution, and use by or on the order of a physician.

This device should be used only by physicians with a thorough understanding of angiography and percutaneous interventional procedures.

CONTENTS

One Rebar™ Micro Catheter with tip shaping mandrel

DESCRIPTION

The Rebar™ Micro Catheter is an endhole, single-lumen catheter designed to be introduced over a steerable guidewire into the vasculature. The proximal end of the catheter incorporates a standard luer adapter to facilitate the attachment of accessories. The catheter has a semi-rigid proximal shaft which transitions into the flexible distal shaft to facilitate the advancement of the catheter in the anatomy. Single or dual radiopaque markers at the distal end facilitate fluoroscopic visualization. The outer surface of the catheter is coated to increase lubricity.

INDICATIONS FOR USE

The Rebar™ Micro Catheter is intended for the controlled selective infusion of physician-specified therapeutic agents or contrast media into the vasculature of the peripheral and neuro anatomy.

CONTRAINDICATIONS

- It is not intended for use in the coronary vasculature.
- The Rebar™ Micro Catheter is contraindicated when, in the medical judgment of the physician, such a procedure may compromise the patient's condition.
- The Rebar™ Micro Catheter is contraindicated for neonatal and pediatric use.

POTENTIAL COMPLICATIONS

Potential complications include, but are not limited to:

- Hematoma at the site of entry
- Vessel perforation
- Vessel spasm
- Hemorrhage
- Pain and tenderness
- Vascular thrombosis
- Thrombotic episodes
- Neurological deficits including stroke and death
- Ischemia

WARNINGS

- Infusion pressure with this device should not exceed pressure as specified in the package insert. Excess pressure may result in catheter rupture, which may result in patient injury.
- If flow through the micro catheter becomes restricted, do not attempt to clear the device by high pressure infusion. Either remove the micro catheter to determine the cause of the obstruction or replace it with a new micro catheter. Excessive pressure may cause catheter rupture, which may result in patient injury.
- Never advance or withdraw an intraluminal device against resistance until the cause of the resistance is determined by fluoroscopy. Excessive force against resistance may result in damage to the device or vessel perforation.
- This device is supplied STERILE for single use only. Do not reprocess or re-sterilize. Reprocessing and re-sterilization increase the risks of patient infection and compromised device performance.

PRECAUTIONS

- Prior to use, carefully examine the Rebar™ Micro Catheter and its packaging to verify that it has not been damaged during shipment.
- Emboloc agents must be used in accordance with manufacturer's instructions.
- Prior to use, all accessory devices and agents should be fully prepared according to the manufacturer's instructions.
- Always monitor infusion rates when using this micro catheter.
- The micro catheter has a lubricious hydrophilic coating on the outside of the catheter. It must be kept hydrated in order to be lubricious. This can be accomplished by attaching a Y-connector to a continuous saline drip.
- When injecting contrast for angiography, ensure that the catheter is not kinked or occluded.
- When the infusion catheter is in the body, it should be manipulated only under fluoroscopy. Do not attempt to move the catheter without observing the resultant tip response.
- Because the Rebar™ may be easily advanced into narrow, selective vasculature, repeatedly verify (for example by slightly withdrawing the catheter) that the Rebar™ has not been advanced so far as to interfere with its removal.

STORAGE

This device should be stored in a dry place at between 50°F (10°C) and 90°F (32°C).

DIRECTIONS FOR USE

NOTE: See package insert for Rebar™ functional characteristics.

- Place the appropriate guide catheter following recommended procedures. Connect a one-way stop-cock to the guide catheter to prevent backflow of blood during insertion of catheter. Connect a hemostatic side arm adapter one way stop-cock to allow for insertion of catheter and to continuously flush the guide catheter with saline.
NOTE: It is recommended that a guide catheter with a minimum internal diameter of 0.053" be used with the Rebar™ Micro Catheter.
- Gently remove the packaging coil for the Rebar™ Micro Catheter from the pouch.

- Flush the coil with heparinized saline through the female luer attached to the coil.
- Flush the ID of the catheter with heparinized saline by attaching a saline filled syringe to the catheter hub.
- Remove the appropriate steerable guidewire from its packaging and inspect for damage. Follow manufacturer's instructions for preparing and using the guidewire.
- Carefully insert guidewire into the hub of the Rebar™ Micro Catheter and advance into the catheter lumen.
- Remove the micro catheter by removing the hub from the clip and gently pulling on the hub.
- Inspect the catheter prior to use to verify that it is undamaged. Retain the packaging coil for storage of the catheter when it is not in use during the procedure.
- Close the one-way stop-cock.
- Loosen the hemostatic valve.
- Introduce the guidewire and micro catheter as a unit through the hemostatic sidearm adapter into the lumen of the guiding catheter. Advance guidewire/catheter assembly to the distal tip of the guiding catheter.
- Tighten the valve around the micro catheter to prevent backflow, but allowing some movement through the valve by the micro catheter.
- Open the one-way stop-cock.
- Alternatively advance the guidewire and micro catheter until the desired site has been reached.
NOTE: To facilitate catheter handling the proximal portion of the micro catheter is uncoated to ensure a non-slip grip.
- When ready to infuse, withdraw the guidewire completely from the Rebar™ Micro Catheter. Connect a syringe containing infusate to the micro catheter hub, and infuse as required.
- After completing the procedure, withdraw the micro catheter and discard.

STEAM SHAPING MANDREL

In order to maintain catheter integrity and dimensional stability of the inner diameter, it is recommended that the user follow these instructions.

WARNINGS

Shaping Mandrel is not intended for use in the human body.

Use only a steam source to shape the catheter tip. Do not use other heat sources.

Prior to use, inspect the catheter tip for any damage that may have resulted from shaping. Do not use a catheter that has been damaged in any way. Damaged catheters may rupture causing vessel trauma or tip detachment during steering maneuvers.

- Remove shaping mandrel from card and insert into distal tip of the catheter.
- Carefully bend catheter tip and shaping mandrel to desired shape. A slight over exaggeration of the shape may be required to accommodate for catheter relaxation.
- Shape the catheter by holding the shaped portion approximately 1 inch (2.5 cm not less than 1 cm) from the steam source for about 20 seconds (DO NOT EXCEED 30 SECONDS.)
- Allow catheter tip to cool in air or saline prior to removing mandrel.
- Remove mandrel from catheter and discard. Multiple shaping is not recommended.
- Inspect the tip for any damage that may have resulted from steam shaping the catheter tip. If any damage is found, do not use the catheter.

FUNCTIONAL CHARACTERISTICS

Flow Rates @ 100 PSI (690 kPa) and 300 PSI (2,070 kPa) – ml/sec

Pressure	Saline		76% Contrast/ Saline (50/50)		60% Ionic Contrast		76% Ionic Contrast	
	100	300	100	300	100	300	100	300
Rebar™-10								
105-5078-153	0.26	0.68	0.10	0.38	0.04	0.16	0.03	0.09
105-5078-170	0.26	0.59	0.10	0.38	0.04	0.15	0.02	0.08
Rebar™-14								
105-5080-143	0.4	0.8	0.19	0.6	0.07	0.30	0.03	0.11
105-5080-153	0.3	0.6	0.12	0.4	0.06	0.18	0.03	0.09
105-5080-170	0.2	0.6	0.08	0.3	0.05	0.15	0.01	0.07
Rebar™-18								
105-5083-153	0.9	1.5	0.4	1.0	0.2	0.6	0.1	0.3
105-5081-153	0.9	1.5	0.4	1.0	0.2	0.6	0.1	0.3
105-5081-130	0.9	1.7	0.5	1.1	0.2	0.7	0.1	0.4
105-5081-110	1.0	1.8	0.6	1.3	0.3	0.7	0.1	0.4
Rebar™-27								
105-5082-145	1.4	2.5	0.6	1.8	0.3	1.3	0.2	0.7
105-5082-130	1.8	3.0	1.1	2.6	0.7	1.7	0.3	0.9
105-5082-110	2.0	3.4	1.3	2.7	0.9	1.9	0.3	0.9

All Rebar™ micro catheters have a maximum allowable dynamic pressure of 700 PSI (4,830 kPa) and static pressure of 300 PSI (2,070 kPa).

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Micro Therapeutics, Inc.