
**FAVOCCATH**

**I.V. CATHETER WITH INJECTION VALVE**

**INSTRUCTIONS FOR USE**

EN

**USM HEALTHCARE**  
**MEDICAL DEVICES FACTORY JSC.**  
 Lot L4b-13, N3 Street, Saigon Hi-Tech Park, Long Thanh  
 My Ward, Thu Duc City, Ho Chi Minh City 700000, Vietnam.  
 Tel: +84 28 730 44 888 · Fax: +84 28 730 43 888  
 Email: info@usm.com.vn · Website: www.usm.com.vn

**CMC Medical Devices & Drugs S.L.**  
 C/Horacio Lengo n16.C. P 29006. Málaga-Spain  
 Email: info@cmcmedicaldevices.com  
 Tel: +34 951 214 054

IF U-0001\_26-v01

**1 PRODUCT DESCRIPTION**

FAVOCCATH I.V. Catheter With Injection Valve has optimal design with 3-beveled needle reduces puncture strength through the skin, flexible material, durable, taper smoothly, kink resistant and non-irritating catheter, injection port, butterfly-wing design makes it easy to fix catheter and flashback chamber with anti-bacterial PTFE filter.

**2 INTENDED USE**

The device provides access to a vein of a patient in order to sample blood or administer fluids for the continuous use up to 96 hours.

**3 ADVANTAGES**

- Catheter consist of radiopaque stripes is composed by biocompatible and non-toxic material - Polyurethane which are flexible, durable and kink resistant. These attributes make FAVOCCATH I to an ideal choice for having access into the peripheral vascular system.
- Flashback chamber is provided with PTFE filter for anti-bacterial action and it enhances success in first-prick of needle.
- The taper smoothly catheter tip and closing needle body minimize the risk of damage to the vessel.

**4 INDICATIONS**

- Patients need to be recuperated from volume in circulatory system in case of diarrhoea, dehydration, severe burns.
- Patients need to be maintained several hours or several days a constant amount of drug in blood.
- Patients can not eat normally: coma, esophageal lesion, gastrointestinal lesion.
- Patients need to detox diuretic.
- The I.V. Catheter connects with the I.V. infusion line or the injection stopper in case of the suspending of infusion.

**5 CONTRAINDICATIONS**

Products should not be used in patients with known hypersensitivity to any of the materials used.

**6 PATIENT POPULATION**

Any patient population with consideration given to patient size, appropriateness for the solution being infused, and duration of therapy

**7 POTENTIAL COMPLICATIONS**

Insertion of a venous cannula is a painful procedure that can lead to anxiety and stress, infection, phlebitis, extravasation, infiltration, air embolism, hemorrhage (bleeding) and formation of a hematoma (bruise) may occur.

**8 INTRUCTION FOR USE**

- Select suitable size of FAVOCCATH and inspect to ascertain that package is intact.
- Confirm all items are accessible for FAVOCCATH insertion/ stabilization.

- Wash hands with antiseptic solution and wear gloves.
- Apply tourniquet.
- Locate a accessible and suitable peripheral vein visually and confirm by palpation.
- Disinfect the venipuncture site properly.

- Get FAVOCCATH out of the package.
- Remove needle cover in a straight outward motion by holding the wing body catheter hub and needle hub of FAVOCCATH.
- Ensure that injection port is upward direction prior to insertion.
- Do not touch the catheter.
- Hold FAVOCCATH by wing and needle body catheter, keeping needle beveled upward.
- Puncture the vein with needle at 35° angle and confirm successful venipuncture by visualizing blood in flashback chamber.
- Lower catheter at 15° angle.
- Push a part of catheter into the vein.
- Continue to push the catheter and simultaneously withdraw the needle body out slowly and continuously until the catheter is adequately into the vein.
- Discard needle into sharps container in trash.
- Remove tourniquet.
- Stabilize wing body catheter hub to prevent catheter dislodgement.
- Connect the catheter's hub with infusion line. If you need to suspend the process, use the injection stopper to close the hub.
- Cover and protect catheter sites by applying I.V. dressing or common methods of hospital.

**9 WARNINGS**

- Consult the instructions before use.
- Do not use after the expiration date.
- We disclaim any responsibility for possible consequences resulting from improper using.
- Re-sterilization may affect the structural integrity of the device and/or lead to damage device.
- Re-use causes contamination include, but is not limited to the transmission of infectious diseases from one patient to another. Contamination of the device may lead to injury, illness or death of the patient.

**10 PRECAUTIONS**

- Do not reinsert the completely or partially withdrawn needle into the catheter.
- Do not bend the catheter while using the product.
- Do not put high viscosity fluids on an I.V.
- It is not recommended to use this catheter for intra-arterial.
- This product should be used by trained medical staff only.
- Do not use if the package seems to be opened or damaged

**11 STORAGE AND TRANSPORTATION**

- Keep dry and any storage temperature 15°C - 30°C, protect against sunlight.
- The product should be transported with standard transport methods used for medical devices.

**12 SHELF LIFE**

The product shall conform to the specifications and functional requirements for a minimum of 5 years from the date of manufacture.

**13 DISPOSAL**

After use, the device should be disposed off in accordance with the local regulations in force.

**14 COMPATIBILITY WITH OTHER MEDICAL DEVICES**

FAVOCCATH is intended for use together with the following medical devices, but not limited to:

- Infusion Set.
- Saline Flush Syringe.
- Injection Stopper.

**15 SYMBOLS GLOSSARY**


Symbols	Explanation
	Manufacturer
	Date of manufacture
	Use-by date
	Catalogue number
	Batch code
	Sterilized using Ethylene oxide
	Do not re-sterilize
	Do not reuse
	Do not use if package is damaged
	Consult instructions for use
	Caution/WARNING/Danger
	Non-pyrogenic
	Keep away from sunlight
	Keep dry
	Fragile, handle with care
	Temperature limit
	CE Mark
	Authorized representative in the European Community
	Prescription Only

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