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
HS[®]
 FUTUREMADEPRESENT

General information:

The ALM SET is a device for the preoperative localisation of breast lesions with non-repositionable markers.

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Instructions for use:

- Remove the device from its packaging and remove the protective sheath.
- Prepare the intervention site.
- Make an incision with a scalpel if necessary.
- Prior to insertion, the marker must be visible just outside the needle's connector, indicating that the tip of the marker is inside the needle.
- Insert the needle into the lesion to be marked using an imaging technique.
- Advance the marker while keeping the needle fixed. The hook is fully positioned once the marks on the marker are in the cone of the needle.
- Make sure that the hooks are in the correct position before removing the needle.
- Medicate, fold, and place a bandage over the part of the marker that protrudes from the skin in order to keep it in place.

Warnings and precautions:

- This device is supplied sterile.
- Sterilised with ethylene oxide.
- Sterility is only guaranteed if the packaging is intact.
- Check the integrity of the packaging. Do not use the device if the packaging is not intact.
- Store in a cool, dry place away from any direct light sources.
- Integrity and sterility are only guaranteed if the specified storage requirements are met.
- Check the size and expiration date before opening.
- Always consider the possibility of allergic reactions.
- This device is intended for single use only. If re-used, the device may not function as expected, and could cause serious harm to the patient and the practitioner. The re-use of this device entails a risk of contamination. Do not re-use or re-sterilise the device.
- The device could pose a potential biohazard after use. The device must be handled and disposed of as infected hospital waste, in accordance with the current medical standards and the applicable local and national laws and regulations.
- The potential complications associated with preoperative localisation procedures are those which are clinically recognised, and include: infection, haemorrhage, damage to the surrounding tissues, and pneumothorax. If complications should arise, the practitioner should act accordingly based on his or her own expertise.
- Any accidents caused by this device must be promptly reported to the manufacturer and to the competent local authorities.
- Not for use in MRI.
- For use exclusively in hospital environments and sterile rooms.
- Never extract without surgical intervention if the marker has been positioned.
- The ALM SET S TYPE version's slim cone allows for insertion through the holes in the mammography machine's compression plate.
- The reinforced versions (identified by the letter S at the end of the product code) have a steel reinforcement near the hook.

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General information:

The PEIT NEEDLE is a device intended for the percutaneous alcoholisation of various organs.

It consists of an AISI 304 steel needle with a closed tip with 3 or 6 lateral holes positioned at 120°, and comes complete with a coupling with a Luer lock connector and an introducer needle. The tip's particular shape prevents unwanted needle deviations.

Careful medical assessments are required when considering PEI procedures on patients undergoing anticoagulant therapy or with blood clotting problems.

This device must only be used by experienced and properly trained physicians.

These instructions are not intended to establish or suggest any particular medical or surgical technique. The physician is responsible for the procedure and for this device's current usage techniques.

The treatments administered after a PEI may vary depending on the patient's condition.

CAUTION: Federal Law (USA) restricts this device to sale by or on the order of a physician.

Instructions for use:

- Remove the device from its packaging and remove the protective sheath.
- Prepare the intervention site.
- Identify the point at which to perform the PEI procedure, preferably under the guidance of ultrasound or another imaging technique.
- Infiltrate with a local anaesthetic, if necessary.
- The needle-guide can be used by positioning it on the ultrasound probe using an appropriate system and passing the needle through it.
- Insert the needle with the syringe attached to the coupling.
- Advance the needle into the target.
- Slowly inject the alcohol, checking its distribution on the ultrasound monitor.
- Extract the needle slowly.

Warnings and precautions:

- This device is supplied sterile.
 - Sterilized with ethylene oxide.
 - Sterility is only guaranteed if the packaging is intact.
 - Check the integrity of the packaging. Do not use the device if the packaging is not intact.
 - Store in a cool, dry place away from any direct light sources.
 - Integrity and sterility are only guaranteed if the specified storage requirements are met.
 - Check the size and expiration date before opening.
 - Always consider the possibility of allergic reactions.
 - This device is intended for single use only. If re-used, the device may not function as expected, and could cause serious harm to the patient and the practitioner. The re-use of this device entails a risk of contamination. Do not re-use or re-sterilise the device.
 - The device could pose a potential biohazard after use. The device must be handled and disposed of as infected hospital waste, in accordance with the current medical standards and the applicable local and national laws and regulations.
 - The potential complications associated with PEI procedures are those which are clinically recognised, and include: haematoma, haemorrhage, haemoperitoneum, haemobilia, chemical cholangitis, cholecystitis, abscesses, intestinal lesions, caval thrombosis, infection, localised pain, peritoneal irritation, hepatic insufficiency, biliary fistula, hepatic infarction, hypotension, renal insufficiency, and seeding. If complications should arise, the practitioner should act accordingly based on his or her own expertise.
 - Any accidents caused by this device must be promptly reported to the manufacturer and to the competent local authorities.
 - Not for use in MRI.
 - For use exclusively in hospital environments and sterile rooms.
-   The device's coupling (not in direct contact with the patient) is made from PVC containing phthalates (DEHP).

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