# Passeo-35

<b>English</b>	Peripheral Dilatation Catheter	Instructions for Use	5
<b>Deutsch</b>	Peripherer Dilatationskatheter	Gebrauchsanweisung	6
<b>Français</b>	Cathéter de dilatation périphérique	Mode d'emploi	7
<b>Italiano</b>	Catetere periferico per dilatazione	Istruzioni per l'uso	8
<b>Español</b>	Catéter de dilatación periférico	Instrucciones de uso	9
<b>Български</b>	Периферен дилатационен катетър	Инструкции за употреба	11
<b>Hrvatski</b>	Periferni dilatacijski kateter	Upute za uporabu	12
<b>Česky</b>	Katetr pro periferní dilataci	Pokyny k použití	13
<b>Dansk</b>	Perifert dilatationskateter	Brugsanvisning	14
<b>Nederlands</b>	Perifere dilatatiekatheter	Gebruiksaanwijzing	15
<b>Suomi</b>	Perifeerinen laajennuskatetri	Käyttöohjeet	16
<b>Ελληνικά</b>	Περιφερικός καθετήρας διαστολής	Οδηγίες χρήσης	17
<b>Magyar</b>	Perifériás tágitókatéter	Használati utasítás	19
<b>Latviešu</b>	Perifērais dilatācijas katetrs	Lietošanas norādījumi	20
<b>Lietuvių</b>	Periferinis plečiamasis kateteris	Naudojimo instrukcija	21
<b>Norsk</b>	Perifert dilatasjonskateter	Bruksanvisning	22
<b>Polski</b>	Obwodowy cewnik poszerzający	Instrukcja użycia	23
<b>Português</b>	Cateter de dilatação periférica	Instruções de utilização	24
<b>Română</b>	Cateter periferic de dilatare	Instrucțiuni de utilizare	25
<b>Русский</b>	Периферический расширяющий катетер	Инструкция по применению	27
<b>Slovenčina</b>	Katéter na periférnu dilatáciu	Návod na použitie	28
<b>Slovenščina</b>	Periferni dilatacijski kateter	Navodila za uporabo	29
<b>Svenska</b>	Perifer dilatationskateter	Bruksanvisning	30
<b>Türkçe</b>	Periferal Dilatasyon Kateteri	Kullanma Talimatı	31

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## Description

The Passeo-35 peripheral dilatation catheter is intended for dilatation of stenotic segments in peripheral vessels and arteriovenous dialysis fistulae. The dilatation balloon is designed to inflate to a known diameter at a specific inflation pressure consistent with the compliance chart on the label. One radiopaque marker is located at each end of the balloon to facilitate fluoroscopic visualization and positioning of the balloon catheter towards and across the lesion. The dilatation catheter includes a soft tapered tip to facilitate advancement of the catheter. The dilatation catheter has two Luer-parts at the proximal end. One port (inflation port) serves for connecting an inflation device to inflate/deflate the balloon. The other port enables flushing of the guidewire lumen. The dilatation catheter has a hydrophobic silicone coating on the shaft outer surface and on the guide wire lumen (inner surface), and a hydrophobic patchwork coating on the balloon. The dilatation catheter is compatible with guide wire and introducer sheath sizes according to the recommendations on the label.

**Caution:** In-vitro testing with the introducer sheaths 5F/6F Cordis Avanti+ (11 cm), 5F/6F Terumo Radifocus (10 cm), 5F/6F Terumo Destination (45 cm), 6F Cook Flexor Check-Flo Balkin (40 cm) and the BIOTRONIK-distributed 5F/6F Fortress from Contract Medical International (45 cm), has demonstrated that the Passeo-35 is compatible with the indicated minimum introducer sheath sizes. If Passeo-35 is used in conjunction with other, especially long and/or braided introducer sheaths, a larger French size than indicated on the label may be necessary to reduce friction.

## How Supplied

Sterile, Non-pyrogenic. Device is sterilized with ethylene oxide. DO NOT use if the package is opened or damaged, or if any information provided is obscured or damaged.

## Contents

- One (1) Passeo-35 peripheral dilatation catheter in a sealed, peel-open pouch.
- One (1) Instructions for Use Manual.

## Storage

Keep away from sunlight and store in a dry location between 10 °C and 40 °C / 50 °F and 104 °F.

## Indications

The Passeo-35 peripheral dilatation catheter is indicated to dilate stenosis in the renal, iliac, femoral, popliteal and intrapopliteal arteries and for the treatment of obstructive lesions of native or synthetic arteriovenous dialysis fistulae.

**Please note:** Passeo-35 peripheral dilatation catheter, does not have TGA (Therapeutic Goods Administration – Australia) approval for use within the renal and common iliac arteries.

## Contraindications

All general contraindications for percutaneous transluminal angioplasty (PTA) are contraindications for this device.

Contraindications for this device and peripheral dilatation catheters in general are:

- Lesions that cannot be reached or treated with the system
  - Large amounts of acute or sub-acute thrombus at the target lesion
  - Perforated vessels
  - Lesion that lies within or adjacent to an aneurysm
  - Uncorrected bleeding disorders
  - Renal insufficiency or an allergy to contrast media
- Furthermore, all procedure-related contraindications as described in the national and international guidelines of the respective medical associations apply.

## Warnings

- The Passeo-35 peripheral dilatation catheter is not indicated for use in coronary, cervical and intracranial arteries.
- This device is designed and intended for single use only. DO NOT resterilize and/or reuse. Reuse of single-use devices creates a potential risk of patient or user infections. Contamination of the device may lead to injury, illness or death of the patient. Cleaning, disinfection and sterilization may compromise essential material and design characteristics leading to device failure. BIOTRONIK will not be responsible for any direct, incidental or consequential damages resulting from resterilization or reuse.
- DO NOT use if the package is opened or damaged, or if any information provided is obscured or damaged.
- Use prior to the "use by" date as specified on the label.
- The inflated diameter of the balloon should never exceed the original diameter of the vessel proximal and distal to the lesion.

- Use only an appropriate balloon inflation medium, (e.g. 50:50 mixture by volume of contrast medium and saline). Never use air or any gaseous medium to inflate the balloon.
- DO NOT expose the catheter to organic solvents, e. g. alcohol.
- DO NOT use the balloon catheter if a vacuum cannot be held, as this indicates a system failure.
- When the dilatation catheter is in the body, it should only be manipulated whilst under high-quality fluoroscopy.
- DO NOT exceed the RBP stated on the label. The use of a pressure monitoring device is mandatory to prevent over-pressurization.

## Precautions

- Only physicians thoroughly trained and educated in the performance of percutaneous transluminal angioplasty (PTA) should use the dilatation catheter.
- Use only guide wires with a 0.035" [0.89 mm] diameter. Use only with introducers of an appropriate size as specified on the label.
- In-vitro testing with the introducer sheaths 5F/6F Cordis Avanti+ (11 cm), 5F/6F Terumo Radifocus (10 cm), 5F/6F Terumo Destination (45 cm), 6F Cook Flexor Check-Flo Balkin (40 cm) and the BIOTRONIK-distributed 5F/6F Fortress from Contract Medical International (45 cm), has demonstrated that the Passeo-35 is compatible with the indicated minimum introducer sheath sizes. If Passeo-35 is used in conjunction with other, especially long and/or braided introducer sheaths, a larger French size than indicated on the label may be necessary to reduce friction.
- Exercise care during handling to reduce the possibility of accidental breakage, bending or kinking of the catheter shaft.
- Prior to commencing the procedure, the dilatation catheter should be examined to verify functionality and ensure that its size is suitable for the specific procedure for which it is to be used.
- Precautions to prevent or reduce clotting should be taken. Flush or rinse all products entering the vascular system with sterile isotonic saline or a similar solution prior to use. The use of systemic heparinization during the procedure is recommended.
- If strong resistance is experienced during manipulation, stop the procedure and determine the cause of the resistance before proceeding.
- It is recommended to maintain vacuum whenever the dilatation catheter is subsequently advanced or withdrawn. Not maintaining the vacuum may result in difficulties pulling back the balloon through the introducer sheath.
- If the balloon has been dilated several times, there may be some resistance as it is pulled back into the introducer sheath. In case of difficulties, pull out the dilatation catheter and the introducer together.

## Adverse effects when used as intended

- Puncture site hematoma or bleeding
- Hemorrhage requiring transfusion or other treatment
- Hematoma (e.g. retroperitoneal)
- Injury to the vessel wall, dissection, perforation, rupture, intimal tear
- Vessel spasm
- Arteriovenous fistula
- Embolization of air, thrombotic or atherosclerotic material
- Pseudoaneurysm formation
- Total occlusion of the vessel
- Restenosis of the dilated vessel
- Thrombosis
- Failure to reach or cross the lesion
- Rupture or pinhole of the balloon
- Difficulties with inflation, deflation or withdrawal of the device
- Embolization of catheter material

- Local or systemic infection
- Allergic reactions to contrast media, antiplatelets, anticoagulants
- Emergency surgery to correct vascular complications
- Tissue necrosis and limb loss
- Death

Furthermore, all procedure-related adverse effects as described in the national and international guidelines of the respective medical associations apply.

## Directions for use

### Dilatation catheter preparation

1. Remove the dilatation catheter protection ring from the package and place it onto a sterile field.
2. Gently withdraw the dilatation catheter from the protection ring.
3. Carefully remove the balloon protector by pulling on the very distal end of the protector.

## Flush guide wire lumen

4. Connect a 10 ml or 20 ml syringe containing sterile saline to the Luer port of the guide wire lumen at the proximal end of the dilatation catheter.
5. Flush the guide wire lumen.
6. Remove the syringe.

## Purge air from catheter inflation lumen

7. Fill a 20 ml capacity inflation device with 10-12 ml of contrast solution.
8. Remove air from the inflation device according to the manufacturer's recommendations and instructions.
9. Connect a 3-way stopcock to the inflation port of the dilatation catheter and attach the inflation device to the stopcock.
10. Open the stopcock so that the fluid path between the catheter and the inflation/deflation device is established.
11. Pull back on the plunger and aspirate for 30 seconds until no bubbles appear inside the barrel of the inflation device during aspiration.

**Warning:** DO NOT use the balloon catheter if a vacuum cannot be held, as this indicates a system failure.

12. Close the stopcock so that the fluid path to the catheter is closed and evacuate all air from the inflation/deflation device through the stopcock.

13. Repeat steps 10-12 several times, if necessary, to ensure air contained in the balloon and inflation lumen are removed.

14. Return to neutral pressure and set the system aside for use.

## Insertion technique

15. Position the guide wire, under fluoroscopic guidance, in accordance with standard percutaneous transluminal angioplasty (PTA) techniques.
16. Insert the distal tip of the dilatation catheter onto the proximal end of the guide wire and advance until the guide wire exits the Luer lock at the proximal end of the catheter.
17. Carefully insert the dilatation catheter through the introducer.
18. Advance the dilatation catheter over the guide wire towards the lesion.
19. Position the balloon across the lesion using the balloon radiopaque markers as reference points.

## Balloon inflation/deflation

20. Ensure that inflation device is filled with 10-12 ml contrast solution.
21. Inflate the balloon with the inflation device to dilate the lesion using standard PTA techniques.

**Caution:** DO NOT exceed the RBP stated on the label.

22. To deflate the balloon completely, pull back the plunger of the inflation device and lock it in this position. Apply vacuum to the balloon under fluoroscopic control for at least 30-60 seconds depending on the balloon size.

**Caution:** It is recommended to maintain vacuum whenever the dilatation catheter is subsequently advanced or withdrawn. Not maintaining the vacuum may result in difficulties pulling back the balloon through the introducer sheath.

## Dilatation catheter removal

23. While maintaining stable guide wire position across the lesion withdraw the dilatation catheter carefully from the lesion and out through the introducer sheath. If the balloon cannot be withdrawn from the lesion easily, slightly advance and retract the device very carefully until it is possible to pull it out.

**Caution:** If the balloon has been dilated several times, there may be some resistance as it is pulled back into the introducer sheath. In case of difficulties, pull out the dilatation catheter and the introducer together.

24. Completely remove the dilatation catheter from the guide wire.
25. After use dispose the product and packaging in accordance with hospital, administrative and/or local government policy.

## Warranty/Liability

This product and each of its components (hereinafter product) were designed, manufactured, tested and packaged with all reasonable care. However, since BIOTRONIK does not have control over the conditions under which the product is used, the contents of this Instructions For Use (IFU) are to be considered as an integral part of this disclaimer for cases when a disturbance of the product's intended function may occur for various reasons.

BIOTRONIK does not guarantee that the following events will not occur:

- Product malfunctions or failures
- Patient's immune response to the product
- Medical complications during the use of the product or as a consequence of the product being in contact with the patient's body

BIOTRONIK does not assume any liability for:

- The use of the product that is not in accordance with the stated intended use/indication, contraindications, warnings, precautions and the directions for use of this IFU
- Modification to the original product
- Events which could not have been foreseen at the time of product delivery using the available levels of science and technology
- Events originating from other BIOTRONIK products or products not from BIOTRONIK
- Force majeure events including, without being limited to, natural disasters.

The above provisions shall be without prejudice to any disclaimer and/or limitation of liability agreed separately with the customer to the extent permitted by applicable laws.

For US only:

BIOTRONIK disclaims all warranties, expressed or implied regarding the product, including but not limited to, any warranty of merchantability or fitness for a particular purpose of the product. BIOTRONIK is not liable for any direct, incidental or consequential damages or expenses caused by any use of the product whether the claim is based on contract, warranty, tort or otherwise.

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