**[**REF**]** CE-919-553

## UF-Fluorocell™ SF

**Identification of the IVD reagent**  
UF-Fluorocell™ SF

**Intended use**

For in vitro diagnostic use only

UF-Fluorocell SF is used to mark formed elements in urine or in body fluids for determination of red blood cells, casts, etc. counts with the Sysmex fully automated urine particle analyzers.

UF-Fluorocell SF is intended for use only in conjunction with Sysmex diluent UF-CELLPACK SF.

**Principles of the examination method**

A sample volume of urine or body fluids is introduced into the analyzer where a portion is automatically diluted 1:4 with the special Sysmex diluent UF-CELLPACK SF and UF-Fluorocell SF. UF-CELLPACK SF reduces the interferences from amorphous salts and mucus and adjusts the condition of the specimen to enable marking of the formed elements. The entire dilution is maintained at a constant temperature for a defined time period in order to allow UF-Fluorocell SF to mark the formed elements in the specimen. The marked specimen is then introduced into the sheath flow detector where light scatter and fluorescent emission are measured, and red blood cells, casts, etc. counts are determined.

Refer to the analyzer's Instructions for Use for the measurement principle.

**Components**

Polymethine dye	0.05 <span> </span> %
Ethylene glycol	99.9 <span> </span> %

**Warnings and precautions**

1. Follow the warnings and precautions for handling described on the reagent container, package box, package insert or the analyzer's Instructions for Use, and use the reagent correctly. The reliability of the analysis values cannot be guaranteed if the reagent is used outside of the written intended use.

2. Do not use expired reagents, as the reliability of the analysis values cannot be guaranteed.

3. Do not use a reagent which you suspect has been frozen.

4. When replacing this reagent, do not refill and use the same container.

5. Handle the reagent with care to prevent air bubbles. If air bubbles are present, the performed analysis may not be correct.

6. Do not remove the RFID (Radio-Frequency Identification Device) tag until disposal. All the product information is managed by the RFID tag on the cartridge.

7. After setting, reset of the package is not recommended. Removing the reagent cartridge from the analyzer may cause deterioration of the reagent by contamination from the already opened sealing film.

8. Always wear protective coat and gloves when carrying out work on or with the reagent.

9. NEVER make contact of this reagent with the human body. Avoid contact with skin and eyes, and avoid ingestion. In case of skin contact, rinse immediately with plenty of water. In case of contact with eyes, rinse immediately with water or normal saline, occasionally lifting upper and lower eyelids until no evidence remains. Obtain medical attention. If swallowed, seek medical advice immediately.

Warning	Harmful if swallowed.	<div><div><div><div><div><span></span></div></div></div><div><div><div><span></span></div></div></div></div></div>
H302	P301+P312	IF SWALLOWED: Call a POISON CENTER/doctor if you feel unwell.
P501		Dispose of contents/container in accordance with local and national regulations.
16CFR1500 (U.S.A.):	WARNING	HARMFUL OR FATAL IF SWALLOWED.

**Additional required equipment**

UF-Fluorocell SF is intended for use exclusively only with the Sysmex reagents. If other reagents are used, the product performance of Sysmex instruments cannot be guaranteed.

Refer to the analyzer's Instructions for Use for further information.

**Examination procedure**

Put a UF-Fluorocell SF cartridge in the prescribed position and then connect the UF-Fluorocell SF. When using UF-Fluorocell SF, do not peel off the sealing film.

Refer to the analyzer's Instructions for Use for further information.

**Storage and shelf life of unopened product**

1. Store unopened UF-Fluorocell SF at 2-35 °C and avoid exposure to sunlight.

2. The expiration date of the unopened product is printed on the box and labels.

**Storage and shelf life after first opening**

Once opened and installed in the instrument, UF-Fluorocell SF is stable for 90 days.

**Performance characteristics**

Performance should be within instrument specifications.

Refer to the analyzer's Instructions for Use for performance characteristics.

**Limitations of the examination procedure**

Always use the reagent within the specified analysis environment for the instrument. Accurate analysis results may not be obtained if analyzed in an environment that does not meet the specified conditions. The results may be inaccurate due to the influences of foreign substances, bacteria growth, etc.

Refer to the analyzer's Instructions for Use for further information.

The use of this reagent, respectively control or calibrator, is validated on specific analyzers to optimize product performance and meet product specifications. Please refer to Instructions For Use of your analyzer whether the use of this reagent, respectively control or calibrator, is authorised by Sysmex. Sysmex cannot take the responsibility for patient results received from the use of Sysmex reagents, controls or calibrators on unauthorised analyzers. It is the responsibility of the user to validate modifications to these instructions or use of the reagent, control or calibrator on analyzers other than those specified by Sysmex.

**Primary sample collection, handling and storage**

UF-Fluorocell SF is intended for use with human urine and body fluids.

See the analyzer's Instructions for Use for further information regarding sample requirements.

**Disposal procedures**

1. Dispose of the cartridge after removing residual reagent solution and tightly seal the spout of the cartridge before disposing to prevent residual reagent solution leakage. You may use tape to secure the spout.

2. Disposal procedures of cartridge, reagent solution and waste effluents from the instrument should meet the requirements of applicable local regulations.

**Manufacturer**

	<b>Systemex Corporation</b> 1-5-1 Wakinohama-Kaigandori, Chuo-ku, Kobe 651-0073, Japan
---	---

**Authorized representatives**

Europe, Middle East and Africa:

	<b>Sysmex Europe GmbH</b> Bornbarch 1, 22848 Norderstedt, Germany
Americas:	<b>Sysmex America, Inc.</b> 577 Aptakistic Road, Lincolnshire, IL 60069, U.S.A.

Asia-Pacific:	<b>Sysmex Asia Pacific Pte Ltd.</b> 9 Tampines Grande #06-18, Singapore 528735
---------------	---

<b>Product information</b>		
UF-Fluorocell SF (UFF-800A)	29 mL x 2 cartridges	

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