QUẢ LỌC CF-609N

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
EN

ULTRA FILTER CF-609N

FILTER FOR CLEAN DIALYSATE AND WATER

Instructions for Use

Please read the instructions carefully before using the CF-609N.

How to use

1. How to connect

To connect the CF-609N, proceed according to the following instructions.


- (1) Connect the coupler onto the supply inlet and the filtered liquid outlet.
- (2) Remove the cap from the flushing port.
- (3) After clamping the filtered liquid outlet, prime the filter with dialysate or water from the supply inlet to the flushing port to release air from the filter housing.
- (4) After confirming that the air has been released from the filter housing, stop the flow of dialysate or water and attach the cap and clasp. Be sure not to let any air into the filter housing. If air does enter the filter housing, release the air according to the "flushing" instructions.
- (5) Remove the clamp on the filtered liquid outlet and pass dialysate or water. Continue the flow of dialysate or water for more than 10 minutes at a rate of 500 ml/minute to clean the filter.

* The CF-609N is designed for use with the external pressure filtration method (method in which filtration is carried out from the outside of the hollow fiber to the inside). Do not use with the vacuum pressure filtration method.

2. Flushing

To flush, proceed according to the following instructions.

- (1) Stop supplying dialysate or water.
- (2) Clamp the filtered liquid outlet and remove the clasp and cap attached to the flushing port.
- (3) Continue running dialysate or water through the unit for several minutes at a rate of 500ml/minute or higher. The CF-609N can be flushed with cleanser design for this system.
- (4) Stop supplying dialysate or water. Attach clasp and cap to the flushing port. Be sure not to let any air into the filter housing.
- (5) Remove clamp on the filtered liquid outlet and pass dialysate or water. If cleanser is used in flushing, rinse very carefully and be careful not to leave any residue of cleanser in the unit.

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Caution

1. Caution for Storage

- * Store at 0~35°C avoiding exposure to direct sunlight, severe vibrations and high humidity levels.
- * Do not unwrap until ready for use. Failure to observe this warning could result in contamination.
- * Be careful to avoid dropping or exposure to shock.

2. Caution in Use

- * The expiry date indicated on the labels and packages stands for the date by which the CF-609N must be used, as opposed to the period of use after installation on the system.
- * Do not use if the package has been damaged or has any pin holes.
- * The CF-609N is designed for use with the external pressure filtration method. Do not use with the vacuum pressure filtration method.
- * The CF-609N is a filter for high ultrafiltration and is designed under the assumption that processed dialysate or water will be permeated by the user. Therefore, if tap water is treated, the efficiency of the permeability will be rapidly reduced. If you have to use unprocessed dialysate or water, use a pre-filter to prevent the reduction of permeability.
- * If air enters the filter housing during use, release air from the flushing port.
- * By installing the CF-609N in the system, the flow rate or pressure will be reduced at the end of the system. Check each part after installation and confirm smooth operation.
- * Since loss of pressure may increase due to extended use, check periodically for irregularities.
- * To avoid contamination, make sure the CF-609N is thoroughly cleaned after installation.

Intended Use

The CF-609N is suitable for purifying dialysate and water by removing endotoxin, bacteria and particles.

Contraindications

There are no special contraindications when using the CF-609N for filtration of dialysate and water.

Guarantee

- * The CF-609N is manufactured under strict quality control and the quality is assured. However, if the CF-609N is defective (broken package, damaged product), it shall be replaced with a new one at our expense upon the return of the broken package or damaged product. We will, however, not be responsible for injury to patients or persons or damage to any object that is attributed to transport, storage and operation in your institution.
- * We will not be responsible for injury or damage if a patient or any person is injured or any object is damaged when using the CF-609N, unless we are clearly identified as being at fault.
- * We will not be responsible for any injury or damage caused by use of the CF-609N after the expiry date stated on the label or package.

Specifications & Performance

Hollow Fiber	Material	Polyethersulfone
	Membrane Area	0.6 m ²
	I. D.	500 μm
	Thickness	150 μm
Housing	Polycarbonate	
Potting Compound	Polyurethane	
O-ring	Silicone Rubber	
Maximum Operating Pressure	1.5 kgf/cm ²	
Pressure Loss(37 °C)	0.07 kgf/cm ² at 500ml/min	
	0.14 kgf/cm ² at 1000ml/min	
Cut-off Point (MW)	6000 (Protein)	
Sterilization	Gamma Ray	

NIPRO
NIPRO CORPORATION
3-9-3, Honjo-Nishi, Kita-Ku, Osaka 531-8510, Japan

EC REP NIPRO EUROPE N.V.
Weihoek 3H, 1930 Zaventem, Belgium

CE 0123

STERILE R

T.T.N.H.H
S