Item	Model Number	Description	Photo
HiRes™ Ultra Reusable Surgical Tool Kit	CI-4509	Kit containing all the reusable surgical templates required for the HiRes Ultra implant. It includes two (2) BTE templates, two (2) recess gauges, and two (2) coil gauges.	
HiRes Ultra BTE Template	CI-4421	The reusable template helps determine placement of the BTE sound processor on the patient's ear, which enables the surgeon to determine the location for the incision and implant placement.	
HiRes Ultra Recess Gauge	CI-4331	Reusable recess gauge used to size the bone bed for proper diameter and depth of the HiRes Ultra implant.	
			

HiRes Ultra Coil Gauge	CI-4341	Reusable metal gauge used for sizing of the skin flap or pocket for the HiRes Ultra implant.	
HiRes Ultra Mock Up and Literature Kit	CI-4426	The single-use HiRes Ultra silicone mock-up, provided sterile in a peel pack, is used to determine the proper fit and placement of the implant and electrode fantail. Please note: one (1) HiRes Ultra silicone mock-up is also included in the implant box.	
HiFocus SlimJ Depth Gauge	CI-1605	The HiFocus SlimJ electrode is a single-use, EO Sterilized depth gauge intended to be used to determine the patency of suspected occluded cochlea prior to insertion of the HiFocus SlimJ electrode.	
HiRes Ultra Magnet Replacement Tool Kit	CI-1418	The kit contains all reusable instruments needed to remove and replace the magnet.	

<p>HiRes Ultra Magnet Insertion Tool</p>	<p>CI-4348</p>	<p>Reusable tool used to remove the magnet from the silicon lip of the implant.</p>	
<p>HiRes Ultra Magnet Pusher Tool</p>	<p>CI-4349</p>	<p>The Magnet Pusher Tool holds the magnet while it is inserted under the silicon lip.</p>	
<p>HiFocus™ Mid-Scala Instrument Kit</p>	<p>CI-4508</p>	<p>The HiFocus Mid-Scala electrode instrument kit contains two (2) HiFocus Mid-Scala cochleostomy gauges and two (2) electrode claw tools.</p>	
<p>HiFocus™ Mid-Scala Cochleostomy Gauge</p>	<p>CI-4347</p>	<p>The cochleostomy gauge can be used to verify the minimum recommended cochleostomy or round window size.</p>	

<p>HiFocus™ Mid-Scala Claw Tools</p>	<p>CI-4254</p>	<p>The reusable Electrode Claw Tool can be used to stabilize the electrode during and/or following insertion or to push the array in during freehand insertion.</p>	
<p>90K™ Magnet Insert Dummy</p>	<p>CI-1410</p>	<p>The single-use insert dummy is used in cases where the HiRes 90K Magnet is removed for an extended period of time to allow for an MRI. The titanium dummy, provided in a sterile peel pack, is used to fill the cavity that houses the magnet. This prevents fibrous tissue growth in the cavity. The Magnetic Resonance Imaging (MRI) for the HiRes 90K brochure is included.</p>	
<p>HiRes™ Ultra Temporary Non Magnetic Plug</p>	<p>CI-1411</p>	<p>The single-use insert dummy is used in cases where the HiRes Ultra Magnet is removed for an extended period of time to allow for an MRI. The titanium dummy, provided in a sterile peel pack, is used to fill the cavity that houses the magnet. This prevents fibrous tissue growth in the cavity.</p>	

<p>HiRes™ Ultra 3D Temporary Non-Magnetic Plug</p>	<p>CI-1420</p>	<p>The single-use insert dummy is used in cases where the HiRes Ultra 3D magnet is removed for an extended period of time to allow for an MRI. The titanium dummy, provided in a sterile peel pack, is used to fill the cavity that houses the magnet. This prevents fibrous tissue growth in the cavity.</p>	
<p>HiRes 90K™ Replacement Magnet</p>	<p>CI-1412</p>	<p>The replacement magnet is provided in a sterile peel pack. It is used when a patient has had the original magnet removed for an MRI. The Magnetic Resonance Imaging (MRI) for the HiRes 90K brochure is included.</p>	
<p>HiRes™ Ultra Replacement Magnet</p>	<p>CI-1413</p>	<p>The replacement magnet is provided in a sterile peel pack. It is used when a patient has had the original magnet removed for an MRI.</p>	
<p>HiRes™ Ultra 3D Replacement Magnet</p>	<p>CI-1419</p>	<p>The replacement magnet is provided in a sterile peel pack. It is used when a patient has had the</p>	

		original magnet removed for an MRI.	
HiFocus Electrode Forceps Kit	CI-4350-02	It is two (2) reusable electrode forceps to aid in the insertion of the HiFocus SlimJ electrode.	
UHP Magnetic Clinic Pack	CI-7516	The UHP Magnet can be used in the internal and external wells of UHP and UHP 2.0, and ONLY the internal well of the AquaMic and Concave UHP. Add magnets as needed per the Increasing Magnet Strength charts using the minimum magnet strength required to maintain retention.	
Harmony™ Listening Check™ Earbuds	CI- 5822	Is a plastic headphone cord, including 2 speakers and 1 plug. Plug headphones into the audio processor	

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Digitally signed by Công Ty TNHH Sonova Việt Nam

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Reason: I am the author of this document

Location: your signing location here

Date: 2022.01.19 13:28:43+07'00'

Foxit PDF Reader Version: 11.1.0

AB

Control  
10/19/2021

# HiRes™ Ultra 3D Replacement Magnet

CI-1419

and

# HiRes™ Ultra 3D Temporary Non-Magnetic Plug CI-1420

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*Caution: Federal law restricts this device to sale, distribution and use by or on the order of a physician. For use in children, Federal law restricts this device to sale, distribution and use by or on the order of a physician who is trained in the pediatric cochlear implantation procedures for Advanced Bionics Cochlear Implants.*

## Symbols Glossary

This table provides an explanation of symbols that may be found on products and packaging

SYMBOL	SYMBOL MEANING OR EXPLANATORY TEXT
	European Community Mark of Conformity. Authorized to affix the CE Mark in 2018
	Caution (refer to "Cautions & Warnings" section for further information)
	See Instructions for Use
	Manufacturer
	Model Number
	Lot Number
	Sterile Lot
	Sterilized using ethylene oxide
	Use By Date
	Do not re-use
	Fragile
	Suitable temperature range for transport and storage
	Date of Manufacture
	MR Conditional: Specific circumstances are required for the cochlear implant recipient to undergo an MRI procedure with this product in place. See the HiRes Ultra 3D cochlear implant Instructions for MRI Safety Information for HiRes Ultra 3D cochlear implant for details.

## Introductory Information

Carefully read all instructions prior to use. Observe all warnings and precautions throughout these instructions.

## Cautions & Warnings

- Electrosurgical instruments (e.g. electrocautery) are capable of producing radio-frequency voltages of such magnitude that a direct coupling might occur between the cautery tip and the electrode. Induced currents may cause damage to the cochlear tissues or permanent damage to implants. Monopolar electrocautery must not be used. For the use of bipolar electrosurgical equipment, the probe tips must not contact the implant and should be kept more than 1 mm (0.04 in.) from the implant.
- Careful handling is required as the silicone pocket can be easily damaged by sharp instruments, forceps and similar instruments.
- If the sterile package is damaged, the product should not be used.

## Intended Use

The Ultra 3D Temporary Non-Magnetic Plug is intended to be installed in the HiRes Ultra 3D cochlear implant magnet pocket, temporarily, to avoid tissue growth in the pocket following the removal of the cochlear implant magnet.

The Ultra 3D Replacement Magnet is intended to be installed in the HiRes Ultra 3D cochlear implant magnet pocket as a replacement magnet following the removal of the Temporary Non-Magnetic Plug.

## Description of the Device

The HiRes Ultra 3D Replacement Magnet and Temporary Non-Magnetic Plug are non-reusable, single-use sterile components of the HiRes Ultra 3D cochlear implant. The Replacement Magnet and Temporary Non-Magnetic Plug fit within the magnet pocket of the HiRes Ultra 3D cochlear implant. Table 1 provides the specifications for both the Replacement Magnet and Temporary Non-Magnetic Plug.

TABLE 1. SPECIFICATIONS

	HiRes Ultra 3D Replacement Magnet (CI-1419)	HiRes Ultra 3D Temporary Non-Magnetic Plug (CI-1420)
		
<b>Material</b>	Titanium	Titanium
<b>Diameter</b>	12.7 mm	13.2 mm
<b>Thickness</b>	3 mm	3 mm
<b>Markings</b>	"3D"	None

The HiRes Ultra 3D Replacement Magnet and Temporary Non-Magnetic Plug are MR Conditional. Unilateral and bilateral recipients implanted with the HiRes Ultra 3D cochlear implant can be safely scanned with the magnet in place in a horizontal closed bore quadrature coil MR system meeting the following conditions:

MRI FIELD STRENGTH	1.5T	3.0T
Maximum Spatial Field Gradient	20 T/m	

The largest artifact at 3.0T is >15 cm when imaged under spin echo and gradient echo sequences. Smaller artifacts are possible if using different scanning parameters or a non-magnetic plug.

***For additional information regarding MRI see the MRI Safety Information for the HiRes™ Ultra 3D Cochlear Implant.***

## Operating Information

The HiRes Ultra 3D Replacement Magnet and Temporary Non-Magnetic Plug are ONLY compatible with the HiRes Ultra 3D cochlear implants (CI-1601) and the Universal Headpiece UHP 3D (CI-5316) and AquaMic™ 3D (CI-5317).

In order to remove and install the HiRes Ultra 3D cochlear implant magnet, Replacement Magnet and Temporary Non-Magnetic Plug you will need the following materials:

## Components Required

- HiRes Ultra 3D Replacement Magnet (CI-1419)
- HiRes Ultra 3D Temporary Non-Magnetic Plug (CI-1420)
- Recipient's Headpiece (CI-5316 or CI-5317)
- Marking Pen

## Magnet Removal / Replacement

The House Ear Elevator, Freer Elevator, Cerumen Curette, and forceps are common OR tools that can be used for magnet removal and replacement (Figure 1). Ensure that these or similar tools are sterilized and available for the magnet removal and replacement procedures.

*Figure 1. The OR tools should be sterilized prior to procedure.*

- Prior to preparing the sterile surgical site place the patient's headpiece over the coil / magnet portion of the implant.
- Verify that the headpiece magnetically attracts to the cochlear implant.
- Decide on location for an incision to gain access to the implant. This incision should be 1 centimeter away from the patient's headpiece either to the side of the implant or behind the implant (Figure 2).
- Once the location is determined draw the incision location with a marking pen.

*Figure 2. BTE and receiver package sites drawn with incision location.*

This location needs to remain visible after preparations for creating a sterile field including any hair removal.

***NOTE:*** *If marked location is not maintained for any reason, the patient's headpiece can be placed in a sterile bag and placed over the coil/magnet site and the incision location can be re-established.*

- Create an incision to gain access to the receiver package. Typically the fibrous capsule is accessed and the fibrous tissue lifted away from the silicone portion of the coil / magnet prior to cutting open.
- Open the fibrous capsule to gain visualization of the magnet pocket location on top of the receiver package coil / magnet region.
- Extend the skin incision and fibrous capsule as necessary to gain sufficient access to remove the magnet.

***NOTE:*** *Use a microscope to access and visualize the fibrous tissue, receiver package, and silicone aspects of the magnet pocket.*

### ***Magnet Removal***

- Ensure the magnet pocket site has sufficient fluid (water) to provide lubrication.

- Slide the distal end of the House Ear Elevator or similar tool under the far side of the magnet pocket lifting the silicone lip (Figure 3).


*Figure 3. Hold House Ear Elevator with angle facing downward.*

- Now slowly flip the magnet out of the pocket by tilting the tool on its side as shown. The magnet is attracted to the instrument so it is automatically retained on the tool (Figure 4).

*Figure 4. House Ear Elevator “flipping” magnet out of magnet pocket*

- Remove the magnet from the instrument by sliding off sideways (less magnetic attraction) (Figure 5).

*Figure 5. Magnet removed.*


### *Magnet Installation*

**NOTE:** The magnet is laser marked as shown with "3D" on the face. The "3D" marking may face in or out.

*Figure 6.  
Ultra 3D Magnet.*

*Figure 7.  
Place Freer Elevator  
under distal silicone lip.*

*Figure 8.  
Freer Elevator tool hold-  
ing magnet in position  
as House Ear Elevator is  
used to push the magnet  
into the silicon pocket.*

*Figure 9.  
Magnet Silicone Pocket  
Inspected – No Damage  
Verified (torn or split  
silicone is obvious).*

## Temporary Non-Magnetic Plug Installation

There are two options in which the magnet can be removed and replaced with a Temporary Non-Magnetic Plug. The first is prior to the implantation of the receiver package and the other is after implantation. Both scenarios are described below:

*NOTE: Since the Temporary Non-Magnetic Plug does not have a magnetic attraction there are no orientation marks on it. Either side can face up.*

### *Temporary Non-Magnetic Plug Installation Procedure Prior to Implantation*

Place the receiver package in the sterile field on a non-magnetic flat surface and ensure the device has fluid over the site of the magnet pocket region. Fluid provides lubrication to aid in magnet removal and Temporary Non-Magnetic Plug installation.

### *Temporary Non-Magnetic Plug Installation Procedure After Implantation*

Remove the magnet as previously described.

- The House Ear Elevator, Freer Elevator, Cerumen Curette, and forceps are common OR tools that can be used for insertion and removal of the Temporary Non-Magnetic Plug (Figure 10). Ensure that these or similar tools are sterilized and available for the Temporary Non-Magnetic Plug removal and replacement procedures.
- Obtain a sterile Temporary Non-Magnetic Plug from the sterile packaging.
- Ensure sterile water or sterile saline is in the magnet pocket and applied to the adjacent area.

*Figure 10. The OR tools should be sterilized prior to procedure*

- Place the Freer Elevator, House Ear Elevator, or Cerumen Curette under the silicone lip (between the Temporary Non-Magnetic Plug and silicone lip) and slide around in a circular fashion until the silicone lip is completely over the Temporary Non-Magnetic Plug (Figure 11).


*Figure 11. Inserting Temporary Non-Magnetic Plug using the Freer Elevator and forceps.*

- Use a gimmick (House Ear Elevator or a Cerumen Curette) to encourage the silicone lip over the remaining Temporary Non-Magnetic Plug (Figure 12). Use short 1-2 mm movements and proceed side to side, if needed.

*Figure 12. Using gimmick to manipulate the silicone lip.*

- Ensure no damage to the silicone lip prior to proceeding to closure (Figure 13).

*Figure 13. Magnet Silicone Pocket Inspected – No Damage Verified (torn or split silicone is obvious).*

## Handling of Sterile Products

The Replacement Magnet and Temporary Non-Magnetic Plug are supplied in ethylene oxide sterile packaging with indicators of sterilization. The sterile package should be inspected carefully to confirm it has not ruptured. Sterility cannot be guaranteed if the sterile package is damaged or opened. If the sterile package is damaged, the product should not be resterilized or used.

## Safe Disposal Instructions

All HiRes Ultra 3D magnets and Replacement Magnets should be returned to Advanced Bionics. All contaminated parts must be returned in the Cochlear Implant Explant Return Kit, in accordance with International Air Transport Association (IATA) rules. Please contact Advanced Bionics for further information.

All contaminated Temporary Non-Magnetic Plugs may be returned to Advanced Bionics in the Cochlear Implant Explant Return Kit, in accordance with International Air Transport Association (IATA) rules or be disposed of in accordance with local regulations for disposing of biohazardous materials.

## User Assistance Information

In the event that you experience any issues with your product, please contact Advanced Bionics. Do not attempt to service or modify your product as doing so may compromise system performance and will void the manufacturer's warranty. Products should be serviced only at Advanced Bionics.

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029-M785-49

029-M785-49 Rev A 2019-02-21

©2019 Advanced Bionics AG and affiliates. All rights reserved.

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029-M252-85

029-M252-85 Rev C

© 2016 Advanced Bionics AG. All Rights Reserved.

Instructions for Use  
English, Français, Español, Deutsch, Nederlands, Português

## **Warnings**


- The color cap could cause choking if swallowed.
- The color cap must be attached when the headpiece is in use.
- Replace the color cap every few months, or sooner, if a change in sound quality is noted.

## ***Intended Use***

The UHP and AquaMic™ Color Caps cover the headpiece magnet well and microphone.

## Placing the Color Cap


- Align the color cap with the headpiece.
- Lightly push the cap onto the headpiece until it snaps into place.


*UHP shown*

## Replacing the Color Cap

Removal Tool  
(CI-7514)

- Gently slide the removal tool between the color cap and the headpiece until the cap disengages.

## Cleaning

If necessary, the color caps can be cleaned with a slightly dampened cloth or tissue. Do not use soap or solvents.

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0°C (32°F) - 45°C (115°F)

### Labeling symbols and their meanings

	European Community Mark of Conformity. Authorized to affix the CE Mark in 2011		
	Model number		Manufacturer
	See Instructions for Use		Date of manufacture
	Lot Number		Warning
	Suitable for exposure to relative humidity between 0-95%		
	Store at temperatures between -20°C (-4°F) and +55°C (131°F)		
	Suitable for atmospheric range between 70kPa and 106kPa, which is equivalent to 3000m above sea level to 380m below sea level.		

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*AdvancedBionics.com*


029-M681-85

029-M681-85 Rev D 2018-06-26  
©2018 Advanced Bionics AG. All Rights Reserved.

AquaMic™ 3D  
Color Caps

CI-7130

instructions for use

English, Français, Español, Deutsch, Nederlands, Português

## Symbols Glossary

This table provides an explanation of symbols that may be found on products and packaging.

SYMBOL	SYMBOL MEANING OR EXPLANATORY TEXT
	European Community Mark of Conformity. Authorized to affix the CE Mark in 2018
	Caution (refer to "Cautions & Warnings" section for further information)
	See Instructions for Use
	Manufacturer
	Model number
	Lot number
	Fragile
	Date of Manufacture
	Suitable temperature range for transport and storage
	Suitable relative humidity range for transport and storage

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### Cautions and Warnings


- **CHOKING HAZARD:** *The color cap could cause choking if swallowed.*
- *Do not use broken, frayed, or cracked color caps.*
- *The color cap must be attached when the headpiece is in use.*

### Intended Use

The AquaMic™ 3D Color Cap is intended to cover the headpiece magnet well and microphone on your Aqua-Mic 3D headpiece.

### Description of the Device

The AquaMic 3D Color Cap is a detachable cover that snaps in place on top of the AquaMic 3D headpiece, securely covering the magnet well and microphone.

*AquaMic 3D Color Cap  
CI-7130*

*AquaMic 3D  
CI-5317*

*3D Color Cap Removal Tool CI-7525*

## Operating Information

The AquaMic 3D Color Cap is **ONLY** compatible with the AquaMic 3D headpiece (CI-5317). For instructions on how to change or replace the color cap on your headpiece, see the following illustration.

### *Attach*

### *Remove*

## Cleaning and Maintenance

Periodically, clean your color cap using a soft cloth. Do not use solvents. Replace the color cap every few months, or sooner, if a change in sound quality is noted.

## Operating and Storage Instructions

### Recommended Operating Conditions

Temperature: 0°C (32°F) - 45°C (113°F)

Relative Humidity: up to 95% RH

### Recommended Storage Conditions

Temperature: -20°C (-4°F) – 55°C (131°F)

Relative Humidity: up to 95% RH

## User Assistance Information

If you experience any issues with your product, please contact your cochlear implant professional or the manufacturer. Do not attempt to service or modify your product as doing so may compromise system performance and will void the manufacturer's warranty. Products should be serviced only at Advanced Bionics.

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*AdvancedBionics.com*


029-M666-85

029-M666-85 Rev D 2018-05-14  
©2018 Advanced Bionics AG. All Rights Reserved.

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## Introductory Information

Carefully read all instructions prior to use. Observe all warnings and precautions noted in these instructions.

### Cautions & Warnings


- **CHOKING HAZARD:** *The color cap could cause choking if swallowed.*
- *Do not use broken, frayed or cracked cables.*
- *The color cap must be attached when the headpiece is in use.*
- *Replace the color cap every few months, or sooner, if a change in sound quality is noted.*

## Intended Use

The UHP 2.0 Color Cap is intended to cover the headpiece internal magnet well and microphone on your UHP 2.0 Headpiece.

## Description of the Device

The UHP 2.0 Color Cap is a detachable cover that snaps in place on top of the UHP 2.0 Headpiece, securely covering the internal magnet well and microphone.

*UHP 2.0 Color Cap*  
*CI-7128*

*UHP 2.0 Headpiece*  
*CI-5315*


*UHP Removal Tool*  
*CI-7514*

## **Operating Information**

The UHP 2.0 Color Cap is **ONLY** compatible with the UHP 2.0 Headpiece (CI-5315). For instructions on how to change or replace the color cap on your headpiece, see the following illustration.

### Attaching and removing color cap:

## Cleaning & Maintenance

Periodically, clean your color cap using a soft cloth. Do not use solvents.

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[AdvancedBionics.com](http://AdvancedBionics.com)

029-M716-85

029-M716-85 Rev D 2018-06-27  
©2018 Advanced Bionics AG. All Rights Reserved.

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SYMBOL	SYMBOL MEANING OR EXPLANATORY TEXT
	European Community Mark of Conformity. Authorized to affix the CE Mark in 2018.
	Manufacturer
	Suitable temperature range for transport and storage
	Suitable relative humidity range for transport and storage
	Model number
	Lot number
	Fragile
	Date of manufacture
	See Instructions for Use
	Caution (refer to "Cautions & Warnings" section for further information)

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### Cautions & Warnings

- **CHOKING HAZARD:** *The color cap could cause choking if swallowed.*
- *Do not use broken or cracked color caps.*
- *The color cap must be attached when the headpiece is in use.*

## Intended Use

The UHP 3D Color Cap is intended to cover the headpiece magnet well and microphone on your UHP 3D headpiece.

## Description of the Device

The UHP 3D Color Cap is a detachable cover that snaps in place on top of the UHP 3D headpiece, securely covering the magnet well and microphone.

*UHP 3D Color Cap  
CI-7129*

*UHP 3D Headpiece  
CI-5316*

*3D Color Cap Removal Tool  
CI-7525*

## Operating Information

The UHP 3D Color Cap is **ONLY** compatible with the UHP 3D headpiece (CI-5316). For instructions on how to change or replace the color cap on your headpiece, see the following illustration.

### *Attach*

*Remove*

## Cleaning & Maintenance

Periodically, clean your color cap using a soft cloth. Do not use solvents. Replace the color cap every few months, or sooner if a change in sound quality is noted.

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Manufactured by:  
Advanced Bionics LLC  
California, U.S.A.  
+1 661 362 1400

[AdvancedBionics.com](http://AdvancedBionics.com)

029-M486-85 Rev A

© 2015 Advanced Bionics AG. All Rights Reserved.

# MRI Antenna Coil Cover

CI-7521

instructions for use

English, Français, Español, Deutsch, Nederlands, Português

## Labeling symbols and their meanings:

 REF	Model number		Date of manufacture
 LOT	Lot Number		Manufacturer
	Do Not Reuse		See Instructions for Use
	MR Conditional: Specific circumstances are required for the patient to undergo an MRI procedure with this device in place. See complete Instructions for Use for details		

## MRI Warnings

*MRI is contraindicated except under the circumstances described below. Do not allow patients with a HiRes 90K or HiRes 90K Advantage cochlear implant to be in the area of an MRI scanner unless the following conditions have been met:*

- The bandaging protocol recommended by Advanced Bionics is followed when the patient undergoes an MRI procedure with the magnet left in place or,*
- the internal magnet is surgically removed and possibly replaced with the Magnet Insert Dummy before the patient undergoes an MRI procedure.*

- *The external sound processor and headpiece are removed before entering a room where an MRI scanner is located.*
- *MRI preparations must be conducted outside of the MRI room.*
- *Verify that the implant, or both implants if bilaterally implanted, are compatible for conducting an MRI before proceeding. Failure to do so can lead to device movement, device damage, magnet movement, patient discomfort, or trauma and pain to the patient.*
- *MRI procedures are contraindicated for CLARION (C1 and CII) Cochlear Implant recipients.*

## Intended Use

The MRI Antenna Coil Cover is intended to be used along with a bandaging protocol to allow the patient to undergo an MRI procedure with the magnet left in place. Instructions for use must be followed, including the recommended bandaging protocol.

# MRI Antenna Coil Cover

## Compression Dressing prior to MRI Procedure with HiRes 90K Users

This document will provide general guidance for applying appropriate compression dressing over the site of the implanted device (and magnet) in order to perform an MRI procedure.

### Required Materials


- Hurt-Free Wrap – Hospital Grade  
Latex Free 5 cm width  
**Johnson and Johnson (or equivalent)**
- Flexible tape  
Latex Free  
2 – 2.5 cm width  
**3M – Nexcare (or equivalent)**
- Patient Headpiece
- MRI Antenna Coil Cover  
CI-7521  
**Advanced Bionics**
- Marking Pen

## Procedure

1. Place patient in sitting position to allow access to the implant site.


**NOTE:** *If the patient is bilaterally implanted, prepare to locate both implant devices using the following steps. Both implant locations must have compression dressing placed over their magnets.*

2. Obtain the patient headpiece (with cable removed) and place over the implant site. Allow magnetic attraction to occur.


3. Once in place, hold the headpiece while drawing a line around it.

4. Obtain the MRI Antenna Coil Cover and place it over the drawn location. Once centered, place flexible tape over the site to maintain the location.

5. Wrap the head using Hurt-Free Wrap sufficiently tight to ensure stabilization of the MRI Antenna Coil Cover.

### After the MRI

Once the MRI procedure has been completed, gently remove the Hurt- Free Wrap, tape, and the MRI Antenna Coil Cover.

Discard all components including the MRI Antenna Coil Cover.

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Manufactured by:  
Advanced Bionics LLC  
California, USA  
+1 661 362 1400

*AdvancedBionics.com*

029-M663-85

029-M663-85 Rev D 2017-11-01  
©2018 Advanced Bionics AG. All Rights Reserved.

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 European Community Mark of Conformity. Authorized to affix the CE Mark in 2018.	
 Model number	 Manufacturer
 See Instructions for Use	 Date of manufacture
 Lot Number	 Fragile
 Suitable temperature range for transport and storage	
 Suitable relative humidity range for trans- port and storage	
 Caution (refer to “Cautions & Warnings” section for further information)	

## *Introductory Information*

Carefully read all instructions prior to use. Observe all cautions and warnings noted throughout these instructions.

### *Cautions & Warnings*

- **CHOKING HAZARD:** The magnet could cause choking if swallowed.
- Seek medical attention if any magnets are swallowed.
- Use the minimum number of magnets required to ensure sufficient magnetic attraction of the headpiece.
- Before placing a magnet in the external well, be sure to remove additional internal magnets to avoid potential skin irritation, pressure, or redness.
- If the user experiences any discomfort or pain they should remove the headpiece immediately and contact their cochlear implant professional.

- If the headpiece magnet is too strong or is in contact with the skin, pressure sores may develop at the headpiece site. If this happens, or if the user experiences any discomfort in this area, they should contact their cochlear implant professional.
- Store spare magnets safely and away from cards that may have a magnetic strip (e.g. credit cards).

## *Intended Use*

These headpiece magnets are intended to provide retention between the headpiece and the implant of an Advanced Bionics Cochlear Implant System.


## Description of the Device

*CI-7522*  
*N52 Magnet*

*CI-7511*  
*UHP Magnet*

*CI-7518*  
*Thin Magnet*

## *UHP Magnet*

The UHP Magnet can be used in the internal and external wells of UHP and UHP 2.0, and **ONLY** the internal well of the AquaMic and Concave UHP.

## *N52 Magnet*

The N52 Headpiece Magnet should **ONLY** be used in the internal well of the headpiece if there is insufficient retention between the user's headpiece and cochlear implant after using four headpiece magnets in the internal well and one headpiece magnet in the external well (as shown in the Increasing Magnet Strengths charts).

## *Thin Magnet*

The Thin Magnet should **ONLY** be used in the external well of the AquaMic and Concave UHP.


## *Operating Information*

The following charts explain how to ensure adequate retention between the headpiece and cochlear implant in each headpiece by gradually adding magnets to increase magnet strength.

## UHP and UHP 2.0:

Increasing Magnet Strength →

### Internal Well


### External Well

## AquaMic and Concave UHP:

Increasing Magnet Strength →

### Internal Well

### External Well

Magnet Key	
	UHP Magnet
	Bottom Cover Patch
	Thin Magnet

## To Change Magnets

1. Remove the color cap with the removal tool.
2. Remove the foam spacer from the internal magnet well (if present).
3. Remove the magnets from the internal well. If the external magnet is present, remove this first with the Magnet Removal Tool.
4. Add magnets as needed per the Increasing Magnet Strength charts using the minimum magnet strength required to maintain retention.

***NOTE: If retention issues remain after using the maximum number of magnets per the Increasing Magnet Strength charts, remove all internal well magnets and replace them with an N52 Magnet. If retention issues still persist following use of the N52 Magnet, contact Advanced Bionics Customer Service.***

5. Add the Bottom Cover Patch if no magnet is used in external well.

6. Place a foam spacer over the magnets in the internal well only if the spacer can remain flush with the top of the well. Otherwise, do not install the foam spacer.
7. Install the color cap on the headpiece.

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Universal Headpiece  
(UHP)  
CI-5305

**instructions for use**

English, French, Spanish, German, Dutch, Portuguese

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[AdvancedBionics.com](http://AdvancedBionics.com)

029-M298-85 Rev A  
©2013 Advanced Bionics AG. All Rights Reserved.


UHP

Spacer

CI-7517


UHP

Magnet

Color Cap/Magnet Removal Tool

Bottom Cover Patch

CI-<sup>1</sup>

7513


 **Warnings**

- The UHP contains small parts that could cause choking if swallowed.
- Do not use the UHP if the housing is cracked.
- Do not use broken, frayed or cracked cables.
- Color Cap must be attached when headpiece is in use.
- Use the minimum number of magnets required to ensure sufficient magnetic attraction of the UHP.
- If irritation, pressure, or redness occurs contact your cochlear implant center.

Generally, additional magnets will not be required. If additional magnets are needed, they should only be inserted under the direction of a cochlear implant professional. The UHP Bottom Cover Patch should only be removed by a cochlear implant professional. If the patch is dislodged, replace it with the extra patch provided, or contact your cochlear implant center.

### Adding Magnets


If needed, additional magnets (\*) may be added one at a time within the UHP (Figure 1). If UHP retention is inadequate with four internal magnets, the Bottom Cover Patch (\*\*) may be removed, and one UHP magnet may be added to the external well (Figure 2).


**⚠ Warning:** Before placing a magnet in the external well, be sure to remove additional internal magnets to avoid potential skin irritation, pressure, or redness.

## Cable\*

To avoid damage to the cable when removing from the UHP, grasp only by the connector.


## Placing the Color Cap\*




\*packaged separately

## Removing the Color Cap

## Replacing the Bottom Cover Patch

Peel the paper backing off of the patch and press into place adhesive side down.

## Cleaning

- Use a soft cloth to clean. Do not use solvents.

**Recommended operating temperature:**

0°C (32°F) - 45°C (115° F)

**Labeling symbols and their meanings:**

 0123	European Community Mark of Conformity. Authorized to affix the CE Mark in 2011		
REF	Model number	SN	Serial number
	See Instructions for Use	 2000-00	Date of manufacture
	Manufacturer		Do not get wet
	Fragile		Warning

	Type of protection: B
	Dispose of in accordance with applicable national and local regulations
	Suitable for exposure to relative humidity between 0-95%
	Store at temperatures between -20°C (-4°F) and +55°C (131°F)
	Suitable for atmospheric range between 70kPa and 106kPa, which is equivalent to 3000m above sea level to 380m below sea level.

## Compatibility

The following tables display the compatibility between products in the HiResolution Bionic Ear System family and/or previous generation product.

Table 1

		Implant Type			
Processor Type		C1	CII	HiRes 90K	HiRes 90K Advantage
	<i>Neptune</i>	-	✓ <sup>1</sup>	✓ <sup>1</sup>	✓ <sup>1</sup>
	<i>Harmony</i>	✓ <sup>2</sup>	✓ <sup>3</sup>	✓ <sup>3</sup>	✓ <sup>1</sup>
	<i>Auria</i>	-	✓	✓	✓ <sup>1</sup>
	<i>Platinum Sound Processor (PSP)</i>	✓ <sup>4</sup>	✓	✓	✓ <sup>1</sup>

<sup>1</sup> Requires SoundWave 2.1 or later

<sup>2</sup> Requires SoundWave 2.0 or later

<sup>3</sup> Requires SoundWave 1.4 or later

<sup>4</sup> Requires SCLin2000


Table 2

		Implant Type			
Software Type/Version		C1	CII	HiRes 90K	HiRes 90K Advantage
	SClin2000	✓	✓	-	-
	SoundWave (versions 1.x)	-	✓	✓	-
	SoundWave 2.0	✓ <sup>1</sup>	✓	✓	-
	SoundWave 2.1 or later	✓ <sup>1</sup>	✓	✓	✓

<sup>1</sup> Only on Harmony

Table 3

		Processor Type			
Headpiece Type		<i>Neptune</i>	<i>Harmony</i>	<i>Auria</i>	<i>PSP</i>
	<i>Universal Headpiece (UHP)</i>	✓	✓	✓	✓
	<i>AquaMic</i>	✓	-	-	-
	<i>HR 90K Auria Headpiece</i>	-	✓	✓	-
	<i>Platinum Headpiece</i>	-	✓	✓	✓

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

029-M717-85

029-M717-85 Rev E 2018-07-18  
©2018 Advanced Bionics AG. All Rights Reserved.

# CI-7527

3D Magnets  
CI-7523

3D Headpiece  
Foam Spacer  
CI-7524

instructions for use

English, Français, Español, Deutsch, Nederlands, Português

## Symbols Glossary

*This table provides an explanation of symbols that may be found on products and packaging.*

SYMBOL	SYMBOL MEANING OR EXPLANATORY TEXT
	European Community Mark of Conformity. Authorized to affix the CE Mark in 2018.
	Caution (refer to "Cautions & Warnings" section for further information)
	See Instructions for Use
	Manufacturer
	Model Number
	Lot Number
	Fragile
	Date of Manufacture
	Suitable temperature range for transport and storage
	Suitable relative humidity range for transport and storage

## *Introductory Information*

Carefully read all instructions prior to use. Observe all warnings and precautions noted throughout these instructions.

### *Cautions & Warnings*

- **CHOKING HAZARD:** The magnet could cause choking if swallowed.
- Seek medical attention if any magnets are swallowed.
- Use the lowest strength magnet required to ensure sufficient magnetic attraction of the headpiece.
- If the magnetic attraction of the headpiece is insufficient the user should contact their cochlear implant professional.
- If the headpiece magnet is too strong, pressure sores may develop at the headpiece site. If this happens, or if the user experiences any discomfort, irritation, pressure, or redness in this area, they should contact their cochlear implant professional.

- Store spare magnets safely and away from cards that may have a magnetic strip (e.g. credit cards).



## *Intended Use*

The 3D Magnets are intended to be used with Advanced Bionics' UHP 3D or AquaMic™ 3D headpieces.

## *Description of the Device*

The 3D Magnet is a removable magnet that is available in four different configurations, each of which offers different strengths to ensure retention between the cochlear implant and the headpiece. A foam spacer is affixed to the magnet to ensure a secure fit inside the headpiece magnet well. Arrows on top of the magnet indicate the magnet strength and help to orient the magnet in the general upward direction within the headpiece. An additional 3D Headpiece Foam Spacer is provided to be used with the smallest strength 3D Magnet (one arrow). The 3D Color Cap Removal Tool includes a magnetic bottom to assist with the magnet placement in the headpiece and a tip to aid in

the removal of the headpiece color cap.

## Operating Information



The 3D Magnets are ONLY compatible with the UHP 3D (CI-5316) and AquaMic 3D (CI-5317) headpieces and the HiRes Ultra 3D cochlear implants (CI-1601). Instructions on how to change or replace the headpiece magnet are provided in this booklet.

## Magnet Strength

The four different configurations of the 3D Magnets offer different retention strength based on the orientation of the magnet. Magnets can be flipped over should you need the spacer closer or farther from the implant to make retention strength stronger or weaker.

*NOTE: If using the smallest strength 3D Magnet (one arrow), add the 3D Headpiece Foam Spacer in the magnet well.*

The following chart details these strengths:

## *Adjusting Magnet Strength*

To adjust the magnet strength of the headpiece you will need to:

1. Remove the color cap using the 3D Color Cap Removal Tool.
2. Remove the magnet from the magnet well.
3. Replace magnet as needed per the Magnet Strength chart using the minimum magnet strength required to maintain retention.
4. Orient the arrow on the magnet in a general upward direction using the 3D Color Cap Removal Tool.
5. Attach the color cap.

If you note retention issues, it may be because the arrow on the magnet is not oriented in the recommended direction. If retention issues persist following use of the highest strength magnet option, contact Advanced Bionics Customer Service.

## ***Cleaning & Maintenance***

Use a soft cloth to clean magnets. Do not use solvents.

## ***Operating & Storage Instructions***

### ***Recommended Operating Conditions***

Temperature: 0°C (32°F) – 45°C (113°F)

Relative Humidity: up to 95% RH

### ***Recommended Storage Conditions***

Temperature: -20°C (-4°F) – 55°C (131°F)

Relative Humidity: up to 95% RH

## ***User Assistance Information***

The 3D magnets may be ordered either by individual model number or as a 3D Headpiece Magnet Kit (CI-7527). If you experience any issues with your product, please contact the manufacturer. Do not attempt to service or modify your product as doing so may compromise system performance and will void the manufacturer's warranty. Products should be serviced only at Advanced Bionics.

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