INSTRUCTIONS FOR USE LEAFLET FOR HUBER NEEDLES

DESCRIPTION

Huber needles are available straight or curved in different lengths and gauges depending on their use. They are packaged in boxes of 50 units.

INDICATIONS

Huber needles are intended for administration or removal of fluids from implantable ports. (IP).

CONTRA-INDICATIONS

On the skin against the septum:

- erythema, exanthema,
- unusual pain,
- weeping, discharge of fluid, fistula, abscess, any sign or suspected local infection,
- any sign or suspected turning around of the reservoir,
- difficulty or inability to locate the septum of the port.

POSSIBLE COMPLICATIONS

The complications most frequently described with the use of Huber needles are:

- pain following puncture or during chemotherapy treatment,
- infection,
- turning around, expulsion with or without extravasation,
- insufficient flow or needle obstruction,
- cutaneous fistula due to repeated punctures at the same skin site,
- distal buckling of the needle (being thrust into the bottom of the reservoir or crushing due to compression under the dressing).

WARNINGS AND PRECAUTIONS

- This product does not protect against leakage of contaminated fluids.
- Single use** medical device. **DO NOT RE-USE**. A damaged needle or attempt to re-introduce needle could cause an accidental needlestick injury with a contaminated needle.
- Sterilised with ethylene oxide. **DO NOT RESTERILISE**.
- Do not use the device if the packaging is damaged, opened or after the expiry date**.
- Implantable port puncture: move slightly up from the previous puncture area in order to avoid leakage or a fistula between the septum and the patient’s skin.
- Pain during the infusion: stop the infusion immediately after seeking medical device and review.
- Check that the needle length is suitable for the port and thickness of the patient’s subcutaneous tissue. If it is too long the needle or the port can be damaged during puncture. If it is too short the needle cannot completely puncture the septum and the drug solution may be administered into surrounding tissue and/or the needle may be blocked.
- Do not press the needle tip too firmly on the bottom of the port (risk of buckling the needle and damaging the septum during removal).

USE

Observe hygiene and aseptic technique when using Huber needles. Prepare the implantation site to insert the sterile needle.

Setting up

- Follow the steps of cutaneous aseptic technique before introducing a sterile needle**, according to the hospital protocol.
- Select the length of needle according to the patient’s morphology: curved lengths: 25-30-38 mm; lengths available for curved needles: 20-25-30-35 mm). The 35 mm and 38 mm lengths are reserved for patients with corpulent body morphology. Selecting a length of needle which is not suitable for the patient’s body morphology may cause complications:
 - too short: partial or total reduction in flow or self expulsion of the needle from the septum which may cause e extravasation (necrosis)
 - too long: rocking effect which may result in the needle becoming removed from the septum (extravasation) or poor fixation of the needle onto the patient’s skin.
- Select the gauge** according to the viscosity of the drug to be infused, desired flow rate and the **thickness of the patient’s** subcutaneous skin. 23 Gauge is reserved for paediatric use.
- Purge the device with physiological saline
- Hold the Huber needle so that the **body of the needle is perpendicular to the surface of the skin**.
- Pass the needle completely through the septum, not crushing the tip against the bottom of the port. The nurse should feel the tip of the needle touching the bottom of the port (an indication that the septum has been completely fully penetrated).
- Once the needle is in place check its correct position and that the system is patent (**checking for blood reflux**).
- Fix the Huber needle firmly with a transparent occlusive dressing preferably with a central non-adhesive POLYFILM® area.
- Warning: Do not press the needle too hard against the bottom of the septum otherwise the tip of the needle may become buckled with a risk of damaging the septum when the needle is removed.
- Monitor regularly throughout the infusion.

Frequency of replacing

This product must be replaced in accordance with the instructions from the Centres for Disease Control or CDC, USA) for administration devices, with specific local and national regulations, standard professional practice and/or the protocols in your establishment.


Pursuant to the centres for disease control (CDC, USA) instructions, in order to prevent infections associated with intravascular catheters (9 August 2002/51(RR10); 1-26), replace the administration devices including the secondary and ancillary devices at intervals of at most of 72 hours , unless catheter-related infection is suspected or found.

Removal

- After a final rinse with physiological saline, remove the needle applying positive pressure (continue to inject physiological saline as the needle is removed)
- Compress the puncture area for one or two minutes, cleaning the skin and applying a dry dressing left in place for 24 hours.

Discarding of used device

Needles and other cutting instruments must be disposed of in appropriate sharps boxes supplied in order to avoid any risk to staff. Needles and syringes must not be re-sheathed before disposal. Used devices must be disposed on in accordance with the CDC (USA) instructions and/or the protocols in your establishment for protection against blood pathogens.

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