## POLYPERF®

### CONNECTING LINES POLYPERF® WITH HUBER NEEDLE

#### Instructions for use

##### DESCRIPTION

POLYPERF® are connecting lines with a curved Huber needle. POLYPERF® are available in different lengths and gauges and with or without lateral injection site. There are provided with luer lock connectors and clamps. POLYPERF® tubing is DEHP free (without di-2 ethylhexyl phtalate). There are packaged in box of 10 units.

##### INDICATIONS

POLYPERF® connecting lines are indicated for the administration into or withdrawal of fluids from implanted ports.

##### CONTRA-INDICATIONS

On the skin, presence of:

- Erythema, exanthema
- Unusual pain
- Oozing, discharge of liquid, fistula, abscess, all sign or suspicion of local infection
- All sign or suspicion of reservoir turning round
- Difficulty or impossibility to locate the septum and the port.

##### POSSIBLE COMPLICATIONS

The complications most frequently described, linked with the use of Huber needles are the following:

- Pain during the puncture or during a chemotherapy treatment
- Infection
- Instability, turning round, expulsion with or without extravasation
- Restriction of flow, or needle obstruction
- Cutaneous fistula due to repeated punctures at the same cutaneous area
- Distal bending of the needle (thrust in the bottom of the reservoir, or crushing due to a compression under the dressing)
- Difficulty or impossibility to perform positive pressure during the needle removal
- Accidental Blood Exposure during the removal or destruction of the contaminated needles.

##### WARNINGS AND PRECAUTIONS

- Single use medical device. **DO NOT REUSE**. The damaging or the reintroduction attempt of the needle could conduct to an accidental pricking with a contaminated needle.
- Sterilized with ethylene oxide. **DO NOT RESTERILIZE**.
- Do not use the device if the packaging is damaged, opened or after the expiry date.**
- Fully tighten Y site injection before each use.
- Failure to attach an appropriate needless access device after removing the injection cap can result in air embolism.
- Implantable port puncture: slightly move up from the last puncture area in order to avoid sluice or fistula between the septum and the patient skin.
- Pain during the treatment: immediately stop the perfusion with medical advice and control.
- Check that needle length suits with the port and patient subcutaneous thickness: if it is too long, the needle or the port can be damaged during the puncture; if it is too short, the septum can not be entirely punctured and the medication can be administered in the surrounding tissues and/or the needle can be blocked.
- Avoid the contact between the Huber needlepoint and the port bottom (risk of septum damage during the removal).
- When administering light sensitive solutions, the tubing must be protected from light.
- POLYPERF® is not intended to be used with high-pressure infusion pumps.

***Caution:** Federal law (FDA) restricts this device to sale by or on the order of a physician*

##### STEP-BY-STEP PROCEDURE

**Respect the hygiene and asepsis measures during the Huber needles use. Prepare the portal site for sterile needle insertion.**

##### Setting up

- Select the gauge** according to the viscosity of the medication to be infused, to the flow rate required, and **the suitable length** according to the patient subcutaneous thickness.
- Flush the device. **Use of luer lock Y site:** if appropriate, remove injection cap from luer lock Y site, attach an appropriate needless access device and prime the Y site. **Or**, insert a needle or blunt cannula into septum in the center of the injection cap and prime the Y site.

##### Before each use of the injection site

- Laboratoires Perouse recommend the following validated disinfecting protocol according to the FDA Guidance on Premarket Notifications For Intravascular Administration Sets (April 15, 2005): Disinfect the surface with cleansing swabs containing 70% Isopropyl Alcohol. Drying.
- Remove the needle guard and insert the needle into the septum of the implanted port. Start injection after checking the correct position and the good permeability of the system.
- Apply of a sterile occlusive dressing allowing a perfect visibility of the puncture area. Regular checking during the injection

##### Frequency of replacing

This product must be changed in accordance with U.S. Centers for Disease Control (CDC, USA) guidelines for administration sets, local or country specific guidelines, professional standards of practice, and /or according to your institution's polly. According to CDC Guidelines (USA) for the Prevention of Intravascular Catheter-Related Infections (August 9, 2002 / 51(RR10);1-26) replace administration sets, including secondary sets and add-on devices, no more frequently than at 72-hour intervals, unless catheter-related infection is suspected or documented

##### Removal

- After a last rinsing with physiological serum, remove the needle with positive pressure (keep on injecting during the needle removal)
- Compression of the puncture area during 1 or 2 minutes, skin cleaning and applying of a dry dressing during 24 hours.

##### Discarding of used device

Needles and other sharps must be disposed in the appropriate needle boxes provided to avoid danger to other personnel. Needles and syringes are not recapped prior to disposal. Used device must be disposed in accordance with CDC guidelines (USA) and / or your institution's blood borne pathogens program.

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