EUROSCREW®NG / EUROSCREW®TCP NG

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EUROSCREW®NG / EUROSCREW®TCP NG

Bioabsorbable interference screw Operator information

Made in France by :
TEKNIMED S.A.S
8, rue du Corps Franc Pommies
65500 VIC EN BIGORRE - FRANCE
T l. (33)5 62 96 88 38
Fax (33)5 62 96 28 72
www.teknimed.com

Date CE marking first applied: 2014

Before using Teknimed products, the user should pay careful attention to the safety recommendations included in the instructions for use, as well as information relating to each product (description, operating technique, brochure, etc.). The appropriate information can be obtained from TEKNIMED.

GENERAL INFORMATION

Teknimed products must be implanted by qualified operators having in-depth knowledge and complete understanding of operating techniques specific to Teknimed products. The operating techniques can be obtained from distributors. The surgeon is responsible for complications or harmful consequences that could result from the wrong indication or operating technique, incorrect use of the equipment and a failure to observe safety recommendations given in the instructions for use. These complications may be attributed neither to the manufacturer nor to the relevant Teknimed representative. The EUROSCREW NG medical device is a bioresorbable interference screw for ligament surgery.

COMPOSITION

EUROSCREW NG:
100% poly(70/30; L/DL)lactide.
EUROSCREW TCP NG :
70% poly(70/30; L/DL)lactide + 30% TCP.

INDICATIONS

EUROSCREW NG screws are recommended for ligament attachment, for example in anterior/posterior cruciate ligament reconstructions, reinsertion of the long-head biceps or ankle ligaments.

PRECAUTIONS FOR USE

- Surgeons must have a thorough knowledge of surgical procedures and, after instruction in a ligament fixing technique using interference screws, abide by this rigorously.
- For the proper use of the EUROSCREW® NG and EUROSCREW® TCP NG interference screws, it is essential to establish beforehand the procedures to follow before and during the operation, and also to have a thorough knowledge of surgical procedures and to establish the type and the location of the implant to be used.
- Maximum guide pin diameter: 1.3 mm.
- Firmly mount the screw onto the screwdriver and ensure that it is fully engaged. It is essential to insert the EUROSCREW® NG screwdriver fully into the screw in order to prevent the screwdriver from slipping and damaging the internal cavity of the screw or breaking the screw during its insertion or

removal.

This is packaged and sterilized, strictly for single-use. Do not reuse, retreat or resterilize. Reusing, retreating or resterilizing the device may compromise its structural integrity and/or lead to the device malfunctioning, which can itself cause lesions, infections or death of the patient. Retreating or resterilizing single-use devices can also create a contamination risk and/or cause infections or cross-infections in the patient, including but not limited to transmitting infectious disease(s) from one patient to another. Contamination of the device can lead to lesions, infections or death of the patient.

CONTRA-INDICATIONS

- These include conditions that might weaken the firm grip of the screw thread within the bone, such as deficient bone quantity or quality due to some tumour-related or osteoporotic illness.
- Inadequate irrigation.
- Previous infection or risk of infection.
- Mental or neurological conditions that might prevent the application of the restrictions and requirements of the healing program, especially during the first two weeks following the operation, e.g., the use of drugs or some mental illness.
- Patients under 16 years old, or wherever bone growth is incomplete.
- Patients exhibiting some risk of allergic reaction to foreign bodies. In such cases prior testing is to be carried out, before the operation, in order to prevent any risk of acquired sensitivity.

SIDE EFFECTS

- Those complications ordinarily encountered with any form of internal fixation, such as fixation defects, the impossibility of regaining limb extension and/or patello-femoral problems, irritation arising from the materials used, graft conflicts and arthro-fibrosis. Complications may also take the form of a failure of the fixation and repositioning of the pin.
- As with any bio-absorbable implant, there is some risk of local inflammation during the period of the material's osseo-integration.

INFORMATION FOR THE PATIENT

The patient must be informed by the physician of the potential consequences of the factors mentioned in the 'safety precautions', 'contra-indications' and 'side-effects' paragraphs, i.e., those issues likely to hinder the success of the operation as well as any possible complications. The patient must also be informed of the measures to be taken in order to reduce the possible consequences of these factors.

PACKAGING – STORAGE

Each product is packaged in sealed double bags and sterilized using ethylene oxide. Sterility is guaranteed until the expiry date for use if the packaging has not been opened or damaged. Resterilization of the product is strictly forbidden. The screw must be stored unopened in its original packaging in a clean, dry place at a temperature less than 37°C.

RECOMMENDATIONS FOR DISPOSAL

The device or its components must be disposed of in accordance with current local waste regulations.

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