HEPAX KIT

- I** - Kit monouso sterile per biopsia epatica secondo la tecnica Menghini.
- GB** - Disposable sterile kit for liver biopsies in accordance with the Menghini technique..
- F** - Kit stérile à usage unique pour biopsie hépatique selon la technique de Menghini.
- D** - Steriles Einwegset zur Leberpunktion nach der Menghini-Technik.
- E** - Kit de un solo uso estéril para biopsia hepática según la técnica de Menghini.
- P** - Conjunto estéril não reutilizável para biópsias ao fígado de acordo com a técnica de Menghini.
- GR** - Αποστειρωμένο κιτ μιας χρήσεως για ηπατική βιοψία, σύμφωνα με την τροποποιημένη τεχνική Menghini.
- NL** - Steriele wegwerkit voor leverbiopsieën overeenkomstig de gewijzigde Menghini-techniek.
- FIN** - Kertakäyttöiset steriilit välineet maksanäytteiden ottoon.
- S** - Steril engångssats för leverbiopsi enligt den modifierade Menghini-tekniken.
- DK** - Engangs sterile sæt for lever-biopsier efter modificeret Menghini teknik..
- PL** - Sterylizowany jednorazowy sprzęt do biopsji wątroby zmodyfikowaną metodą Menghini.
- CZ** - Sterilní sada na jedno použití pro jaterní biopsii podle modifikované techniky Menghiniho.
- SK** - Sterilná súprava na jedno použitie pre pečňovú biopsiu podľa modifikovanej techniky Menghiniho.
- SLO** - Sterilni komplet za enkratno uporabo za jetrno biopsijo s prilagojeno tehniko Menghini.
- H** - Egyszer használatos, módosított Menghini-technikával végzett májbiopsziához való készlet.
- RU**

8451410000/02 del 20/07/2012

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LIMITATIONS AND WARNINGS:

- Remove the syringe from the packaging and check that the plunger is fully pushed in.
- Screw the injection needle to the luer cone on the syringe and draw up approximately 3 ml of the physiological saline thereby preventing air bubbles from being introduced into the syringe, then remove the needle.
- Tighten the biopsy needle to the syringe and remove the needle protective cap.
- Perform an incision in the skin with the scalpel (following a local anaesthetic) and penetrate the intercostal region to reach the peritoneal cavity.
- Inject approximately 2 ml of physiological saline: a certain resistance during injection indicates that the cannula of the needle is still blocked by intercostal muscles; in this case penetrate further.
- Withdraw the plunger till it stops to depress the syringe.
- Instruct the patient to hold his/her breath and quickly proceed to the liver to perform the biopsy.
- Remove the needle attached to the syringe, release the stop lever and discharge the collected sample together with the remaining physiological saline.
- At the end of the procedure, remove all devices and accessories observing the necessary precautions for potentially contaminated materials.
- Throw away after use and dispose with the hospital waste according to current legislation.

WARNINGS:

- For use only by qualified physicians. MEDAX shall not be held responsible for patient injury in case of improper or incorrect use of the device, incorrect procedures, or in case of use by non qualified personnel.
- The procedure described is for guidance only: physicians should use the device on the basis of their own experience.
- Non-pyrogenic, sterile device: not to be reused, sterilised or reconditioned in any way.
- Disposable single-use/single-patient: re-use of the disposable may cause cross infection to patient and/or user
- Check that the packaging is intact and the product is valid before use.
- Use immediately after opening the outer packaging.
- The product must not be used if the packaging is unsealed, damaged or wet.
- Not suitable for use in magnetic resonance imaging.
- Store in a cool and dry place; keep away from light and heat.

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