## OXUS- NEO OXUS

- I** - Ago per biopsia ossea.
- GB** - Bone marrow biopsy needle.
- F** - Aiguille à biopsie osseuse.
- D** - Nadel für Knochenbiopsie.
- E** - Aguja para biopsia ósea.
- P** - Agulha para biópsia à medula óssea.
- GR** - Βελόνα βιοψίας μυελού των οστών.
- NL** - Naald voor beenmergbiopsieën.
- FIN** - Neula luuydinbiopsiaan.
- S** - Nål för benmärgsbiopsi.
- DK** - Knoglemarvsbiopsinål.
- PL** - Igła biopsyjna do szpiku kostnego.
- CZ** - Bioptická jehla na kostní dřevě.
- SK** - Ihla na biopsiu kostnej drene.
- SLO** - Igla za kostno biopsijo.
- H** - Csontvelő biopszia tű.
- RU** - игла для ТРЕПАН-БИОПСИИ КОСТНОЙ ТКАНИ

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**INSTRUCTIONS FOR USE:**

- Insert the needle with the stylet perpendicular to the bone surface.
- Slowly turn and push the needle forward until the bone marrow cavity has been reached.
- Remove the stylet.
- Push the needle for 2 or 3 cm with rotating movement.
- Slowly insert the graduated stylet to verify the length of the biotpic sample. Remove the stylet.
- Withdraw the needle for a few millimeters and then insert it again with a different angle: this movement allows to obtain the bone marrow sample.
- Close the handle with the cap.
- Withdraw the needle slowly and completely with a rotating movement.
- Remove the bone marrow sample obtained pushing it out from the proximal end of the needle (handle side) inserting the remover guide from the distal end (tip side) with the help of safety cap.
- At the end, remove all devices and accessories paying attention to the necessary precautions for potentially contaminated materials.
- Throw away after the use and discharge together with the hospital waste following, the laws in force.

Note: If marrow aspiration, or explant or perfusion is to be performed, insert an appropriate syringe into the luer cone and aspirate.

**LIMITATIONS AND WARNINGS:**

- Usage allowed only to qualified physicians. MEDAX declines every responsibility for damages caused to the patient when the device is used in an improper way, with wrong manoeuvres or if handled by non qualified personnel.
- The procedure described is for guidance only: the use of the device must be based on the clinical training of each physician.
- Don't move sideways the needle into the narrow without the madrel into the cannula
- To avoid possible break of the steel cannula, avoid strong and fast movements backward and forward , with traction angles too elevated during sampling.
- Non-pyrogenic, sterile device: not to be reused, sterilised or reconditioned in any way.
- Disposable single-use/single-patient: re-use of the disposable may cause cross infection to patient and/or user
- Verify the integrity of the packaging and the validity of the product before use.
- The product must not be used if the package is unsealed, damaged or wet.
- Keep in a fresh and dry place, avoid the exposure to light and high temperature.

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