BIOFEATHER

- I** - Ago per biopsia ai tessuti molli con funzionamento semiautomatico e doppia profondità.
- GB** - Soft tissue semi-automatic biopsy needle with adjustable notch size.
- F** - Aiguille semi-automatique pour biopsie des tissus mous avec encoche ajustable.
- D** - Nadel für Biopsie für weiche Gewebe mit halbautomatischer Funktion und doppelter Tiefe.
- E** - Aguja para biopsia de tejidos blandos con funcionamiento semi-automático y doble profundidad.
- P** - Agulha para biópsia semiautomática para tecidos moles com boca de tamanho ajustável.
- GR** - Ημιαυτόματη βελόνα βιοψίας μαλακού ιστού με ρυθμιζόμενο μέγεθος εντομής.
- NL** - Biopsienaald voor zacht weefsel met halfautomatische werking en dubbele diepte.
- FIN** - Puoliautomaattinen biopsianeula pehmytkudokseen ja kaksinkertainen syvyys.
- S** - Semi-automatisk biopsinål för mjuk vävnad med justerbar skära.
- DK** - Halvautomatisk nål til blødvævsbiopsi med justerbar beholderstørrelse.
- PL** - Półautomatyczna igła biopsyjna do tkanek miękkich z regulowanym rozmiarem wycięcia.
- CZ** - Poloautomatická bioptická jehla na měkké tkáně s nastavitelnou velikostí žlábků.
- SK** - Poloautomatická ihla na biopsiu mäkkého tkaniva s nastaviteľnou veľkosťou priestoru na vzorku.
- SLO** - Igla za biopsijo mehkih tkiv s polavtomatskim delovanjem in dvojno globino.
- H** - Félautomata lágyrész biopszia tű állítható mintatartó réssel.
- RU** - Полуавтоматическая игла для биопсии мягких тканей с возможностью регулировки глубины забора биоптата.

845154000/02 del 20/07/2012

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INSTRUCTIONS FOR USE:

- Cock biopsy needle by pulling back hard on the plunger (the first click sets the notch to 15 mm, the second click sets the notch to 22 mm).
- With stylet fully retracted (specimen notch completely covered by cannula), slowly insert needle proximal to area to be sampled. Be careful not to press plunger until needle is inserted to desired location.
- Slowly move plunger with thumb to first stop to advance the stylet and expose the specimen notch. Be careful not to push plunger past first stop during this step. Move the needle back and forth a few millimeters to fill the specimen notch. Press completely the plunger with thumb, automatically triggering the cannula to close, trapping a specimen in the sample notch.
- Draw needle out. Pull back on plunger to re-cock the spring. Push plunger to first stop exposing the specimen notch. Be careful not to push plunger past first stop during this step.
- Remove tissue specimen from notch.
- At the end of procedure, remove all devices and accessories paying attention to the necessary precautions for potentially contaminated materials.
- Throw away after the use and discharge together with the hospital waste following, the laws in force.

LIMITATIONS AND WARNINGS:

- Usage allowed only to qualified physicians. MEDAX declines every responsibility for damages caused to the patient when the device is used in an improper way, with wrong manoeuvres or if handled by non qualified personnel.
- The procedure described is for guidance only: the use of the device must be based on the clinical training of each physician.
- Before the use, make sure that the mandrel slides freely and smoothly in the cannula; forward motion must open the notch and backwards motion must cover it completely without friction. Should this not be the case, do not use the needle.
- non-pyrogenic, sterile device: not to be reused, sterilised or reconditioned in any way.
- Disposable single-use/single-patient: re-use of the disposable may cause cross infection to patient and/or user.
- Verify the integrity of the packaging and the validity of the product before use.
- The product must not be used if the package is unsealed, damaged or wet.
- Keep in a fresh and dry place, avoid the exposure to light and high temperature.

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