**vi**

HAPTOGLOBIN  
9D91  
G06934R04  
B9D9BV

**ARCHITECT**

# HAPTOGLOBIN

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- Burtis CA, Ashwood ER, editors. *Tietz Textbook of Clinical Chemistry*, 2nd ed. Philadelphia, PA: WB Saunders; 1994:707–8.
- Tietz NW, editor. *Clinical Guide to Laboratory Tests*, 3rd ed. Philadelphia, PA: WB Saunders; 1995:306–9.
- Henry JB. *Clinical Diagnosis and Management by Laboratory Methods*, 18th ed. Philadelphia, PA: WB Saunders; 1991:225–7.
- US Department of Labor, Occupational Safety and Health Administration. 29 CFR Part 1910.1030. *Bloodborne Pathogens*.
- US Department of Health and Human Services. *Biosafety in Microbiological and Biomedical Laboratories*, 5th ed. Washington, DC: US Government Printing Office, December 2009.
- World Health Organization. *Laboratory Biosafety Manual*, 3rd ed. Geneva: World Health Organization, 2004.
- Clinical and Laboratory Standards Institute (CLSI). *Protection of Laboratory Workers From Occupationally Acquired Infections; Approved Guideline—Fourth Edition*. CLSI Document M29-A4. Wayne, PA: CLSI; 2014.
- Young DS. *Effects of Preanalytical Variables on Clinical Laboratory Tests*, 2nd ed. Washington, DC: AACC Press; 1997:4-497.
- US Pharmacopeial Convention, Inc. General notices. In: *US Pharmacopeia National Formulary*, 1995 ed (USP 23/NF 18). Rockville, MD: The US Pharmacopeial Convention, Inc; 1994:11.
- Ledue TB, Collins MF, Ritchie RF. Development of immunoturbidimetric assays for fourteen human serum proteins on the Hitachi 912. *Clin Chem Lab Med* 2002;40(5):520–8.
- Ritchie RF, editor. *Serum Proteins in Clinical Medicine*, Vol 1. AACC, 1996:7.04-5.
- Sasse EA, Aziz KJ, Harris EK, et al. *How to Define and Determine Reference Intervals in the Clinical Laboratory; Approved Guideline (C28-A)*. Villanova, PA: The National Committee for Clinical Laboratory Standards, 1995.
- Powers DM, Boyd JC, Glick MR, et al. *Interference Testing in Clinical Chemistry; Proposed Guideline (EP7-P)*. Villanova, PA: The National Committee for Clinical Laboratory Standards, 1986.
- Kennedy JW, Carey RN, Coolen RB, et al. *Evaluation of Precision Performance of Clinical Chemistry Devices—Second Edition; Tentative Guideline (EP5-T2)*. Villanova, PA: The National Committee for Clinical Laboratory Standards, 1992.
- Kennedy JW, Carey RN, Coolen RB, et al. *Evaluation of Precision Performance of Clinical Chemistry Devices; Approved Guideline (EP5-A)*. Wayne, PA: The National Committee for Clinical Laboratory Standards, 1999.
- Kennedy JW, Carey RN, Coolen RB, et al. *Method Comparison and Bias Estimation Using Patient Samples; Approved Guideline (EP9-A)*. Wayne, PA: The National Committee for Clinical Laboratory Standards, 1995.

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Abbott Laboratories  
Diagnostics Division  
Abbott Park, IL 60064 USA

ABBOTT  
Max-Planck-Ring 2  
65205 Wiesbaden  
Germany  
+49-6122-580

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**MANUFACTURED FOR** Abbott Laboratories

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Configure assay parameters — General			
<input checked="" type="radio"/> General	<input type="radio"/> Calibration	<input type="radio"/> SmartWash	<input type="radio"/> Results <input type="radio"/> Interpretation
Assay: <b>Hapt</b>	Type: <b>Photometric</b>	Version: †	
Number: <b>1060</b>	Run controls for onboard reagents by: <b>Lot</b>		
<input checked="" type="radio"/> Reaction definition	<input type="radio"/> Reagent / Sample	<input type="radio"/> Validity checks	
Reaction mode: <b>End up</b>	Primary	Secondary	Read times
Wavelength: <b>604</b> / <b>None</b>	Main: <b>31 – 33</b>		
Last required read: <b>33</b>			
Absorbance range: ___ - ___	Color correction: ___ - ___		
Sample blank type: <b>Self</b>	Blank: <b>14 – 16</b>		

Configure assay parameters — Calibration			
<input type="radio"/> Reaction definition	<input checked="" type="radio"/> Reagent / Sample	<input type="radio"/> Validity checks	
Reagent: <b>HAPTO</b>	Reagent volume: <b>167</b>	R1	R2
Diluent: <b>Saline</b>	Water volume: ___		
Diluent dispense mode: <b>Type 0</b>	Dispense mode: <b>Type 0</b>	<b>Type 0</b>	<b>Type 0</b>
Dilution name	Sample	Diluted sample	Dilution factor
<b>STANDARD</b>	<b>2.0</b>	___	<b>1:1.00</b>
<b>1:2</b>	<b>25.0</b>	<b>4.0</b>	<b>1:2.02</b>
<b>1:4</b>	<b>25.0</b>	<b>2.0</b>	<b>1:4.00</b>

Configure assay parameters — Calibration			
<input type="radio"/> Reaction definition	<input type="radio"/> Reagent / Sample	<input checked="" type="radio"/> Validity checks	
Reaction check: <b>None</b>			
Maximum absorbance variation: ___			

Configure assay parameters — Calibration			
<input type="radio"/> General	<input checked="" type="radio"/> Calibration	<input type="radio"/> SmartWash	<input type="radio"/> Results <input type="radio"/> Interpretation
Assay: <b>Hapt</b>	Calibration method: <b>Spline</b>		
<input checked="" type="radio"/> Calibrators	<input type="radio"/> Volumes	<input type="radio"/> Intervals	<input type="radio"/> Validity checks
Calibrator set: <b>SP</b>	Blank: <b>Water</b>	Calibrator level: <b>0†††</b>	Concentration: ___
Replicates: <b>3</b> [Range 1 – 3]	Cal 1: <b>SP1</b>	Cal 2: <b>SP2</b>	Cal 3: <b>SP3</b>
	Cal 4: <b>SP4</b>	Cal 5: <b>SP5</b>	

Configure assay parameters — Calibration			
<input type="radio"/> Calibrators	<input checked="" type="radio"/> Volumes	<input type="radio"/> Intervals	<input type="radio"/> Validity checks
Calibrator: <b>SP</b>	Calibrator level	Sample	Diluted sample
Blank: <b>Water</b>	<b>2.0</b>	___	___
Cal 1: <b>SP1</b>	<b>2.0</b>	___	___
Cal 2: <b>SP2</b>	<b>2.0</b>	___	___
Cal 3: <b>SP3</b>	<b>2.0</b>	___	___
Cal 4: <b>SP4</b>	<b>2.0</b>	___	___
Cal 5: <b>SP5</b>	<b>2.0</b>	___	___

Configure assay parameters — Calibration			
<input type="radio"/> Calibrators	<input type="radio"/> Volumes	<input checked="" type="radio"/> Intervals	<input type="radio"/> Validity checks
Calibration intervals:			
Full interval: <b>1368</b> (hours)			
Calibration type:			
Adjust type: <b>None</b>			

Configure assay parameters — Calibration			
<input type="radio"/> Calibrators	<input type="radio"/> Volumes	<input type="radio"/> Intervals	<input checked="" type="radio"/> Validity checks
Blank absorbance range: ___ - ___			
Span: <b>Blank</b> - <b>Blank</b>			
Span absorbance range: ___ - ___			
Expected cal factor: <b>0.00</b>			
Expected cal factor tolerance %: <b>0</b>			

Configure assay parameters — SmartWash				
<input type="radio"/> General	<input type="radio"/> Calibration	<input checked="" type="radio"/> SmartWash	<input type="radio"/> Results	<input type="radio"/> Interpretation
Assay: <b>Hapt</b>				
COMPONENT	REAGENT / ASSAY	WASH	Volume	Replicates
R1	DIG00	Detergent A	345	1
R1	AMIK9	Detergent A	345	1
R1	VANCO	Detergent A	345	1
R1	GENT9	Detergent A	345	1
R1	TOBRA	Detergent A	345	1
R1	DGT0B	Detergent A	345	1
R2	DIG00	Detergent A	345	1
R2	AMIK9	Detergent A	345	1
R2	VANCO	Detergent A	345	1
R2	GENT9	Detergent A	345	1
R2	TOBRA	Detergent A	345	1
R2	DGT0B	Detergent A	345	1

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Configure assay parameters — Results			
<input type="radio"/> General	<input type="radio"/> Calibration	<input type="radio"/> SmartWash	<input checked="" type="radio"/> Results <input type="radio"/> Interpretation
Assay: <b>Hapt</b>		Assay number: <b>1060</b>	
Dilution default range:		Result units: <b>mg/dL</b>	
Low-Linearity: <b>8</b>		High-Linearity: <b>250<sup>††</sup></b>	
Gender and age specific ranges:*			
GENDER	AGE (UNITS)	NORMAL	EXTREME
<b>Male</b>	<b>12 – 60 (Y)</b>	<b>14 – 258</b>	
<b>Female</b>	<b>12 – 60 (Y)</b>	<b>35 – 250</b>	
<b>Male</b>	<b>61 – 130 (Y)</b>	<b>40 – 268</b>	
<b>Female</b>	<b>61 – 130 (Y)</b>	<b>63 – 273</b>	

Configure result units	
Assay: <b>Hapt</b>	Version: †
Result units: <b>mg/dL</b>	Decimal places: <b>0</b> [Range 0 – 4]
Correlation factor: <b>1.0000</b>	Intercept: <b>0.0000</b>

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Configure assay parameters — Results			
<input type="radio"/> General	<input type="radio"/> Calibration	<input type="radio"/> SmartWash	<input checked="" type="radio"/> Results <input type="radio"/> Interpretation
Assay: <b>Hapt</b>		Assay number: <b>1060</b>	
Dilution default range:		Result units: <b>g/L</b>	
Low-Linearity: <b>0.08</b>		High-Linearity: <b>2.50<sup>††</sup></b>	
Gender and age specific ranges:*			
GENDER	AGE (UNITS)	NORMAL	EXTREME
<b>Male</b>	<b>12 – 60 (Y)</b>	<b>0.14 – 2.58</b>	
<b>Female</b>	<b>12 – 60 (Y)</b>	<b>0.35 – 2.50</b>	
<b>Male</b>	<b>61 – 130 (Y)</b>	<b>0.40 – 2.68</b>	
<b>Female</b>	<b>61 – 130 (Y)</b>	<b>0.63 – 2.73</b>	

Configure result units	
Assay: <b>Hapt</b>	Version: †
Result units: <b>g/L</b>	Decimal places: <b>2</b> [Range 0 – 4]
Correlation factor: <b>1.0000</b>	Intercept: <b>0.0000</b>

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©2009, 2015 Abbott Laboratories