



05402255001V11.0

# CREJ2

**cobas**<sup>®</sup>**Creatinine Jaffé Gen.2 (compensated)****Order information**

| REF                                    | CONTENT  | Analyzer(s) on which kit(s) can be used |
|--|--|---|
| 05401755 190                           | Creatinine Jaffé Gen.2 (compensated) (4 × 100 tests) | cobas c 111                             |
| Materials required (but not provided): |  |   |
| 10759350 190                           | Calibrator f.a.s. (12 × 3 mL)                        | Code 401                                |
| 10759350 360                           | Calibrator f.a.s. (12 × 3 mL, for USA)               | Code 401                                |
| 12149435 122                           | Precinorm U plus (10 × 3 mL)                         | Code 300                                |
| 12149435 160                           | Precinorm U plus (10 × 3 mL, for USA)                | Code 300                                |
| 12149443 122                           | Precipath U plus (10 × 3 mL)                         | Code 301                                |
| 12149443 160                           | Precipath U plus (10 × 3 mL, for USA)                | Code 301                                |
| 05117003 190                           | PreciControl ClinChem Multi 1 (20 × 5 mL)            | Code 391                                |
| 05947626 190                           | PreciControl ClinChem Multi 1 (4 × 5 mL)             | Code 391                                |
| 05947626 160                           | PreciControl ClinChem Multi 1 (4 × 5 mL, for USA)    | Code 391                                |
| 05117216 190                           | PreciControl ClinChem Multi 2 (20 × 5 mL)            | Code 392                                |
| 05947774 190                           | PreciControl ClinChem Multi 2 (4 × 5 mL)             | Code 392                                |
| 05947774 160                           | PreciControl ClinChem Multi 2 (4 × 5 mL, for USA)    | Code 392                                |
| 03121313 122                           | Precinorm PUC (4 × 3 mL)                             | Code 240                                |
| 03121291 122                           | Precipath PUC (4 × 3 mL)                             | Code 241                                |

**English****System information****CREJ2:** ACN 690**CRJ2U:** ACN 691**Intended use**

In vitro test for the quantitative determination of creatinine in human serum, plasma and urine on the **cobas c 111** system.

**Summary**<sup>1,2,3,4,5</sup>

Chronic kidney disease is a worldwide problem that carries a substantial risk for cardiovascular morbidity and death. Current guidelines define chronic kidney disease as kidney damage or glomerular filtration rate (GFR) less than 60 mL/min per 1.73 m<sup>2</sup> for three months or more, regardless of cause.

The assay of creatinine in serum or plasma is the most commonly used test to assess renal function. Creatinine is a break-down product of creatine phosphate in muscle, and is usually produced at a fairly constant rate by the body (depending on muscle mass). It is freely filtered by the glomeruli and, under normal conditions, is not re-absorbed by the tubules to any appreciable extent. A small but significant amount is also actively secreted.

Since a rise in blood creatinine is observed only with marked damage of the nephrons, it is not suited to detect early stage kidney disease. A considerably more sensitive test and better estimation of glomerular filtration rate (GFR) is given by the creatinine clearance test based on creatinine's concentration in urine and serum or plasma, and urine flow rate. For this test a precisely timed urine collection (usually 24 hours) and a blood sample are needed. However, since this test is prone to error due to the inconvenient collection of timed urine, mathematical attempts to estimate GFR based only on the creatinine concentration in serum or plasma have been made. Among the various approaches suggested, two have found wide recognition: that of Cockcroft and Gault and that based on the results of the MDRD trial. While the first equation was derived from data obtained with the conventional Jaffé method, a newer version of the second is usable for IDMS-traceable creatinine methods. Both are applicable for adults. In children, the Bedside Schwartz formula should be used.<sup>6,7,8,9</sup>

In addition to the diagnosis and treatment of renal disease, the monitoring of renal dialysis, creatinine measurements are used for the calculation of the fractional excretion of other urine analytes (e. g., albumin, α-amylase). Numerous methods were described for determining creatinine. Automated assays established in the routine laboratory include the Jaffé alkaline picrate method in various modifications, as well as enzymatic tests.

**Test principle**<sup>10,11,12</sup>

This kinetic colorimetric assay is based on the Jaffé method. In alkaline solution, creatinine forms a yellow-red complex with picrate. The rate of dye formation is proportional to the creatinine concentration in the specimen. The assay uses "rate-blanking" to minimize interference by bilirubin. To correct for non-specific reaction caused by serum/plasma pseudo-creatinine chromogens, including proteins and ketones, the results for serum or plasma are corrected by -18 μmol/L (-0.2 mg/dL).

**Reagents - working solutions**

**R1** Potassium hydroxide: 900 mmol/L; phosphate: 135 mmol/L; pH ≥ 13.5; preservative; stabilizer

**SR** Picric acid: 38 mmol/L; pH 6.5; non reactive buffer

**Precautions and warnings**

For in vitro diagnostic use for health care professionals. Exercise the normal precautions required for handling all laboratory reagents.

**Infectious or microbial waste:**

Warning: handle waste as potentially biohazardous material. Dispose of waste according to accepted laboratory instructions and procedures.

**Environmental hazards:**

Apply all relevant local disposal regulations to determine the safe disposal.

Safety data sheet available for professional user on request.

For USA: Caution: Federal law restricts this device to sale by or on the order of a physician.

This kit contains components classified as follows in accordance with the Regulation (EC) No. 1272/2008:

**Danger**

H314 Causes severe skin burns and eye damage.

EUH 001 Explosive when dry.

**Prevention:**

P280 Wear protective gloves/ protective clothing/ eye protection/ face protection.

**Response:**

P301 + P330 IF SWALLOWED: Rinse mouth. Do NOT induce vomiting.  
+ P331

P303 + P361 IF ON SKIN (or hair): Take off immediately all contaminated clothing. Rinse skin with water.  
+ P353

P304 + P340 IF INHALED: Remove victim to fresh air and keep at rest in a position comfortable for breathing.

P305 + P351 IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.  
+ P338

P310 Immediately call a POISON CENTER /doctor.

**Disposal:**

P501 Dispose of contents/container to an approved waste disposal plant.

Product safety labeling follows EU GHS guidance.

Contact phone: all countries: +49-621-7590, USA: 1-800-428-2336

**Reagent handling**

Ready for use

Reagent bottle 1 contains excess reagent to reduce the effect of CO<sub>2</sub> uptake. Any remaining reagent should be discarded after the recommended period of on-board use (see below).

**Storage and stability**

Shelf life at 15-25 °C: See expiration date on reagent

On-board in use and refrigerated on the analyzer: 4 weeks

**Specimen collection and preparation<sup>13</sup>**

For specimen collection and preparation only use suitable tubes or collection containers.

Only the specimens listed below were tested and found acceptable.

Serum

Plasma: Li-heparin and K<sub>3</sub>-EDTA plasma

The sample types listed were tested with a selection of sample collection tubes that were commercially available at the time of testing, i.e. not all available tubes of all manufacturers were tested. Sample collection systems from various manufacturers may contain differing materials which could affect the test results in some cases. When processing samples in primary tubes (sample collection systems), follow the instructions of the tube manufacturer.

Urine

Collect urine without using additives. If urine must be collected with a preservative for other analytes, only hydrochloric acid (14-47 mmol/L urine, e.g. 5 mL 10 % HCl or 5 mL 30 % HCl per liter urine) or boric acid (81 mmol/L, e.g. 5 g per liter urine) may be used.

Urine samples are automatically prediluted 1:25 (1+24) with water by the instrument.

Centrifuge samples containing precipitates before performing the assay.

See the limitations and interferences section for details about possible sample interferences.

Stability in serum/plasma:<sup>14</sup> 7 days at 15-25 °C  
7 days at 2-8 °C  
3 months at (-15)-(-25) °C

Stability in urine (without preservative):<sup>14</sup> 2 days at 15-25 °C  
6 days at 2-8 °C  
6 months at (-15)-(-25) °C

Stability in urine (with preservative): 3 days at 15-25 °C  
8 days at 2-8 °C  
3 weeks at (-15)-(-25) °C

Sample stability claims were established by experimental data by the manufacturer or based on reference literature and only for the temperatures/time frames as stated in the method sheet. It is the responsibility of the individual laboratory to use all available references and/or its own studies to determine specific stability criteria for its laboratory.

**Materials provided**

See "Reagents – working solutions" section for reagents.

**Materials required (but not provided)**

See "Order information" section

General laboratory equipment

**Assay**

For optimum performance of the assay follow the directions given in this document for the analyzer concerned. Refer to the appropriate operator's manual for analyzer-specific assay instructions.

The performance of applications not validated by Roche is not warranted and must be defined by the user.

**Application for serum, plasma and urine****cobas c 111 test definition**

|                       |                         |
|-----------------------|-------------------------|
| Measuring mode        | Absorbance              |
| Abs. calculation mode | Kinetic                 |
| Reaction direction    | Increase                |
| Wavelength A/B        | 512/583 nm              |
| Calc. first/last      | 21/26                   |
| <i>Serum/plasma</i>   |                         |
| Compensation          | -18 µmol/L (-0.2 mg/dL) |
| Unit                  | µmol/L                  |
| Reaction mode         | R1-S-SR                 |

*Urine*

|               |         |
|---------------|---------|
| Unit          | mmol/L  |
| Reaction mode | R1-S-SR |
| Predilution   | 25      |

**Pipetting parameters**

| <i>Serum/plasma</i> |        | Diluent (H <sub>2</sub> O) |
|---------------------|--------|----------------------------|
| R1                  | 13 µL  | 71 µL                      |
| Sample              | 10 µL  | 20 µL                      |
| SR                  | 17 µL  | 16 µL                      |
| Total volume        | 147 µL |                            |

*Urine*

|              |        | Diluent (H <sub>2</sub> O) |
|--------------|--------|----------------------------|
| R1           | 13 µL  | 71 µL                      |
| Sample       | 10 µL  | 20 µL                      |
| SR           | 17 µL  | 16 µL                      |
| Total volume | 147 µL |                            |

**Calibration**

|                  |   |
|------------------|---|
| Calibrator       | Calibrator f.a.s.   |
|                  | Deionized water is used automatically by the instrument as the zero calibrator. |
| Calibration mode | Linear regression   |



|                       |  |
|-----------------------|--|
| Calibration replicate | Duplicate recommended  |
| Calibration interval  | Each lot, every 7 days, and as required following quality control procedures |

Calibration interval may be extended based on acceptable verification of calibration by the laboratory.

Traceability: This method has been standardized against ID/MS.<sup>a)</sup>  
For the USA, this method has been standardized against a primary reference material (SRM<sup>b)</sup> 914 and SRM 967 (ID/MS)).

a) Isotope Dilution Mass Spectrometry

b) Standard Reference Material

### Quality control

For quality control, use control materials as listed in the "Order information" section. In addition, other suitable control material can be used.

The control intervals and limits should be adapted to each laboratory's individual requirements. Values obtained should fall within the defined limits. Each laboratory should establish corrective measures to be taken if values fall outside the defined limits.

Follow the applicable government regulations and local guidelines for quality control.

### Calculation

The **cobas c 111** analyzer automatically calculates the analyte concentration of each sample.

|                     |  |
|---------------------|--|
| Conversion factors: | $\mu\text{mol/L} \times 0.0113 = \text{mg/dL}$ |
|                     | $\text{mmol/L} \times 11.336 = \text{mg/dL}$   |

### Limitations - interference

Criterion: Recovery within  $\pm 10\%$  of initial value at a creatinine concentration of  $80 \mu\text{mol/L}$  ( $0.90 \text{ mg/dL}$ ) in serum and  $2500 \mu\text{mol/L}$  ( $28.3 \text{ mg/dL}$ ) in urine.

#### Serum/plasma

Icterus:<sup>15</sup> No significant interference up to an I index of 5 for conjugated and unconjugated bilirubin (approximate conjugated and unconjugated bilirubin concentration:  $86 \mu\text{mol/L}$  or  $5 \text{ mg/dL}$ ).

Hemolysis:<sup>15</sup> No significant interference up to an H index of 400 (approximate hemoglobin concentration:  $248 \mu\text{mol/L}$  or  $400 \text{ mg/dL}$ ).

Lipemia (Intralipid):<sup>15</sup> No significant interference up to an L index of 250. There is poor correlation between the L index (corresponds to turbidity) and triglycerides concentration.

Pyruvate: No significant interference from pyruvate up to a concentration of  $0.4 \text{ mmol/L}$  ( $3.5 \text{ mg/dL}$ ).

Ascorbic acid: No significant interference from ascorbic acid up to a concentration of  $4 \text{ mmol/L}$  ( $70 \text{ mg/dL}$ ).

Drugs: No interference was found at therapeutic concentrations using common drug panels.<sup>16,17</sup>

Exception: Cefoxitin causes artificially high creatinine results at the therapeutic drug level.<sup>18</sup> Antibiotics containing cephalosporin lead to significant false-positive values.<sup>18,19</sup>

Cyanokit (Hydroxocobalamin) may cause interference with results.

Values  $< 18 \mu\text{mol/L}$  ( $< 0.2 \text{ mg/dL}$ ) or negative results are reported in rare cases in children  $< 3$  years and elderly patients. In such cases use the Creatinine plus test to assay the sample.

Do not use Creatinine Jaffé for the testing of creatinine in hemolyzed samples from neonates, infants or adults with HbF levels  $\geq 30 \text{ mg/dL}$ .<sup>20</sup> In such cases, use the Creatinine plus test ( $\leq 600 \text{ mg/dL HbF}$ ) to assay the sample.

The presence of ketone bodies can cause artificially high results in serum and plasma.

Estimation of the Glomerular Filtration Rate (GFR) on the basis of the Schwartz Formula can lead to an overestimation.<sup>21</sup>

In very rare cases, gammopathy, in particular type IgM (Waldenström's macroglobulinemia), may cause unreliable results.<sup>22</sup>

#### Urine

Icterus: No significant interference up to a conjugated bilirubin concentration of  $854 \mu\text{mol/L}$  or  $50 \text{ mg/dL}$ .

Hemolysis: No significant interference up to a hemoglobin concentration of  $683 \mu\text{mol/L}$  or  $1100 \text{ mg/dL}$ .

Glucose: No significant interference from glucose up to a concentration of  $117 \text{ mmol/L}$  ( $2100 \text{ mg/dL}$ ).

Urobilinogen: No significant interference from urobilinogen up to a concentration of  $676 \mu\text{mol/L}$  ( $40 \text{ mg/dL}$ ).

Urea: No significant interference from urea up to a concentration of  $2100 \text{ mmol/L}$  ( $12612 \text{ mg/dL}$ ).

Drugs: No interference was found at therapeutic concentrations using common drug panels.<sup>17</sup>

High homogentisic acid concentrations in urine samples lead to false results.

Cyanokit (Hydroxocobalamin) may cause interference with results.

The presence of ketone bodies can cause artificially high results in urine.

For diagnostic purposes, the results should always be assessed in conjunction with the patient's medical history, clinical examination and other findings.

### ACTION REQUIRED

**Special Wash Programming:** The use of special wash steps is mandatory when certain test combinations are run together on the **cobas c 111** analyzer. For information about test combinations requiring special wash steps, please refer to the latest version of the carry over evasion list found with the CLEAN Method Sheet and the operator's manual for further instructions.

**Where required, special wash/carry-over evasion programming must be implemented prior to reporting results with this test.**

### Limits and ranges

#### Measuring range

##### Serum/plasma

$18\text{-}1100 \mu\text{mol/L}$  ( $0.2\text{-}12.4 \text{ mg/dL}$ )

Determine samples having higher concentrations via the rerun function. Dilution of samples via the rerun function is a 1:10 dilution. Results from samples diluted using the rerun function are automatically multiplied by a factor of 10.

##### Urine

$0.027\text{-}32.5 \text{ mmol/L}$  ( $0.31\text{-}368 \text{ mg/dL}$ )

Determine samples having higher concentrations via the rerun function. Dilution of samples via the rerun function is a 1:4 dilution. Results from samples diluted using the rerun function are automatically multiplied by a factor of 4.

### Lower limits of measurement

#### Serum/plasma

Lower detection limit of the test:

$18 \mu\text{mol/L}$  ( $0.2 \text{ mg/dL}$ )

The lower detection limit represents the lowest measurable analyte level that can be distinguished from zero. It is calculated as the value lying 3 standard deviations above that of the lowest standard (standard 1 + 3 SD, repeatability,  $n = 21$ ).

#### Urine

Lower detection limit of the test:

$0.027 \text{ mmol/L}$  ( $0.31 \text{ mg/dL}$ )

The lower detection limit represents the lowest measurable analyte level that can be distinguished from zero. It is calculated as the value lying 3 standard deviations above that of the lowest standard (standard 1 + 3 SD, repeatability,  $n = 21$ ).

### Expected values

#### Serum/plasma

##### Adults<sup>23</sup>

Females  $44\text{-}80 \mu\text{mol/L}$  ( $0.50\text{-}0.90 \text{ mg/dL}$ )

Males  $62\text{-}106 \mu\text{mol/L}$  ( $0.70\text{-}1.20 \text{ mg/dL}$ )

##### Children<sup>24</sup>

Neonates (premature)  $25\text{-}91 \mu\text{mol/L}$  ( $0.28\text{-}1.03 \text{ mg/dL}$ )

Neonates (full term)  $21\text{-}75 \mu\text{mol/L}$  ( $0.24\text{-}0.85 \text{ mg/dL}$ )



|            |              |                   |
|------------|--------------|-------------------|
| 2-12 m     | 15-37 µmol/L | (0.17-0.42 mg/dL) |
| 1- < 3 y   | 21-36 µmol/L | (0.24-0.41 mg/dL) |
| 3- < 5 y   | 27-42 µmol/L | (0.31-0.47 mg/dL) |
| 5- < 7 y   | 28-52 µmol/L | (0.32-0.59 mg/dL) |
| 7- < 9 y   | 35-53 µmol/L | (0.40-0.60 mg/dL) |
| 9- < 11 y  | 34-65 µmol/L | (0.38-0.73 mg/dL) |
| 11- < 13 y | 46-70 µmol/L | (0.52-0.79 mg/dL) |
| 13- < 15 y | 50-77 µmol/L | (0.57-0.87 mg/dL) |

**Urine****1<sup>st</sup> morning urine<sup>23</sup>**

|         |                  |                |
|---------|------------------|----------------|
| Females | 2.47-19.2 mmol/L | (28-217 mg/dL) |
| Males   | 3.45-22.9 mmol/L | (39-259 mg/dL) |

**24-h urine<sup>25</sup>**

|         |                |                     |
|---------|----------------|---------------------|
| Females | 7-14 mmol/24 h | (740-1570 mg/24 h)  |
| Males   | 9-21 mmol/24 h | (1040-2350 mg/24 h) |

Creatinine clearance for adults<sup>25,26</sup> 71-151 mL/minRefer to reference 26 for a prospective study on creatinine clearance in children.<sup>27</sup>

Roche has not evaluated reference ranges in a pediatric population.

Each laboratory should investigate the transferability of the expected values to its own patient population and if necessary determine its own reference ranges.

**Specific performance data**Representative performance data on the **cobas c 111** analyzer are given below. Results obtained in individual laboratories may differ.**Precision**

Precision was determined using human samples and controls in an internal protocol with repeatability (n = 21) and intermediate precision (3 aliquots per run, 1 run per day, 10 days). The following results were obtained:

**Serum/plasma**

| Repeatability | Mean<br>µmol/L (mg/dL) | SD<br>µmol/L (mg/dL) | CV<br>% |
|---------------|------------------------|----------------------|---------|
| Precinorm U   | 98.2 (1.11)            | 2.7 (0.03)           | 2.8     |
| Precipath U   | 353 (3.99)             | 3 (0.04)             | 0.9     |
| Human serum 1 | 66.5 (0.751)           | 2.6 (0.030)          | 4.0     |
| Human serum 2 | 548 (6.19)             | 5 (0.05)             | 0.8     |

| Intermediate precision | Mean<br>µmol/L (mg/dL) | SD<br>µmol/L (mg/dL) | CV<br>% |
|------------------------|------------------------|----------------------|---------|
| Precinorm U            | 94.8 (1.07)            | 3.5 (0.04)           | 3.7     |
| Precipath U            | 335 (3.79)             | 7 (0.08)             | 2.1     |
| Human serum 3          | 56.0 (0.633)           | 3.1 (0.035)          | 5.5     |
| Human serum 4          | 584 (6.60)             | 8 (0.09)             | 1.4     |

**Urine**

| Repeatability | Mean<br>µmol/L (mg/dL) | SD<br>µmol/L (mg/dL) | CV<br>% |
|---------------|------------------------|----------------------|---------|
| Precinorm PUC | 8.87 (101)             | 0.06 (1)             | 0.7     |
| Precipath PUC | 4.43 (50.2)            | 0.07 (0.8)           | 1.5     |
| Human urine 1 | 1.71 (19.4)            | 0.06 (0.7)           | 3.4     |
| Human urine 2 | 13.4 (152)             | 0.09 (1)             | 0.7     |

| Intermediate precision | Mean<br>µmol/L (mg/dL) | SD<br>µmol/L (mg/dL) | CV<br>% |
|------------------------|------------------------|----------------------|---------|
| Precinorm PUC          | 8.86 (100)             | 0.16 (2)             | 1.8     |
| Precipath PUC          | 4.48 (50.8)            | 0.12 (1.4)           | 2.7     |
| Human urine 3          | 1.82 (20.6)            | 0.10 (1.1)           | 5.4     |
| Human urine 4          | 13.8 (156)             | 0.4 (4)              | 2.7     |

**Method comparison**Creatinine values for human serum, plasma and urine samples obtained on the **cobas c 111** analyzer (y) were compared with those determined using the corresponding reagent on a COBAS INTEGRA 400 analyzer (x).**Serum/plasma**

Sample size (n) = 80

Passing/Bablok<sup>28</sup>

y = 0.996x + 3.276 µmol/L

r = 0.973

Linear regression

y = 0.996x + 3.680 µmol/L

r = 1.000

The sample concentrations were between 44.9 and 1057 µmol/L (0.507 and 11.9 mg/dL).

**Urine**

Sample size (n) = 50

Passing/Bablok<sup>28</sup>

y = 1.004x - 0.073 mmol/L

r = 0.977

Linear regression

y = 1.008x - 0.118 mmol/L

r = 1.000

The sample concentrations were between 1.58 and 31.3 mmol/L (17.9 and 355 mg/dL).

**References**

- 1 Thomas C, Thomas L. Labordiagnostik von Erkrankungen der Nieren und ableitenden Harnwege. In: Thomas L, ed. Labor und Diagnose, 6th ed. Frankfurt/Main: TH-Books 2005;520-585.
- 2 Lamb E, Newman DJ, Price CP. Kidney function tests In: Burtis CA, Ashwood ER, Bruns DE. Tietz textbook of clinical chemistry and molecular diagnostics. 4th ed. St.Louis, MO: Elsevier Saunders 2006;797-835.
- 3 <http://www.kidney.org/>
- 4 <http://www.nkdep.nih.gov/>
- 5 Lamb EJ, Tomson CRV, Roderick PJ. Estimating kidney function in adults using formulae. Ann Clin Biochem 2005;42:321-345.
- 6 Miller WG. Editorial on Estimating glomerular filtration rate. Clin Chem Lab Med 2009;47(9):1017-1019.
- 7 Schwartz GJ, Muñoz A, Schneider MF, et al. New Equations to Estimate GFR in Children with CKD. J Am Soc Nephrol 2009;20:629-637.
- 8 Schwartz GJ, Work DF. Measurement and Estimation of GFR in Children and Adolescents. Clin J Am Soc Nephrol 2009;4:1832-1843.
- 9 Staples A, LeBlond R, Watkins S, et al. Validation of the revised Schwartz estimating equation in a predominantly non-CKD population. Pediatr Nephrol 2010 Jul 22;25:2321-2326.
- 10 Jaffé M. Ueber den Niederschlag, welchen Pikrinsäure in normalem Harn erzeugt und über eine neue Reaktion des Kreatinins. Z Physiol Chem 1886;10:391-400.
- 11 Fabiny DL, Ertlinghausen G. Automated reaction-rate method for determination of serum creatinine with the CentrifChem Clin Chem. 1971;17:696-700.
- 12 Bartels H, Böhmer M. Micro-determination of creatinine. Clin Chim Acta 1971;32:81-85.
- 13 Guder WG, Narayanan S, Wisser H, et al. List of Analytes; Preanalytical Variables. Brochure in: Samples: From the Patient to the Laboratory. Darmstadt: GIT-Verlag 1996.
- 14 Guder WG, da Fonseca-Wollheim F, Heil W, et al. Die Qualität Diagnostischer Proben, 6. Aufl. Heidelberg: BD Diagnostics, 2009.



- 15 Glick MR, Ryder KW, Jackson SA. Graphical Comparisons of Interferences in Clinical Chemistry Instrumentation. Clin Chem 1986;32:470-475.
- 16 Breuer J. Report on the Symposium "Drug effects in Clinical Chemistry Methods". Eur J Clin Chem Clin Biochem 1996;34:385-386.
- 17 Sonntag O, Scholer A. Drug interference in clinical chemistry: recommendation of drugs and their concentrations to be used in drug interference studies. Ann Clin Biochem 2001;38:376-385.
- 18 Kroll MH. Some observations on the reaction mechanism of Cefoxitin and Cephalothin with picrate. Microchem J 1990;42:241-249.
- 19 Ducharme MP, Smythe M, Strohs G. Drug-induced alterations in serum creatinine concentrations. Ann Pharmacotherapy 1993;27:622-633.
- 20 Mazzachi BC, Phillips JW, Peake MJ. Is the Jaffe creatinine assay suitable for neonates? Clin Biochem Revs 1998;19:82.
- 21 Filler G, Priem F, Lepage N, et al.  $\beta$ -Trace Protein, Cystatin C,  $\beta$ 2-Microglobulin, and Creatinine Compared for Detecting Impaired Glomerular Filtration Rates in Children. Clin Chem 2002;48:729-736.
- 22 Bakker AJ, Mücke M. Gammopathy interference in clinical chemistry assays: mechanisms, detection and prevention. Clin Chem Lab Med 2007;45(9):1240-1243.
- 23 Mazzachi BC, Peake MJ, Ehrhardt V. Reference Range and Method Comparison Studies for Enzymatic and Jaffé Creatinine Assays in Plasma and Serum and Early Morning Urine. Clin Lab 2000;53-55.
- 24 Schlebush H, Liappis N, Klein G. Creatinine and ultrasensitive CRP: Reference Intervals from Infancy to Childhood. Clin Chem Lab Med 2001;39 Special Supplement PO-T042;1-448.
- 25 Junge W, Wilke B, Halabi A, et al. Determination of reference intervals for serum creatinine, creatinine excretion and creatinine clearance with an enzymatic and a modified Jaffé method. Clin Chim Acta 2004;344:137-148.
- 26 Zawta B, Delanghe J, Taes Y, et al. Arithmetic Compensation for Pseudo-Creatinine Interferences of the Creatinine Jaffé Method and its Effect on Creatinine Clearance Results. Clin Chem Part 2, Suppl S June 2001;46(6):487.
- 27 Wuyts B, Bernard D, van den Noortgate N, et al. Reevaluation of Formulas for Predicting Creatinine Clearance in Adults and Children Using Compensated Creatinine Methods. Clin Chem 2003;49:1011-1014.
- 28 Bablok W, Passing H, Bender R, et al. A general regression procedure for method transformation. Application of linear regression procedures for method comparison studies in clinical chemistry, Part III. J Clin Chem Clin Biochem 1988 Nov;26(11):783-790.

A point (period/stop) is always used in this Method Sheet as the decimal separator to mark the border between the integral and the fractional parts of a decimal numeral. Separators for thousands are not used.

Any serious incident that has occurred in relation to the device shall be reported to the manufacturer and the competent authority of the Member State in which the user and/or the patient is established.

### Symbols

Roche Diagnostics uses the following symbols and signs in addition to those listed in the ISO 15223-1 standard (for USA: see dialog.roche.com for definition of symbols used):

|   |                                       |
|---|---------------------------------------|
|  | Contents of kit                       |
|  | Reagent                               |
|  | Volume after reconstitution or mixing |
|  | Global Trade Item Number              |

### FOR US CUSTOMERS ONLY: LIMITED WARRANTY

Roche Diagnostics warrants that this product will meet the specifications stated in the labeling when used in accordance with such labeling and will be free from defects in material and workmanship until the expiration date printed on the label. THIS LIMITED WARRANTY IS IN LIEU OF ANY OTHER WARRANTY, EXPRESS OR IMPLIED, INCLUDING ANY IMPLIED WARRANTY OF MERCHANTABILITY OR FITNESS FOR PARTICULAR PURPOSE. IN NO EVENT SHALL ROCHE DIAGNOSTICS BE LIABLE FOR INCIDENTAL, INDIRECT, SPECIAL OR CONSEQUENTIAL DAMAGES.

COBAS, COBAS C, COBAS INTEGRA, PRECINORM, PRECIPATH and PRECICONTROL are trademarks of Roche.

All other product names and trademarks are the property of their respective owners.

Additions, deletions or changes are indicated by a change bar in the margin.

© 2020, Roche Diagnostics



Roche Diagnostics GmbH, Sandhofer Strasse 116, D-68305 Mannheim  
www.roche.com

+800 5505 6606



Distribution in USA by:  
Roche Diagnostics, Indianapolis, IN  
US Customer Technical Support 1-800-428-2336

