



Document Information:

Document Author: Jaeschke, Irmgard {DOMC~Mannheim}, Maschat, Cordula {DOMC~Mannheim}
Business Area / Unit: Roche Professional Diagnostics
Confidentiality: Confidential
Document Class: Device Master Record
Document Type: Device Master Record
Document Creator: Jaeschke, Irmgard {DOMC~Mannheim}
Document Lifecycle Status: Signed
Valid From: 16-Nov-2018 00:07:22 (UTC)
Valid To:
Document Title: C_0201_RSP_05401755000
Document Number: 0000000000001008018000925
Document Version: 03
Template: No

Global Group: Regulatory Affairs
Global SubGroup: Specification Prod. (Registr.)
Local Group: Reg. Specification
Language: English
Site: RDG Germany
Department: C_DOMCC
Document Applies To: C_QC Centralized Assays
Document Description: Creatinine Jaffe Gen.2 (CREJ2) - "Only for use in Regulatory Affairs"
Change Request eCCP PLM: 400030307

Electronic Signatures:

Signed By: maschac1 (Cordula Maschat {DOMCSD..6164})
Role: Reviewer
Signature Differentiation: Quality Control
Signed Date: 15-Nov-2018 08:48:00 (UTC)

Signed By: loewh (Heiner Loew {DOMC....6164})
Role: Approver
Signature Differentiation: Quality Control
Signed Date: 15-Nov-2018 09:55:56 (UTC)



Specification

Product: **Creatinine Jaffé Gen.2
(CREJ2)**

Catalog No. 05401755190

Intended use:

In vitro test for the quantitative determination of creatinine in human serum, plasma and urine on the cobas c 111 system.

Test principle:

This kinetic colorimetric assay is based on the Jaffé method. In alkaline solution, creatinine forms a yellow-orange complex with picrate.

The rate of dye formation is proportional to the creatinine concentration in the specimen.

The assay uses "rate-blanking" to minimize interference by bilirubin.

To correct for non-specific reaction caused by serum/plasma pseudo-creatinine chromogens, including proteins and ketones, the results for serum or plasma are corrected by -18 µmol/L (-0.2 mg/dL).

Creatinine + picric acid $\xrightarrow{\text{AlkalinepH}}$ yellow-orange complex

Roche Diagnostics GmbH
Roche Diagnostics Global Operations
Quality Control Mannheim

D-68305 Mannheim
Phone: +49-621-759-0
Fax: +49-621-759-2890

1. Reagents-working solutions:

1.1 Alkaline buffer solution (R1):

1.1.1 pH-value (25 °C): ≥ 13.5

1.2 Picric acid solution (R2):

1.2.1 pH-value (25 °C): 6.00 - 7.00

2. Performance test:

2.1 Recovery on cobas c 111 analyzer:

PreciControl ClinChem Multi 1 (Cat. No. 05947626190): 90 - 110 %

PreciControl ClinChem Multi 2 (Cat. No. 05947774190): 90 - 110 %

Also required:

Calibrator f.a.s. (Cat. No. 10759350190)