OneGuide Metal Sleeve User Manual

(2020.10 V1.0)

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(2020.10V1.0)

1. Description

- OneGuide Metal Sleeve is a product that is connected to OneGuide Template (for Metal Sleeve) to prevent the template from being split from the rotation of the drill during drilling, and protects the fixture surface from the inside of the template hole when placing the fixture.
- The OneGuide Metal Sleeve is a product for implant treatment and is optimized for the Osstem Implant system. Do not use with other similar products.

2. Purpose for Use

A tool used to guide the path, location, and surgical site indication of an implant or device.

3. Usage and Procedures

- 1) Correct dimensions for each hole size of OneGuide Template printed by 3D printer.
- 2) Connect the Metal Sleeve to the Press-in Jig and insert it by aligning the Hex direction of the OneGuide Template.
- 3) Using a syringe, inject the adhesive (3D printing material) into both adhesive holes and cure it sufficiently.
- 4) When assembly is complete, proceed with the procedure in the same way as general drilling. (Not available for Sleeve alone)

4. Precautions

- For the safe and effective use of Osstem Implant, dental implant procedures require a highly specialized and complex procedure, so it is highly recommended to perform the procedure after completing a professional training course.
- As this product is a medical device, it is recommended to use this product for surgery after soaking it in 70% ethanol solution.

[Use of Autoclave and hydrogen peroxide is prohibited.]

5. Precautions during Use

- It must be used by a medical professional, and use for other purposes is prohibited.

- If there is a problem with the exterior or if the product is dropped on the floor and received a severe shock, it must be carefully observed and replaced with another product if a problem occurs.
- To prevent damage to OneGuide Template, be careful not to apply excessive force to the drilling process.

6. How to Store after Use

This product is a disposable medical device. So discard after use and it is prohibited to reuse.

7. Side Effects

After placing an implant, problems such as loss of implant stability, damage of prosthesis, etc. may occur. These problems may be caused by a range of factors, including qualitative or quantitative deficiencies of the remaining bones, infection, allergic reactions, poor oral hygiene or non-coordination of the patient, movement of the implants, local tissue degeneration, inappropriate location and arrangement of implants, etc.

8. Warning

This product should be treated only with the directions prescribed by Osstem Implant, and it is forbidden to modify and use in any way.

9. Contraindicated Symptom

Insufficient osseointegration or implant fail may occur to the patients with the symptoms of insufficient bone, poor bony tissue, excessive smoking, blood disorder, or uncontrolled diabetes.

10. Product Information

The Package label includes the product name, specifications, manufactured date, product Lot no., etc.

※ For more information, please see the catalog and procedure manual provided by your sales representative.