

Lấy mẫu

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Blood Sampler for Blood Gas

Only for blood collection – DO NOT USE FOR INJECTION!

Product description:

Primed disposable syringe with cap

Fill volume: printed on product packaging

Anticoagulant: Quantity in IU: printed on product packaging

For blood collection for pH/blood gas, oxymetry, electrolyte and metabolite analysis

Material: Plastic

Syringe Tip: LUER

Storage conditions: Room temperature / store in a dry place

Intended use

For in vitro diagnostic use for blood collection with the purpose of analysis of pH/blood gas, oxymetry, electrolyte and metabolite on a blood gas/ISE instrument.

Notes and limitations

- The syringe may only be used by trained medical personnel.
- Observe the required hygiene laws and regulations when handling biological samples. Dangerous causative agents may be present. Protective gloves and protective clothing must be worn. Danger of infection!
- After blood collection attach cap and mix sample thoroughly! (see point A – Sampling or pictograms on reverse)
- The analysis should be performed immediately after taking the sample. If the sample is stored too long in the syringe before performing the measurement, it may lead to sedimentation and changes in the measured values. This occurs particularly in samples with high hematocrit, extreme measured values, increased self-metabolism as well as increased body temperature! (see point B – Storage)
- Maximum storage period of the sample at room temperature is 15 minutes. The sample may lose its characteristics when stored beyond 15 minutes, even if kept in ice-water slush. (see IFCC study under Literature)
- Mix sample thoroughly before measurement! (see point C – sample preparation or pictograms on reverse)

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A) Sampling

1. Open sterile packing
2. Pull plunger once with the cap still in place
3. Sampling to be done by trained medical personnel only! Remove the syringe cap. Sampling by:
 - a) puncture
 - b) 3-way connector
 - c) in-dwelling catheterPlease check the fill volume required by the sampler. Other mixtures can influence the measurement!
The decision for the sampling procedure and type of sample has to be taken by a physician.
4. Put the plunger cap back on
5. Roll syringe at least five times
6. Turn syringe at least five times

B) Storage / Important

Maximum storage period at room temperature is 15 minutes.

The usability of the sample will be restricted if it is stored for longer than 15 minutes, even if kept in ice-water slush (see IFCC study under Literature).

C) Sample preparation and measurement

1. Roll syringe at least five times
2. Turn syringe at least five times
3. Remove the plunger cap, discard the first sample drop to avoid clogging the analyzer with a blood clot
4. If necessary, use clot catcher
5. Please aspirate the sample as recommended by the instrument manufacturer (see Users Manual)

Disposal

The syringe is for single use only. It has to be disposed of according to the legal and laboratory regulations.

It is biologically contaminated and has to be treated as hazardous material.

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