CÔNG TY TNHH ASIA ACTUAL VIETNAM
MST: 0316765196

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
F&P 850 System

F&P Optiflow+

MicroCell Technology™
Heated Inspiratory Adult
Breathing Circuit and
Optiflow™+ Nasal Cannula

USER INSTRUCTIONS

REF RT332

14








CE 0123

Rx only

www.fphcare.com

**Fisher & Paykel
HEALTHCARE**




Manufacturer Fisher & Paykel Healthcare Ltd, 15 Maurice Paykel Place, East Tamaki, Auckland 2013, PO Box 14 348 Panmure, Auckland 1741, New Zealand Tel: +64 9 574 0100 Fax: +64 9 574 0158 Email: info@fphcare.com Web: www.fphcare.com Importer/ Distributor **Australia** (Sponsor) Fisher & Paykel Healthcare Pty Ltd, 19-31 King Street, Nunawading, Melbourne, Victoria 3131. Tel: +61 3 9871 4900 Fax: +61 3 9871 4998 **Austria** Tel: 0800 29 31 23 Fax: 0800 29 31 22 **Benelux** Tel: +31 40 216 3555 Fax: +31 40 216 3554 **Brazil** Fisher & Paykel do Brasil, Rua Sampaio Viana, 277 cj 21, Paraisópolis, 04004-000, São Paulo - SP, Brazil Tel: +55 11 2548 7002 **China** 代理人/售后服务机构: 费雪派克医疗保健 (广州) 有限公司, 广州高新技术产业开发区科学城科丰路31号G12栋301号 电话: +86 20 32053486 传真: +86 20 32052132 **Denmark** Tel: +45 70 26 37 70 Fax: +46 83 66 310 **Finland** Tel: +358 94 1590 355 Fax: +46 83 66 310 **France** **EC REP** Fisher & Paykel Healthcare SAS, 10 Av. du Québec, Bât F5, BP 512, Villebon-sur-Yvette, 91946 Courtaboeuf Cedex, France Tel: +33 1 6446 5201 Fax: +33 1 6446 5221 Email: c.s@fphcare.fr **Germany** Fisher & Paykel Healthcare GmbH & Co. KG, Deutschland, Österreich, Schweiz, Wiesenstrasse 49, D 73614 Schorndorf, Germany Tel: +49 7181 98599 0 Fax: +49 7181 98599 66 **Hong Kong** Tel: +852 2116 0032 Fax: +852 2116 0085 **India** Tel: +91 80 2309 6400 **Ireland** Tel: 1800 409 011 Fax: +44 1628 626 146 **Italy** Tel: +39 06 7839 2939 Fax: +39 06 7814 7709 **Japan** Tel: +81 3 5117 7110 Fax: +81 3 5117 7115 **Korea** Tel: +82 2 6205 6900 Fax: +82 2 6309 6901 **Norway** Tel: +47 21 60 13 53 Fax: +47 22 99 60 10 **Russia** Tel. and Fax: +7 495 782 21 50 **Spain** Tel: +34 902 013 346 Fax: +34 902 013 379 **Sweden** Tel: +46 8 564 76 680 Fax: +46 8 36 63 10 **Switzerland** Tel: 0800 83 47 63 Fax: 0800 83 47 54 **Taiwan** Tel: +886 2 8751 1739 Fax: +886 2 8751 5625 **Turkey** İthalatçı Firma: Fisher Paykel Sağlık Ürünleri Ticaret Limited Şirketi, İletişim Bilgileri: Ostim Mahallesi 1249. Cadde No:6, Yenimahalle, Ankara, Türkiye 06374, Tel: +90 312 354 34 12 Fax: +90 312 354 31 01 **UK** Fisher & Paykel Healthcare Ltd, Unit 16, Cordwallis Park, Clivernock Road, Maidenhead, Berkshire SL6 7BU, UK Tel: 0800 132 189 Fax: +44 1628 626 146 **USA/Canada** Tel: 1800 446 3908 or +1 949 453 4000 Fax: +1 949 453 4001

For patent information, see www.fphcare.com/ip
F&P 850 System, F&P Optiflow and MicroCell Technology are trademarks of Fisher & Paykel Healthcare Ltd.

REF 900532 REV A 2019-04 © 2019 Fisher & Paykel Healthcare Limited

English (en)[ⓘ]

Symbols

		
This product is intended to be used for a maximum of 14 days	Storage temperature range	Single Use Only
		
Use by	Incorrect water level, replace MR290 Chamber	Correct water level in the MR290 chamber
		
Authorised representative in the European community	Consult instructions for use.	Type BF applied part.
		
Conforms with medical device directive 93/42/EEC	Date of Manufacture	Manufacturer
		
Catalogue number	Batch code	Prescription only

Intended Use

The Fisher & Paykel Healthcare RT332 is a single use consumable kit intended for use with the MR850 to deliver heated and humidified nasal high flow therapy to spontaneously breathing adult patients who require breathing support. The product is designed for use in the hospital environment and must be prescribed by a physician.

Technical Specifications

Inspiratory Limb Length	1.5m
Compressible volume	1L
Flow Range	<60 L/min
Connection to Flow Source	ISO 5356-1 Conical Connector
Minimum Internal Diameter	14mm
Flow Resistance @ 30 L/min BTPS (including 0.02 cmH ₂ O measurement uncertainty)	0.85 ± 0.07 cmH ₂ O
Maximum tube surface temperature	44°C
Maximum chamber operating pressure	8 kPa
This product was not made with natural rubber latex	

Setup

MR850 Humidifier in invasive mode. RT332 kit with 22mm heated inspiratory tube and chamber, ventilator in high flow mode.

- Select appropriate size. Prongs must not create a seal in the nares. A clear gap must be visible around each prong.
- Adjust headstrap to fit. Do not over-tighten.
- Ensure headstrap clip is attached, to prevent cannula from being pulled out of nares.
- Cannula can become unattached if not used with the headstrap clip.
- Attach tubing clip to clothing/bedding to prevent cannula from pulling off face.
- Attach the tubing clip to breathing circuit but ensure probe cable is not crushed by the tubing clip.

Attention

- Check all connections are tight before use.
- Check that the heater wire is evenly distributed along the circuit and not bunched or kinked.
- Avoid prolonged contact with patient’s skin.
- Set appropriate ventilator alarms.
- Check breathing circuit and interface for condensation every 6 hours. Drain as required to prevent loss of therapy.
- Before connecting the interface, check for adequate gas flow and ensure that the system has warmed up.
- Use USP sterile water for inhalation or equivalent

Warning

Failure to comply with the following warnings may impair performance of the device or compromise safety (including potentially cause serious harm):




- Visually inspect components and accessories for damage before use and replace if damaged. Use of damaged components or accessories may impair performance of this device or compromise safety.
- The use of breathing circuit/chamber/interface combinations not recommended by Fisher & Paykel Healthcare may result in poor humidification performance, ventilator malfunction and harm to the patient/user.
- For single patient use only. Reuse may result in transmission of infectious substances, interruption to treatment, or serious harm. Attempting to reprocess will result in degradation of materials and render the product defective.
- Nasal delivery of respiratory gases generates flow-dependent positive airway pressure (PAP). This must be taken into account where PAP could have adverse effects on a patient.
- Appropriate patient monitoring must be used at all times. Failure to monitor the patient may result in loss of therapy, serious injury or death.
- When mounting a humidifier adjacent to a patient ensure that the humidifier is always positioned lower than the patient.
- Ensure there is a water supply connected to the chamber and that water is present within the chamber.
- Do not use if packaging is not sealed.
- Do not use heated wire breathing circuits without gas flow. If gas flow is interrupted turn the humidifier off.
- Do not use the chamber if the seals are not intact when received, or if it has been dropped.
- Do not spike the water source until the blue caps have been removed.
- Do not reinsert the blue cap into the chamber after removal. This could damage components, impair the performance of this device or compromise safety.
- Do not use the chamber if the water level rises above the maximum water level line.
- Do not operate the chamber at an angle in excess of 10°.
- Do not fill the chamber with water in excess of 37 °C.
- Do not touch the heater plate or chamber base. Surfaces may exceed 85 °C.
- Do not use near a naked flame, to avoid fires.
- Do not crush or stretch tube, to prevent loss of therapy (cannula).
- Do not stretch or milk the tubing (breathing circuit).
- Do not cover circuit with materials such as blankets, towels or bed linen as this can affect the quality of the therapy or injure the patient.
- Do not soak, wash, sterilize, or reuse this product. Avoid contact with chemicals, cleaning agents, or hand sanitizers.
- California residents please be advised of the following, pursuant to Proposition 65: This product contains chemicals known to the State of California to cause cancer, birth defects and other reproductive harm. For more information, please visit: www.fphcare.com/prop65.

Notes

- Discard according to appropriate hospital protocol. User may be exposed to breathing tract fluids during disposal.
- This breathing set is validated for use with MR850-system humidifiers only.
- The responsible organisation is accountable for the compatibility of the humidifier and all of the parts and accessories used to connect to the patient and other equipment before use.

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