Vision
■ ZIRKON

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WOHLWEND AG

Einführung

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Press over Zircon

It is possible to press over all commercially Zirconia systems (green, sintered or hiped). The grinded or sintered Zirconia framework has to be prepared with special grinders or rubberwheels to manufacturer's specifications. Here, the manufacturer's minimum strengths of copings and thickness of connections have to be followed.

First step: Liner

1. Apply the special universal liner in a thin and watery consistence to the framework and fire as shown in the firing chart. This serves to reduce the brightness value of the zircon framework material, and to get a better bonding.

Attention: Different firing temperature.

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Second step: Wax up


Use 100% ash free waxes only! Put the framework without pressure on the isolated model. It is possible to wax up in full contour, or to reach a higher esthetic work in a reduced way (cut back).

In such cases apply from our Vision Zircon incisals and transpa powders. There is no loss of dimension during pressing.

Third step: Spruing

Use 6 or 8 gauge ash free sprues with a length of 4 to max. 7 mm. Sprues are connected to the coping in a straight line so as not to restrict the flow of the molten ingot. Sprues should be angled slightly outward with a 3 mm distance between copings.

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Wichtig:

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Fourth step: Investing

Wax weight up to 0.6 g = 1 Ingot
Wax weight 0.6–1.3 g = 2 Ingots

Optional:




Paper ring former or silicon ring former.

Remove strip from paper, form a ring to the appropriate mark. Slip onto base and adapt stabilizing ring. Pour special investment to just below the stabilizing ring.

Carefully remove the stabilizing ring and fit the ring former, twisting it carefully. Investment must leak from the opening. After bench set, remove paper / Silicon ring, base and ring former, then clean the bottom surface with a plaster knife. This must be absolutely flat, so that the ring stands perfectly upright in the pressing furnace.

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Sixth step: Devasting

Mark the length using spare plunger. Cut along the line, using separating disk. Attention, use mask and dust suction. Break segments apart.

Carefully devaste the objects by sandblasting. Use glass beads only.

Start with 60 psi pressure, as soon as objects appear turn pressure down to 35 psi proceed blasting carefully. DO NOT USE ALUMINIUM OXIDE!

Carefully separate the sprue from the objects with a sharp diamond disk.

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Neunter Schritt: Glaser

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Seventh step: Layering (cut back)

After working out the press channels, layer up with Vision Zircon. Complete the incisal part and fire according to the chart. With dentin you can layer supplemented portions. If necessary make a correction firing. Before glazing you can customize using Vision Zircon stains. Fire on ceramic pins. It is important to ensure that the ceramic shoulder has no contact with the pin. It can be glazed or polished.

Eighth step: Staining technique

Individual nuances can now be added with Vision Zircon Stains.

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16 Dentine powder A1–D4 20 g
2 Modellingliquid 25 ml

ASSORTMENT S2

Incisal: 57, 58, 59, 60, 20 g
Opalincisal: Neutral, Blue 20 g; Intesiveincisal: White, Yellow 20 g
Transpa Clear 20 g, Glaze 20 g, Liner: 1 to 3 and Pressliner 20 g,
1 Modellingliquid 25 ml, 1 Glaze & Stainliquid 25 ml

ASSORTMENT S3

Transpa: Neutral, Yellow, Blue, Grey 20 g
Opaquedentine: White, Honey 20 g,
Gingival: Pink and Softpink 20 g
Shouldermasses: Creme, Yellow, Redbrown 20 g
1 Modellingliquid 25ml, 1 Shoulderliquid 25 ml

ASSORTMENT S4

Pressables each 5 pc: A1, A2, A3, B1, B2, B3, C1, C2, D2, Bleach,
TA1, TA2, TA3, TB1, TB2, TB3, TC1, TC2, TD2

STAIN-ASSORTMENT

White, Honey, Olive, Redbrown, Chestnut, Grey, Blue 3 g,
Glaze, A, B, C, D, 3 g, Glaze & Stainliquid 25 ml

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Accessory

- Porcelain croco
- Brushes
- Mixing trays
- Ceramic isolation pen
- Firing paste
- Plunger investment
- Plunger Aluminium oxide
- Silicon ring former
- Investment

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**Wohlwend AG