Vision
■ CLASSIC

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TECHNICAL MANUAL

WOHLWEND AG

Einführung

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Introduction

Since its introduction at the 1995 IDS in Cologne, Germany, Vision Classic has established itself as an easy to use ceramic with consistent shade accuracy. As a clinically proven ceramic, Vision Classic is now used in over 40 countries. Vision Classic is a special fine-structured ceramic developed for all alloys in a CTE range of 13,8–15,2 (25–600°C). Vision Classic also has a natural light transmission and reflexion. By using fluorescent and opalescent additional powders very individual and esthetic results are achieved.

Compatible with alloys having a CTE of 13.8 – 15.2 $\mu\text{m}/\text{mK}$ (25–600°C)

- Precious alloys
- Semi-precious
- Non-precious
- Non greening in combination with alloys containing silver

Wohlwend's unique manufacturing process gives an ideal grain size and pigment distribution. Benefits of this unique process include:

- Consistently accurate shade reproduction – V shades A1–D4
- Easy handling
- No slumping during build up
- Less shrinkage during firing

The vital shade reproduction and natural fluorescence of Vision porcelains correspond to natural teeth. This replication of nature can be obtained by either using a basic two layer buildup or by incorporating additional dentine and incisal modifiers into the restoration.

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Alloy Parameters

- Compatible with alloys having a CTE of 13.8–15.2 (25–600°C).
- Ideally at 14.5 (25–600°C).
- If the CTE is higher 14.7 (25–600°C):
Use slow cooling cycle!
- If below 14.2 (25–600°C):
Remove restoration from lift as soon as furnace has fully opened.
- These parameters may vary between different alloys.
- In case of uncertainty we recommend to use Vision Bonding.

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Material Properties

We confirm that all porcelain powders stated on this list have been tested and conform to EN ISO 9693:2000 and 6872:1998.


Produktname Product name		Vision Classic Massen (Metallkeramik) Typ 1 Porcelain for metal typ 1			
Massen Powders		Biegefestigkeit Flexural Strength	Löslichkeit Solubility	mittlerer WAK thermal expansion	Glastransformationspkt. Glass transition temperature
		MPa	µg/cm ³	25–500°C	TG°C ± 10
Pulveropaquer Powderopaque	Op	> 100 MPa	50 µg/cm ³	13,8	600°C
Paste m. Aerosil Pasteopaque	PO	> 100 MPa	50 µg/cm ³	13,8	630°C
Dentin + Schneide Dentin + Incisal	D+S	> 80 MPa	25 µg/cm ³	13,3	590°C
Opak Dentin Opaque dentin	OD	> 80 MPa	25 µg/cm ³	13,3	590°C
Int. Opaquer Intensive opaque	IO	> 100 MPa	50 µg/cm ³	13,8	600°C
Transpa Transparent	T	> 75 MPa	25 µg/cm ³	13,3	590°C
Schultermasse Shoulder	SP	> 80 MPa	25 µg/cm ³	13,8	590°C
Effekt Schneiden Effect Incisal	EF	> 80 MPa	25 µg/cm ³	13,3	590°C
Int. Dentin Intensive dentin	ID	> 80 MPa	25 µg/cm ³	13,3	590°C
Gingival Gingiva	ZF	> 80 MPa	25 µg/cm ³	13,3	590°C

Farbe Shade	A1	A2	A3	A3.5	A4	B1	B2	B3	B4	C1	C2	C3	C4	D2	D3	D4
Opaker Opaque	O	O	O	O	O	O	O	O	O	O	O	O	O	O	O	O
Opakdentin Opacious Body	OD	OD	OD	OD	OD	OD	OD	OD	OD	OD	OD	OD	OD	OD	OD	OD
Dentin Body	D	D	D	D	D	D	D	D	D	D	D	D	D	D	D	D
Schneide Incisal	S 58	S 58	S 59	S 59	S 60	S 57	S 59	S 59	S 59	S 60	S 59	S 59	S 60	S 60	S 59	S 59
Transpamassen Transpa materials	TR: SCL – CL – NT – OP – RT – BL – GR															

Erweitert Farbe Modified Shade	A1	A2	A3	A3.5	A4	B1	B2	B3	B4	C1	C2	C3	C4	D2	D3	D4
Schultermassen shoulder material	SP neutral + BA					SP neutral + BB				SP neutral + BC				SP neutral + BD		
Effekt Schulter Effect shoulder	ivory, straw, yellow-gold, brown, yellow-grey, brown-grey															
Opaker zur Chromasteigerung Opaque to intensify chroma	O BA				O BB				O BC				O BD			
Intensivopaker zum charakterisieren Intensive opaques for characterization	10%	20%	30%	35%	40%	10%	20%	30%	40%	10%	20%	30%	40%	20%	30%	40%
	I O W GB OR OK BR V															
Dentin zur Chromasteigerung Body to intensify chroma	D BA				D BB				D BC				D BD			
Intensivdentin zum charakterisieren Intensive dentine for characterization	10%	20%	30%	35%	40%	10%	20%	30%	40%	10%	20%	30%	40%	20%	30%	40%
	I D EB ST HG BR MG BL NEON															
Intensivschneide zum individualisieren Intensive incisal for individualization	I S 7 8 9 10															
Effektmassen Marmelons Effect masses for marmelons	ivory, straw, honeyyellow, khaki															

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


Metal processing prior to porcelain application



1. Process alloys according to manufacturers instructions. Avoid sharp edges.

Unless specified otherwise, the alloy should be sandblasted with min. 110 µm aluminum oxide at 29 psi.

A uniform oxide layer is desired. In case of NP, the oxides are to be removed using min. 110 µm aluminum oxide. (Refer to manufacturers manual!). A clean and grease-free surface must be obtained.

2. Apply Vision bonding paste in an uniform and thin layer. Firing cycle according to the bonder manual.
3. The result after firing is also replacing a washfiring.

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First step: 1st Opaque

1. Powder / Liquid:
Mix the opaque to a creamy texture. Apply using brush or glass instrument with a 70% coverage.
Firing temperature 1778°F.
2. After first firing.

Vis Opaque - Paste

Apply using brush or glass instrument. Firing temperature 1778°F. Can be reduced to degas temperature of alloys. Prolong hold time at final temp. 2–3 min. Secure well after usage to avoid a dry out of the material.

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Second step: 2nd Opaque

1. Apply the opaque with a full covering layer.
2. Egg shell appearance after firing.


Options

To achieve internal effects, characterizations can be made by using opaque-modifiers "IO".

As an example:
Layer IO-orange in the central Fossa or IO-brown in the cervical area. To reach a higher chroma from the depth, apply some Chroma opaque to the equatorial area or admix as per shadechart.

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Schulterkombination Basic Vision Classic:

Third step: Shoulder

Our shoulder powders in basicshades A, B, C, D and Neutral are giving you the chance to create a shoulder for every shade situation. Please follow the mixing ratio in the chart below.

Optionally you can reach a larger color palette by using effectshoulder masses.

Shoulder combination Basic Vision Classic:

Proportionales Mengenverhältnis proportional factor					
Basic	2 Teile portions	1 Teil portion	3 Teile portions	4–5 Teile portions	5 Teile portions
Neutral	3 Teile portions	1 Teil portion	2 Teile portions	2 Teile portions	1 Teil portion
Farbe Shade					
A	A1	A2	A3	A3.5	A4
B	B1	B2	B3		B4
C	C1	C2	C3		C4
D		D2	D3		D4
Kann durch individuelle Farbschemata verändert werden. Can be changed by using an individual color scheme.					

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



Erster Dentinbrand.




Fourth step: Build up

1. In case of less space apply a thin layer of opacious dentin (occlusal, cervical, etc.).
2. Build up with dentin (shade that is given) in full contour.
3. Do a small cut back from incisal and approximal. Apply the borders and mamelons with a thin layer of transpa neutral.
4. Complete the anatomic form of the tooth with enamel according to the shade. A light oversize of the build up, brings a perfect result, because of the small shrinkage during the firing cycle.

Use some opacious dentine to fill up interdental space.

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Fifth step: Correction

1. After the first dentine firing.
2. If needed correct the form with Transpa Super Clear.
3. Proceed with second dentine firing.

Sixth step: Finishing

4. When the contact points and the occlusion are worked out, you have got two possibilities to glaze.

Glaze firing without glaze powder and subsequent a mechanical polish (special rubber wheels, brush, Dia glaze and so on). This is only recommended if the situation is well known.

Glaze firing with glaze powder. Mix it with glaze liquid and apply a thin layer. For individual characterization use Vision Classic stains.

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Multi layering technique

The use of Vision's extended range of modifiers in conjunction with simple layering techniques allow a natural looking restorations. These can be used straight or mixed into other materials.

Shade example: Orange neck, Body between A3 and A3.5, Incisal A2.



To raise the chroma mix 30% ID BA with 70% dentine A3.



Neck area Intensive body ID OR, Upper third opacious body OD A3.

Build up with dentin A3. Reduce dentin slightly more than usual.

Incisal S58 and transpa TR Opal are alternatly layered to achieve the shade but also getting natural opaleszenz.

TR Opal can be mixed with all transpas, depending on which level of translucency is to be achieved in correction.

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Multi layering technique

Mix the highly fluorescent ID Neon with the intensive incal IS 7, to achieve higher value and apply in a thin layer.

For mammelons our effect masses will be used very carefully and in a thin layer.
Effectmasses (EF): Ivory, Straw, Honeyyellow, Khaki

Neon can also be mixed with appropriate stain and be applied onto wet surface.
(What you see is what you get!)

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Classic Firing chart

Recommended firing parameters!

Temperatures can vary due to different furnace performances.

* For the Vis-Opak pasteopaque extend the dry time for 4 minutes.

Follow alloy manufacturers instructions. Horizontal furnace systems (Programat®) may require reducing the start temperatures to 50°C / 122°F.

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16 Dentine powders A1–D4 20 g
2 Modellingliquid 25 ml

ASSORTMENT S3

Incisal: 57, 58, 59, 60 20 g; Intensiveincisal: 7, 8, 9, 10 20 g
Transpa: Clear, Superclear, Neutral, Opal, Red, Blue, Grey 20 g
Glaze 10 g; 1 Modellingliquid 25 ml; 1 Glazeliq 25 ml

ASSORTMENT S4

16 Opaque-Dentine powders A1–D4 20 g
2 Modellingliquid 25 ml

ASSORTMENT S5

Intensivopaque: White, Yellow, Orange, Ocher, Brown, Violet 20 g
Intensive dentin: Ivory, Straw, Honeyyellow, Brown, Grey,
Blue, Neon 20 g; Gingiva: Softpink, Darkpink 20 g
1 Opaqueliquid 25 ml; 1 Modellingliquid 25 ml

ASSORTMENT S6

Shouldermasses: Ivory, Straw, Yellowgold, Brown, Yellowgrey,
Brown 20 g; Effectmasses: Ivory, Straw, Honeyyellow, Khaki 20 g
1 Shoulderliquid 25 ml; 1 Modellingliquid 25 ml

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STAINASSORTMENT

01, 02, 03, 04, 05, 06, 07, 08, 09, 10, 11, 12, 13, 14, 15, Glaze
 Additional grey and violet; 1 Brush; Glaze & Stainliquid 25 ml

VIS-OPAK-ASSORTMENT

16 Paste opaque A1–D4 4 g
 White, Ocher, Lila, Pink 3 g; 1 Brush

Zubehör

- Porzellan Krokodil
- Pinsel
- Anmischplatte
- Keramik Isolierstift
- Brennpaste
- Bonder

Accessory

- Porcelain croco
- Brushes
- Mixing trays
- Ceramic isolation pen
- Firing paste
- Bonding

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**Wohlwend AG