## CORMAY MULTICALIBRATOR LEVEL 1 CORMAY MULTICALIBRATOR LEVEL 2

(PL)

### CORMAY MULTICALIBRATOR LEVEL 1

Nr kat. 5-174 10 x 5 ml

Nr kat. 5-176 4 x 5 ml

### CORMAY MULTICALIBRATOR LEVEL 2

Nr kat. 5-175 10 x 5 ml

Nr kat. 5-177 4 x 5 ml

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Cat. No 5-177 4 x 5 ml

#### INTENDED USE

CORMAY MULTICALIBRATOR LEVEL 1 / LEVEL 2 is designed for use as a calibrator in clinical chemistry assays.

The reagents must be used only for *in vitro* diagnostic, by suitably qualified laboratory personnel, only for the intended purpose, under appropriate laboratory conditions.

#### CHARACTERISTICS

CORMAY MULTICALIBRATOR LEVEL 1 / LEVEL 2 is based on lyophilised human serum. The concentrations of organic and inorganic components and enzymes activities are suitable for calibration of clinical chemistry assays on a wide range of automatic analysers. Constituent concentrations are available at 2 levels.

#### INSTRUCTION FOR USE

- Open the vial carefully and add **exactly 5ml** of distilled water at 20-25°C.
- Replace the rubber stopper and leave to stand for **30 minutes** (in case of the ALP determinations stand for 1 hour) out of bright light.
- Swirl **gently** several times during reconstitution period to ensure that the contents are completely dissolved.
- Prior to use, mix the contents by inverting the vial several times. **Do not shake! Avoid the foam formation!**

#### STABILITY AND STORAGE

- Lyophilised serum remains stable when stored at 2-8°C until expiry date given on the product label.
- The constituents in reconstituted control serum are stable for **(except of parameters mentioned in notes)**:  
Temp. 2-8°C 7 days  
Temp. -20°C 30 days

#### NOTES

- The **activity of alkaline phosphatase** in reconstituted calibrator will rise over the stability period. It is recommended that the reconstituted or thawed calibrator be allowed to stand for 1 hour at 20-25°C before calibration. ALP is stable for 6 days at 2-8°C.
- Bilirubin** as well as other organic compounds is light sensitive and it is recommended that the serum be stored in the dark. Stored in the dark bilirubin is stable for 2 days at 2-8°C. Do not store at room temperature and do not freeze the serum assigned for bilirubin level control.
- MG** is stable for 2 days at 2-8°C and 14 days at -20°C.
- Glucose, TG, AST, CK** are stable for 4 days at 2-8°C.
- dtIBC** is stable for 4 days at 2-8°C and 7 days at -20°C.
- In case of microbial contamination or excessive turbidity of the reconstituted calibrator do not use for assays calibration.

- Only one freezing of the reconstituted calibrator is recommended! For freezing prepare portions of the calibrator no less than 500 µl in small (1.5 ml), tightly-closed vials. Thaw at room temperature for about 30 minutes (in case of the ALP determinations stand for 1 hour). Swirl gently several times during unfreezing.
- This calibrator has been tested for HBsAg and the HIV and HCV antibody, and found to be non-reactive. However this material should be handled as thought capable of transmitting infectious disease.

#### WASTE MANAGEMENT

Please refer to local legal requirements.

#### LITERATURE

- Wallhofer H., Schmidt E., Schmidt U.F. W.: Synopsis Der Leberkrankheiten. G. Thieme Verlag, Stuttgart (1974).
- Thefeld W., et al: Dtsch. Med. Wschr. 99, 343 (1974).
- Bergmeyer H.U., Horder M., Rej R.: J. Clin. Chem. Clin. Biochem. 24, 481-495 (1986).
- Henry R.J. Cannon D.C. Winkerman J. W.: Clinical Chemistry Principles and Technics, 2nd ed. Hagerstown MD: Harper and Row, 815, 888 (1974).
- Burtis C.A., Ashwood E.R., ed. Tietz Textbook of Clinical Chemistry, 3<sup>rd</sup> ed. Philadelphia, PA: Moss D. W., Henderson A. R., 652 (1999).
- Dembńska-Kieć A., Naskalski J.W.: Diagnostyka laboratoryjna z elementami biochemii klinicznej, Volumes, 776, (1998).

Date of issue: 05.2019