

HƯỚNG DẪN SỬ DỤNG TIẾNG ANH

Tài liệu được xác nhận bằng chữ ký số

Hà Nội, ngày 07 tháng 05 năm 2022

Người đại diện hợp pháp của cơ sở

GIÁM ĐỐC


Nguyễn Thị Minh Phương



DIAPATH

CytoPath® - Gyn Filters

Filters for thin layer preparation of gyn specimens



Manufacturer: Diapath S.p.A. via Savoldini 71 – 24057 Martinengo (BG) – Italy
Tel. +390363-986411 Fax +390363-948000 www.diapath.com info@diapath.com



Code	Description	Packaging
CP102	CytoPath – GYN filter (Pink)	100 pcs

Intended use

Disposable Filter made of inert material dedicated for the preparation of liquid based citological samples.

Composition

Pink colored cylindric filters made of acrylonitril-styrene copolymer
PET Membrane

Characteristics

Filter's dimension (height)	5.5 cm
Membrane dimension	2 cm
Membrane's pore dimension	7-8 µm

Warning: Do not touch the membrane with hands. Do not leave the filters package open in order to protect them from dust exposure.

Instruction for use

Gyn Filter are intended for the disposable with automatic or semiautomatic devices for thin layer preparation of liquid based cytological samples.

The Cytological sample, collected and preserved into a dedicated vial filled with fixative solution, is filtered with the above mentioned equipments in order to collect the cells inside it.

To avoid mistakes, qualified and trained staff should use the product. Professional use product. The guidelines concerning safety on the workplace must be applied according to current regulations. The tools used for diagnosis must be suitable for diagnostic use in laboratory. The diagnosis should be performed only by authorized, trained and competent staff.

Quality Control

The products and the raw materials are entered and constantly monitored by computer systems that allow traceability between batch number of each single product and batches of their raw materials.

Storage

Store the product according to the specifications listed on the label. The product, if opportunely stored and integrally packed, is stable up to the expiry date reported on the label. Do not use after expiration date.

If the product is not stored as recommended, its performance may change and must be validated by the user. After the opening the reagent is stable till expiration date only if stored inside its container and according to all recommendations shown on the label. Tightly close the container after the use.





DIA PATH

Disposal instructions

The expired and/or unused product must be disposed according to local waste regulations, based on danger classification on the label and after possible contaminations evaluation. In some cases, it may be necessary an analytical evaluation to determine the correct waste classification and the danger features.

Labeling legend



Production



Manufacturer



Storage temperature



Product code



Expiry date



In vitro diagnostic medical device



Photosensitive

For more information, see the MSDS

