TERUMOBCT

11220

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Bruksanvisning

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Intended Use

The Spectra Optia® Single-Needle Connector is intended for use with the Spectra Optia Exchange Set.

Please refer to the Spectra Optia Apheresis System Operator's Manual for a complete listing of contraindications, adverse effects, warnings, cautions and operating instructions.

Warnings

1. Do not disconnect the patient access from the single-needle connector during a single-needle procedure until rinseback is complete.
2. **Do Not Reuse/Not for Reuse:** Terumo BCT, Inc. products bearing the "Do Not Reuse" symbol are intended for single use only and are not intended to be reused or re-sterilized in any manner. Terumo BCT cannot ensure the functionality or sterility of the product if it is reused or re-sterilized.

Reuse of a single-use product could result in:

- Product performance issues due to a loss of product integrity, including but not limited to the following:
 - Fluid leaks
 - Parts that are warped or deformed
 - Plastics that are brittle and discolored
 - Filters that have reduced filtration capabilities
- Exposure to excessive ethylene oxide (EO) residuals
- Viral infections such as hepatitis or human immunodeficiency virus (HIV)
- Bacterial infections
- Cross-contamination

Any of these risks could result in serious injury or death. These risks are shared by product users, donors, patients, and recipients of end products of the device.

Cautions

1. Terumo BCT has validated the single-needle procedure using the Spectra Optia Single-Needle Connector. Terumo BCT is not responsible for the performance of other connectors.
2. Do not use the single-needle connector if the end caps are not in place.

Storage Conditions

Long-term storage range: 0 °C to 35 °C (32 °F to 95 °F)

Storage humidity: 0% to 75% relative humidity (RH), non-condensing

Permitted excursions:

- -29 °C to 0 °C (-20 °F to 32 °F) for up to 72 hours
- 35 °C to 50 °C (95 °F to 122 °F) for up to 6 weeks
- Up to 85% RH ± 5% for up to 72 hours

Procedures for Use

Refer to the Spectra Optia Apheresis System Operator's Manual.

Symbols and Certifications

Symbol	Definition
	Indicates the product quantity when the quantity is placed in the square.
	Indicates the date after which the medical device is not to be used.
	Indicates a non-pyrogenic fluid pathway.
	Indicates the presence of a sterile fluid path. The method of sterilization is ethylene oxide (EO).
	Indicates that the user should consult the instructions for use.
	Indicates the date when the medical device was manufactured (or sterilization date, if the product is sterile).
	Indicates the manufacturer's catalog number so that the medical device can be identified.
	Indicates a medical device that is intended for one use, or for use on a single patient during a single procedure.
	Indicates that the use of product is by prescription only.
	Indicates that the product was manufactured in accordance with Annex II of the European Council Directive, 93/42/EEC, as amended.
	Indicates the manufacturer's batch code so that the batch or lot can be identified.
	Indicates the medical device manufacturer.
	Indicates a medical device that needs protection from light sources.
	Indicates a medical device that needs to be protected from moisture.

Symbol	Definition
	Indicates that the product packaging complies with European Directive 94/62/EC for packaging and packaging waste.
	Indicates a medical device that is not to be resterilized.
	Indicates the range of temperatures to which the device can be safely exposed.

Return of Used Product

If for any reason this product must be returned to Terumo BCT, Inc., a returned goods authorization (an RGA number) is required from Terumo BCT prior to shipping.

Instructions for cleaning and materials, including appropriate shipping containers, proper labeling, and an RGA number, may be obtained from the Terumo BCT Quality Assurance Department.

IT IS THE RESPONSIBILITY OF THE HEALTH CARE INSTITUTION TO ADEQUATELY PREPARE AND IDENTIFY THE PRODUCT FOR RETURN SHIPMENT.

Please contact your local representative for information regarding returned goods and product complaints.

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©2019 Terumo BCT, Inc.