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Method Comparison

Regression Statistics^h

Comparative Method	Slope	Intercept mg/dL [g/L]	Correlation Coefficient	n
Beckman Array 360	0.95	0.6 [0.006]	0.991	127 ⁱ

h. Model equation for Passing-Bablok regression statistics is: [Result of Dimension® system] = [Slope x comparable method result] + Intercept.

i. The range of transferrin values in the correlation study was 89 – 426 mg/dL [0.89 – 4.26 g/L].

Specificity

HIL Interference

The TRNF method was evaluated for interference from hemolysis, icterus and lipemia according to NCCLS EP7-P. Bias, defined as the difference between the control sample (does not contain interferent) and the test sample (contains the interferent), is shown in the table below. Bias exceeding 10% is considered "interference".

Substance Tested	Test Concentration SI Units	TRNF Concentration mg/dL [g/L]	Bias ^j %
Hemoglobin (hemolysate)	1000 mg/dL	191 [1.91]	<10
	[0.62 mmol/L] (monomer)		
Bilirubin (unconjugated)	80 mg/dL	187 [1.87]	<10
	[1368 µmol/L]		
Lipemia (Intralipid®)	600 mg/dL	198 [1.98]	<10
	[6.78 mmol/L]		

j. Analyte results should not be corrected based on this bias.

Non-Interfering Substances

The following substances at the concentrations indicated have no significant effect (less than 10%) on the TRNF method when added to a serum pool containing 257 mg/dL [2.57 g/L] transferrin.

Substance	Test Concentration	SI Units
Acetaminophen	20 mg/dL	1323 µmol/L
Albumin	3.6 g/dL	36 g/L
Amikacin	15 mg/dL	256 µmol/L
Ascorbic Acid	3 mg/dL	170 µmol/L
Caffeine	10 mg/dL	515 µmol/L
Carbamazepine	12 mg/dL	508 µmol/L
Chloramphenicol	25 mg/dL	774 µmol/L
Chlordiazepoxide	2 mg/dL	67 µmol/L
Chlorpromazine	5 mg/dL	157 µmol/L
Cholesterol Ester	500 mg/dL	12.9 mmol/L
Cimetidine	10 mg/dL	396 µmol/L
Creatinine	30 mg/dL	2652 µmol/L
Dextran 75	2500 mg/dL	333 µmol/L
Diazepam	2 mg/dL	70 µmol/L
Digoxin	5 ng/mL	6.4 nmol/L
Erythromycin	20 mg/dL	273 µmol/L
Ethanol	350 mg/dL	76 mmol/L
Ethosuximide	30 mg/dL	2125 µmol/L
Ferrous Sulfate	6.3 mg/dL	225 µmol/L
Furosemide	2 mg/dL	61 µmol/L
Gentamicin	12 mg/dL	251 µmol/L
Heparin (Sodium)	8 U/mL	8000 U/L
Hydralazine	1 µg/mL	6.2 µmol/L
Ibuprofen	40 mg/dL	1939 µmol/L
Iron Dextran	1.4 mg/dL	190 µmol/L
Lidocaine	6 mg/dL	256 µmol/L
Lithium Chloride	3.5 mg/dL	5.07 mmol/L
Nicotine	2 mg/dL	123 µmol/L
Penicillin G	25 U/mL	25000 U/L
Pentobarbital	10 mg/dL	442 µmol/L
Phenobarbital	15 mg/dL	646 µmol/L
Phenytoin	10 mg/dL	396 µmol/L
Primidone	10 mg/dL	458 µmol/L
Procainamide	10 µg/mL	36.8 µmol/L
Propoxyphene	0.4 mg/dL	12 µmol/L
Rheumatoid Factors	1040 IU/mL	1040 IU/mL
Salicylic Acid	50 mg/dL	3.62 mmol/L
Theophylline	25 mg/dL	1388 µmol/L
Urea	500 mg/dL	83.3 mmol/L
Uric Acid	20 mg/dL	1.2 mmol/L
Valproic Acid	50 mg/dL	3467 µmol/L

Recovery

Recovery of reference material CRM470/RPPHS⁸ ranged from 92 – 93%, with a mean recovery of 92%.

Analytical Sensitivity: ≤40 mg/dL [0.40 g/L]

The analytical sensitivity represents the low assay range of TRNF.

Symbols Key: See adjacent panel.

Dimension® and Flex® are trademarks of Siemens Healthcare Diagnostics.

©2008 Siemens Healthcare Diagnostics
All rights reserved.

Dimension® clinical chemistry system

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©2008 Siemens Healthcare Diagnostics
Reservados todos los derechos.

Bibliography/Literatur/Bibliographie/Bibliografia/Bibliografía:

1. Silverman LM, et al. In: Tietz NW, ed. Textbook of Clinical Chemistry. Philadelphia: WB Saunders Company, 1986: 519.
2. Jeppson J.-O., Transferrin, Serum Proteins in Clinical Medicine Volume I Laboratory Section, First Edition, Foundation for Blood Research, Scarborough, ME, 1996.
3. Thomas L. Clinical Laboratory Diagnostics, T-H Books Verlagsgesellschaft, Frankfurt, 1998, First English Edition, pp 275-277 (Transferrin Saturation).
4. Clinical and Laboratory Standards Institute/NCCLS. *Procedures for the Collection of Diagnostic Blood Specimens by Venipuncture; Approved Standard – Fifth Edition*. CLSI/NCCLS document H3-A5 [ISBN 1-56238-515-1]. CLSI, 940 West Valley Road, Suite 1400, Wayne, PA 19087-1898 USA, 2003.
5. Clinical and Laboratory Standards Institute/NCCLS. *Tubes and Additives for Venous Blood Specimen Collection; Approved Standard – Fifth Edition*. CLSI/NCCLS document H1-A5 [ISBN 1-56238-519-4]. CLSI, 940 West Valley Road, Suite 1400, Wayne, PA 19087-1898 USA, 2003.
6. Clinical and Laboratory Standards Institute/NCCLS. *Procedures for the Handling and Processing of Blood Specimens; Approved Guideline – Third Edition*. CLSI/NCCLS document H18-A3 [ISBN 1-56238-555-0]. CLSI, 940 West Valley Road, Suite 1400, Wayne, PA 19087-1898 USA, 2004.
7. National Committee for Clinical Laboratory Standards. How to Define, Determine, and Utilize Reference Intervals in the Clinical Laboratory. CLSI/NCCLS Document C 28-A, 1992, Wayne, PA
8. Whicher JT, et al. New International Reference Preparation for Proteins in Human Serum (RPPHS). Clin. Chem. 1994; 40: 934.

Symbols Key
Symbolschlüssel
Explication des Symboles
Interpretazione simboli
Clave de los Símbolos

	Do not reuse / Nicht zur Wiederverwendung / Ne pas réutiliser / Non riutilizzare / No reutilizar
	Use By / Verwendbar bis / Utiliser jusque / Utilizzare entro / Fecha de caducidad
	Batch Code / Chargenbezeichnung / Code du lot / Codice del lotto / Código de lote
	Catalogue Number / Bestellnummer / Référence du catalogue / Numero di catalogo / Número de catálogo
	Caution, consult accompanying documents / Achtung, Begleitdokumente beachten / Attention voir notice d'instructions / Attenzione, vedere le istruzioni per l'uso / Atención, ver instrucciones de uso
	Manufacturer / Hersteller / Fabricant / Fabbicante / Fabricante
	Authorized Representative in the European Community / Bevollmächtigter in der Europäischen Gemeinschaft / Mandataire dans la Communauté européenne / Mandatario nella Comunità Europea / Representante autorizado en la Comunidad Europea
	Contains sufficient for <n> tests / Inhalt ausreichend für <n> Tests / Contenu suffisant pour "n" tests / Contenido suficiente per "n" saggi / Contenido suficiente para <n> ensayos
	In Vitro Diagnostic Medical Device / In-Vitro-Diagnostikum / Dispositif médical de diagnostic in vitro / Dispositivo medico-diagnostico in vitro / Producto sanitario para diagnóstico in vitro
	Temperature Limitation / Temperaturbegrenzung / Limites de température / Limiti di temperatura / Limite de temperatura
	Consult Instructions for Use / Gebrauchsanweisung beachten / Consulter les instructions d'utilisation / Consultare le istruzioni per l'uso / Consulte las instrucciones de uso
	Non-sterile / Nicht steril / Non stérile / Non sterile / No estéril
	CE Mark / CE Zeichen / Marquage CE / Marchio CE / Marca CE
	Contents / Inhalt / Contenu / Contenuto / Contenido
	Reconstitution Volume / Rekonstitutionsvolumen / Volume de reconstitution / Volume di ricostituzione / Volumen de reconstitución
	Level / Konzentration / Niveau / Livello / Nivel

2014-10_EFR6S

Siemens Healthcare Diagnostics Inc.
 500 GBC Drive
 Newark, DE 19714 USA

Global Siemens
 Headquarters
 Siemens AG
 Wittelsbacherplatz 2
 80333 Muenchen
 Germany

Global Siemens
 Healthcare Headquarters
 Siemens AG
 Healthcare Sector
 Henkestrasse 127
 91052 Erlangen
 Germany
 Phone: +49 9131 84-0
 siemens.com/healthcare

Global Division
 Siemens Healthcare
 Diagnostics Inc.
 511 Benedict Avenue
 Tarrytown, NY 10591
 USA
 siemens.com/healthcare