HbA1c 2R CALIBRATOR SET

Lot: 0000000

Cat. No.	Pack Name	Packaging (Content)
XSYS0097	HbA1c 2R CAL SET	5 × 0.5 ml

EN

INTENDED USE

HbA1c 2R CALIBRATOR SET is a set of five lyophilized calibrators with different levels based on human blood material (erythrocytes).

STORAGE AND STABILITY

The calibrators both unopened and opened must be stored at 2–8 °C, protected from light and heat.

Unopened calibrators are stable until the date of expiry stated at the label.

After opening and reconstitution are calibrators stable for one week if stored under 2–8 °C and if contamination and evaporation are avoided after having opened the vials.

Proper storage and handling must be observed.

WARNINGS AND PRECAUTIONS

Reagents of the kit are not classified like dangerous.

Each individual blood donation used for production of HbA1c 2R CAL SET was found to be non-reactive when tested with approved methods for HBs-Ag, HIV-Ag•Ab and HCV-Ab HTLV-1-Ab. As there is no possibility to exclude definitely that products derived from human blood transmit infectious agents, it is recommended to handle the calibrator with the same precautions used for patient specimens.

PREPARATION

Before use perform reconstitution of all five calibrators by adding 500 µl of Hemolyzing solution (R3) into each vial, mix gently and allow to stand for 5 minutes. After reconstitution are the calibrators ready to use. Don't prepare hemolysate.

PROCEDURE

Please refer to the reagent package insert for instructions for use.

ASSAY VALUES

The calibrator values have been made traceable to the approved IFCC reference method. Values according to DCCT/NGSP in % have been derived from the values according to IFCC by calculation. Calibrator values listed below are specific for this lot number of calibrator only.

WASTE MANAGEMENT

Please refer to local legal requirements.

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LITERATURA

1. The Diabetes Control and Complications Trial Research Group. The effect of intensive treatment of diabetes in the development and progression of long-term complications in insulin-dependent diabetes mellitus. N Engl J Med. 1993; 329:977-86.

2. Little RR, Rohlfing CI, Wiedmeyer HM, Myers GL et al. The national Glycohemoglobin Standardization Program: A five-Year Progress Report. Clin Chem 2001; 47:1985-92.

3. Jeppsson JO, Kobold U, Barr J, Finke A et al. Approved IFCC Reference method for the measurement of HbA1c in human blood. Clin Chem Lab Med 2002; 40:78- 89.

4. Hoelzel W, Weykamp C et al. IFCC Reference system for measurement of HbA1c in Human Blood and the national standardization Schemes in the United States, Japan and Sweden: A Method Comparison Study. Clin Chem 2004; 50:1:1666-74.

5. Biosafety in Microbiological and Biomedical Laboratories. U.S. Department of Health and Human Services, Washington 1993 (HHS Publication No. [CDC] 93-8395).

ATESTOVANÉ HODNOTY

HbA1c hladina	Šarže	Datum expirace	Hodnota kalibrátoru	
			Podle DCCT/NGSP (%)	Podle IFCC (mmol/mol)
1	xxxxxxx	mm/rrrr	0	0
2	xxxxxxx	mm/rrrr	xx,x	xx,x
3	xxxxxxx	mm/rrrr	xx,x	xx,x
4	xxxxxxx	mm/rrrr	xx,x	xx,x
5	xxxxxxx	mm/rrrr	xx,x	xx,x