

MAL CONTROL

Lot: 0000000

Cat. No..	Pack name	Packaging (Content)
BLT20033	MAL CON	1 x 1 ml



INTENDED USE

Accuracy control for quantitative immunochemical determination of Microalbumin in urine by turbidimetry and nephelometry.

COMPOSITION

A dilution of defibrinated human plasma with phosphate buffered saline, liquid stabilised and filtered through 0.2 µm.

Preservative: 0.095 g % sodium azide.

Ready for use.

STORAGE AND STABILITY

The expiry date of the product at 2-8 °C is listed on the label.

After first opening the container, the control can be used for 6 weeks if stored tightly closed at 2-8 °C after use.

Do not freeze.

PRECAUTIONS AND WARNING

1. For *in Vitro* Diagnostic use.
2. Each individual donation intended for use in manufacture of protein control serum was tested for hepatitis B surface antigen (HBsAg), anti-hepatitis C virus (anti-HCV) and anti-HIV1 and HIV 2 by FDA required test. Only donations with negative findings were used for its manufacture. Nevertheless every product obtained from human body fluids should be handled with appropriate care in accordance with recommended procedures for biohazardous materials since absence of infectious agents can never be proven.
3. Reagents containing sodium azide must be handled with due caution: Do not ingest or allow to contact skin or mucous membranes! Sodium azide can form explosive azides when contacting heavy metals such as copper or lead.

ASSIGNED VALUES

The value for CRP has been assigned to ERM-DA474/IFCC

Microalbumin: xyz mg/l

USED SYMBOLS



Catalogue Number



Manufacturer



See Instruction for Use



Lot Number



In vitro Diagnostics



Storage Temperature



Expiry Date



Content

QUALITY SYSTEM CERTIFIED
ISO 9001 ISO 13485



Erba Lachema s.r.o., Karásek 2219/1d, 621 00 Brno, CZ
e-mail: diagnostics@erbamannheim.com, www.erbamannheim.com