



HƯỚNG DẪN SỬ DỤNG
DÂY TRUYỀN MÁU NIPRO BLOOD TRANSFUSION SET

HƯỚNG DẪN SỬ DỤNG

1. Đóng khóa kẹp chỉnh giọt
2. Mở nắp đậy đầu nhựa và cắm hết chiều dài của nó vào lỗ ra của túi máu/ chai máu
3. Treo túi máu, bóp bầu đếm giọt có màng lọc và thả ra nhẹ nhàng để lấy máu vào bầu, lấy khoảng nửa bầu
4. Mở nắp đậy kim (kim đặt vào tĩnh mạch)
5. Mở khóa chỉnh giọt từ từ và mỗi dây truyền
6. Đặt kim vào tĩnh mạch
7. Mở khóa kẹp và chỉnh giọt

CHÚ Ý:

1. Dùng lại hoặc tái sử dụng thiết bị dùng một lần có thể gây nhiễm khuẩn và chức năng của thiết bị hoặc tính năng thống nhất của thiết bị có thể tổn hại
2. Không dùng nếu bao bì bị hỏng
3. Dùng xong bỏ vào thùng chứa nguy hại sinh học



250mm

100mm

IS2101TH-G_NST



BLOOD TRANSFUSION SET

1mL=20drops (1drop = 0.050mL) FILTER 200 μm
Sterilized using ethylene oxide Do not re-use

NEEDLE
18G x 1 1/2"

250mm

100mm

CONTRAINDICATION
Do not reuse.

INSTRUCTIONS FOR USE

1. Close control clamp.
2. Remove protector from plastic needle and insert full length into outlet port of container.
3. Suspend container and gently squeeze and release filter chamber to fill with blood, then drip chamber until half full by same manner.
4. Remove protector from vein needle.
5. Open control clamp gradually to allow blood flow to vein needle.
6. Perform venipuncture.
7. Regulate flow rate with control clamp.

CAUTIONS

- Reuse or reprocessing of a single use device may lead to contamination and compromised device function or structural integrity.
- Do not use if package is damaged.
- Dispose of the device in an approved biohazard container as per facility protocol.



Amsino International, Inc.
708 Corporate Center Drive,
Pomona CA 91768, USA
<http://www.amsino.com>

Emergo Europe B.V.
Prinsessegracht 20
2514 AP The Hague
THE NETHERLANDS

Distributed by:
NIPRO SALES (THAILAND) CO., LTD
89/170, Nipro Building, Viphavadee Rangsit Road,
Talat Bangkhen, Laksi, Bangkok 10210, THAILAND
Made in China

LOT
LOT

Date of manufacture

Use-by date

Color: Pantone 171

AMSINO

	NAME	DATE	TITLE
DRAWN			IS2101TH-G_NST pouch label
CHKD.			
APPR.			LABEL#: L4685 REV. C