Separating Abrasive Strips

INSTRUCTIONS FOR USE

Contraindications:

- Separating abrasive strips contain nickel and should not be used for individuals with known allergic sensitivity to this metal as it may cause hypersensitivity.
- This product contains nickel, a chemical known to the state of California to cause cancer, birth defects or other reproductive harm.

Material Used:

Stainless steel- 304 grade.

Indications:

Beveling, finishing and polishing of interproximal areas.

Instruction for Use:

- The total length of the strip (150mm) allows for optimal finger support on the adjacent teeth during work.
- To get past the contact area, exert a slight horizontal movement. Do not twist, bend or pull strip too tightly during use.
- After introduction, work with the less abrasive strip (Ref# 160-220/160-223), in order to facilitate the access of the more abrasive strip (Ref#106-221/106-224). Then begin the finishing procedure using a back and-forth movement.
- Moving ahead utilize the coarse separating strip with cautious horizontal movements to model the proximal contours. After contouring, switch again to the medium and fine separating strip to complete the procedure.

Maintenance and Sterilization of all mechanical instruments:

- Instruments are delivered non-sterile. Sterilization must be carried out according to the validated procedure before use.
- Sterilization process is referred to ISO 17644. Sterilization is conducted by steam autoclave with the following parameters:
 - Fractioned pre-vacuum (3 times)
 - Temperature: 134°C
 - Pressure: 1 bar
 - Holding time: 2.5 min. (Half cycle)
 - Drying time: 10 min.

Note:

Our instruments are not suitable for a chemiclave or hot air sterilization.

Warning:

- Avoid bending and levering: use of forces greater than 3N could cause damage/breakage of the instrument.

Disposal of Discarded Product and Packaging:

Ensure safe and proper disposal of used/discarded product and packaging to avoid adverse effect to environment. The disposal should be in accordance to local laws in the country of use.

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Ortho Organizers, Inc.
mdi Europa GmbH