Pipelle® (1111000)
Pipelle® Mark II (1111200)

Composition

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Pipelle® (1111000)
Pipelle® Mark II (1111200)

Composition

Polyethylene cannula with a closed distal end and one distal side eye, fitted with a plunger.
The surface of the Pipelle is completely smooth. The surface of the Pipelle Mark II is roughened at its distal end.

Indication

Endometrial biopsy for histological analysis.
The Pipelle Mark II also allows cells to be sampled from the surface of the endometrium in order to combine a cytological analysis with the histological analysis.

Contraindications

This product is contraindicated if the patient is pregnant.

Warnings

Contents are sterile as long as original packaging is intact. Do not use if the packaging is damaged.
Do not use after the expiry date.
Class I disposable sterile device.
Do not re-use: single-use only. Re-use can result in cross-infection independent of the method of cleaning or re-sterilization used.
Biological hazard: after use, this product presents a potential biological risk. It should be handled and disposed of in compliance with current local and national guidelines.
Users must take precautions to protect themselves from the risk of AEs (accidental exposure to blood or body fluids).
Only specialist practitioners should use this device.

Instructions for use

Pipelle and Pipelle Mark II

1. Place the patient in the dorsal decubitus position with the feet in the stirrups.
2. Insert a vaginal speculum to expose the cervix.
3. The Pipelle de Cornier can be bent prior to removing it from the sterile packaging. Some elastic resistance allows the device to maintain an anterior curve or reverse curvature.
4. Disinfect the cervix according to current institutional practices.
5. Gently guide the Pipelle while threading it through the cervix to the uterine fundus. Graduated mark 4 indicates the start of the uterine cavity. Graduated mark 7 is usually reached, except if there is a cervical obstruction.
6. Create a vacuum by pulling the plunger out to the end.
7. Gently sweep several times from the back of the uterus to the internal os of the cervix using regular and repeated back-and-forth movements at the same time as rotations in order to sample the entire uterine cavity.
8. Wait until you can see fragments of uterine mucosa, which have gradually become detached, passing through the cannula. The procedure generally takes 30 seconds.
9. Withdraw the device completely.
10. To collect the cytology sample from the Pipelle Mark II, rinse the roughened end in a bottle containing fixative.
11. To collect the histology sample, cut the tip off the cannula and then expel all of the contents using the plunger into a bottle containing fixative.

Sterilization

Ethylene oxide-sterilized device.

Storage conditions

Store in a dry environment protected from light.
Product shelf-life: 5 years.

For US only :

Caution: Federal law restricts this device to sale by or on the order of a physician.

PRODIMED
4 avenue de l'Europe - 60530
Neuilly en Thelle FRANCE

0459
(1987)

STERILE EO

PRODIMED
4 avenue de l'Europe - 60530
Neuilly en Thelle FRANCE

CE
0459
(1987)

STERILE EO