IFU-FCSSWAF-INTER-2.0

CE 0197

**Circular Stapler and Reloading Unit for Circular Stapler**

## Circular Stapler EN

FCSSWAF

**READ THE FOLLOWING INFORMATION THOROUGHLY BEFORE USING PRODUCT.**

### IMPORTANT NOTES

This booklet is NOT a reference to surgical techniques but a supporting material for the use of this product, which was designed, tested and produced merely for single patient use. It is NOT suggested that this device be reused or reprocessed or resterilized due to the risks of failure, patient injury, contamination as well as patient infection.

### DESCRIPTION

The FCSSWAF™ stapler places circular, double staggered rows of titanium staples in the tissue and resects the excess tissue, thus creating a circular anastomosis. The instrument can be activated when the handle is squeezed firmly to its limit. There are 5 stapler diameters in common use, 24 mm, 26 mm, 29 mm, 32 mm, or 34 mm.

### INDICATIONS

The FCSSWAF™ staplers can be applied throughout the alimentary tract to create end-to-end, end-to-side and side-to-side anastomoses.

### CONTRAINDICATIONS

1. This device is not designed, sold, or intended for use except as indicated.
2. The FCSSWAF™ stapler should NOT be used on tissues that would not tolerate conventional suture materials or conventional closure techniques if so judged by the surgeon.
3. It is NOT suitable to use the FCSSWAF™ stapler on any tissue that cannot compress comfortably or by less than the specified compression requirements. These instructions have to be followed; otherwise, closure failure, tissue trauma, dehiscence, tissue tearing, and displacement may occur, and/or hemostasis may not be obtained.
4. The instrument should NOT be used if unusual effort is required when turning the wing-nut for the visualization of the red line in the middle of the green bar in the indicator window.
5. If the cartridge is too large for the diameter of the structure and the tissue is stretched or thinned out, the FCSSWAF™ stapler should NOT be used. Leakage and narrowness may occur in the anastomosis. (Refer to the staple specification chart.)
6. Do NOT use the FCSSWAF™ stapler unless there is enough tissue to allow for proper inversion of tissue edges and secure placement for staples in the inverted edges. Purse-string suture has to be tied securely around the purse-string notch on the anvil.
7. Do NOT use the stapler if hemostasis cannot be visually proved after firing.

### WARNINGS AND PRECAUTIONS

1. Preoperative radiotherapy may result in tissue thickness to exceed the indicated range

for the selected staple size. Pre-surgical treatment over the patient deserves high attention. The staple size should be carefully selected.

2. Purse-string sutures must be placed less than 2.5 mm from the cut edge of the tissue, to avoid excessive tissue within the closed anvil and cartridge, which may result in staple malformation or leakage.

3. When piercing tissue with the instrument trocar, the point and puncture site have to be visualized all the time to prevent inadvertent damage to surrounding structures.

4. It has to be verified before attaching the anvil to the instrument that the Orange Band is fully visible. If the Orange Band cannot be seen, proper assembly of the instrument and anvil may not be well obtained.

5. It has to be made certain that the space between the cartridge and anvil is closed appropriately and the tissue is compressed. Before firing the instrument, Inspection must be made to ensure that the red line is in the middle of the green bar in indicator window.

6. To fire the instrument, the red safety below the handle has to be released first, and then the handle must be firmly squeezed to its limit. When firing the instrument, make sure that a reduced trigger pressure is felt and a "crunch" is heard as the instrument completes the firing cycle. Partial or incomplete handle squeezes may result in unacceptable staple formation and/or incomplete knife cut. This may result in intraoperative leaks.

**7. When opening the stapler prior to removal, open the instrument four to four and a half revolutions.**

8. After the removal of stapler, always inspect the staple line for hemostasis and the anastomosis must always be checked for integrity. Minor bleeding can be controlled by means of electrocautery or manual sutures.


9. Dispose of the instrument after tissue specimens are taken from the shaft.

10. The tissue specimens ("donuts") should be inspected to verify that all tissue layers have been incorporated in the anastomosis. If the specimens are not complete, a leak could result and cause narrowing.

11. Do not reuse any components of disposable instruments.

12. Application of an improperly matched instrument and anvil combination will lead to staple malformation or failure of proper cut. Malformed staples may compromise the integrity of the staple line, which may result in leakage or disruption.

13. This device is provided STERILE. Each instrument is designed for use in a SINGLE procedure only. DISCARD AFTER USE. DO NOT RESTERILIZATION.


## INSTRUCTIONS FOR USE

1. Take the stapler out of the package. Open the instrument by turning the wing-nut counterclockwise until the orange band is visible. Remove the protector.

2. Place purse-string sutures (A) in the organ to be anastomosed. The surgeon can alternatively use a closed lumen technique (double or triple stapling technique) (B) based on his/her experience and judgment.

3. If a double stapling technique is employed, open the instrument by turning the wing-nut until the orange band is visible.

4. Remove the detachable anvil to expose the trocar.

5. Turn the wing-nut clockwise to retract the trocar until a stop is reached.

6. Check the trocar to verify that it is fully retracted before the next step. **Note:** The instrument may also be inserted without removing the detachable anvil assembly if the preferred application is a sutured purse-string technique. In this case, however, prior to insertion, the detachable anvil assembly and the staple housing must be closed by rotating the wing-nut clockwise.

7. Insert the detachable anvil into the lumen and secure the purse-string onto the anvil shaft above the purse-string notch.

8. Insert the instrument up to the closed lumen with the trocar fully retracted. Extend the trocar fully and rotate the wing-nut counterclockwise to pierce tissue. Push the tissue down until the orange band is visible.

**Caution: Keep the trocar visible at all times to prevent personal injury or inadvertent trauma to adjacent structures.**

9. Reattach the detachable anvil by sliding the anvil shaft over the trocar and pushing until the detachable anvil snaps into its fully seated position.

**Caution: Do NOT clamp across or grip on the locking springs while attempting to reattach the detachable anvil to avoid inclusion of tissue within the anvil shaft.**

10. While closing the instrument, keep the organ segments in proper orientation. Inspect to ensure **no** extraneous tissue. Turn the wing-nut clockwise to close the instrument.

11. As the final adjusting rotation is being performed, the red line moves into the middle of the green bar of the indicator window. Applying gentle counter traction on the distal bowel during approximation may minimize excessive tissue incorporated into the barrel of the stapler.

**Caution: If unusual effort is required to turn the wing-nut for visualization of the red line in the middle of the green bar in the indicator window, DO NOT use the instrument.**

12. **Prefire checklist** : The Red line is in the middle of the green bar. The anvil is securely attached. To fire the instrument, pull the red safety back toward the wing-nut until it seats into the body of the instrument.

13. Once released, **squeeze the firing handle using firm, steady pressure**. The surgeon will feel reduced trigger pressure and hear a “crunch” as the instrument completes the firing cycle.

**Caution: Failure to squeeze the handle fully during firing may result in unacceptable staple formation and/or incomplete knife cut. This may result in intraoperative leaks. Ensure that the handle is compressed to the full extent of its travel.**


14. After firing, release the firing handle, and allow it to return to its original position. After the handle is released, return the safety to the locked position.

15. Turn the wing-nut counterclockwise 4 FULL to 4.5 TURNS to create a space between the cartridge and anvil to remove the instrument. Opening the wing-nut 4 full turns will allow

the anvil to tilt on its side for easy removal through the anastomosis. Pull the instrument straight

out of the new anastomosis, and gently remove it. Then inspect the staple line for hemostasis.

16. Remove the detachable anvil, washer (if present), and donuts from within the circular knife for donuts inspection. Examine the integrity of the donuts which should be intact with all tissue layers included. Incomplete donuts indicate risk of leakage in the anastomosis, and appropriate repairs should be made.

#### STAPLE SPECIFICATION CHART






Q) SIZE OF STAPLE

R) APPROXIMATE CLOSED STAPLE HEIGHT (HEIGHT OF "B") OF STAPLE

S) SHELL OUTSIDE DIAMETE

T) NUMBER OF STAPLES

U) TISSUE COMPRESSION REQUIREMENT

	Q	R	S	T	U
24	 4,5 mm	 1,8 mm	24 mm	18	1,8 mm
26	 4,8 mm	 2,0 mm	26 mm	20	2,0 mm
29	 4,8 mm	 2,0 mm	29 mm	24	2,0 mm
32	 5,0 mm	 2,2 mm	32 mm	30	2,2 mm
34	 5,0 mm	 2,2 mm	34 mm	32	2,2 mm

STORE AT ROOM TEMPERATURE.

AVOID PROLONGED EXPOSURE TO ELEVATED TEMPERATURES.

DO NOT EXPOSE TO TEMPERATURES ABOVE 130° F (54° C).

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