- de** Gebrauchsanweisung
- en** User's Manual
- fr** Mode d'emploi
- it** Istruzioni per l'uso
- es** Instrucciones de empleo
- pt** Manual de operação
- ru** Руководство по применению
- zh** 使用说明书

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QM-1-020L

**Macintosh**

**Miller**

**Foregger**

**KaWe MEGALIGHT®  
Macintosh**

**KaWe FLAPLIGHT®  
Macintosh**

**Economy Macintosh**

**Economy Miller**

<b>de</b>	Fiber Optik (F.O.) Laryngoskope.....	4
<b>en</b>	Fibre Optic Laryngoscopes (F.O.).....	12
<b>fr</b>	Laryngoscopes à fibre optique (F.O.).....	20
<b>it</b>	Laringoscopi a fibre ottiche (F.O.).....	28
<b>es</b>	Laringoscopios de fibra óptica (F.O.).....	36
<b>pt</b>	Laringoscópios de fibras ópticas (F.O.).....	44
<b>ru</b>	Ларингоскопы волоконная оптика (фиброоптика).....	52
<b>zh</b>	光纤(F.O.)喉镜 .....	60

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
## User's Manual Fibre-Optic (F.O.) Laryngoscopes

Dear Customer, we thank you for choosing a KaWe product. Our products are known for their high quality and long life. This KaWe products meets EC standards 93/42/EEC (standards for medical products).

**Please read these instructions thoroughly and carefully before attempting to use this products.**

Before attempting to use this product heed the given instructions.

**Utilization:** The laryngoscopes are only to be used by authorized personnel.


**Purpose:** Used for direct inspection of the larynx. These laryngoscopes are mainly used in cases of medical emergency for oral intubation of accident victims or for intubation before inducing anaesthesia during operations. This instrument

is used in connection with endotracheal tube in order to remove acute respiration disturbances and before operations in hospitals to ensure a clear respiratory tract. This instrument consists of 3 principal parts: the handle, the light source and the blade.

**Unsuited use/contraindication of the products:** Any use other than that described here is not in accordance with the intended use of the instrument. The manufacturer is not liable for any resulting damages. The user alone bears the risk. For the light source please exclusively use the KaWe laryngoscope handles.

**Risks and warning of possible dangers:** **1.** Danger of injury of the sensitive mucosa and the epiglottis during the intubation. **2.** During intubation, ensure that the lips, the teeth, the mucosa, the larynx and the vocal cords are not in danger of being harmed. In most cases, injuries to the lips and the tongue are most often caused when they are pinched between the patient's teeth and the handle of the laryngoscope. **3.** After use, the laryngoscope blade must

be disinfected and sterilized before further use. ( Disinfection in an autoclave is to be limited in order to avoid a shortened service life of the light guide due to the thermal and mechanical stress involved). **4.** After use, please clean the handles with an alcoholic disinfectant solution in order to avoid possible infection.

### Replacing the fibre-optic light guide:

**Please note that the KaWe MEGALIGHT® blades with integrated fibre-optics cannot be replaced!!**

### Replacing the lamp in the fibre-optic handle:

Hint: The satisfactory operation of the laryngoscopes can only be guaranteed when suitable xenon bulbs and/or LED lamps are used. Remove the head from the handle sleeve. Unscrew the inner enclosure from the inside part of the head. Next, unscrew and remove the lamp. If necessary, clean the glass bulb of the new lamp with alcohol. The glass bulb must be properly clean and free from fingerprints (free of grease). The new lamp must be screwed completely into the socket.

**Further information, maintenance, storage:** **1.** Familiarize yourself thoroughly with the assembly of the laryngoscopes before use. **2.** In the case of constant use and cleaning of the fibre-optic blades with exchangeable light guides it is recommended that the screw on the light guide be readjusted at regular intervals.

When used and stored properly, this instrument will serve you many years. Pay attention to the following instruction: Products with chromium-plated surface or stainless steel: Stained laryngoscope charging handles are the result of insufficient care.

**Reasons for stained laryngoscopes may be:** **1.** Insufficient mechanical or manual cleaning. **2.** Unsuitable cleaning materials or disinfectants. **3.** Non-observance of the dosing instructions with regard to cleaning materials, disinfectants or by residues (entrainment). **4.** Influences of foreign ions in the water, such as iron or silicate. **5.** Residues from medications, pens or chemical indicators. **6.** Errors of procedure (e.g. the new surgical instrument has not been cleaned before being sterilized). **7.** Instruments may not be stored in physiological salt solution in order to prevent crevice corrosion. **8.** Avoid contact with chemicals containing chloride. Risk of corrosion. **9.** Sterilization in a hot-air apparatus is not allowed.

**Attention: 3.5V battery handles that can be recharged in the socket (REF 03.41000.811), have a separate user's manual!**

**Guarantee:** When used under normal circumstances and with attention to these instructions (User's Manual), we guarantee this product for two years from the date of purchase (excluding wear parts such as batteries/accumulators and lamps). Should you have any questions or should your device require repair, please contact your dealer!

**Common Accessories:** Further information about accessories for these products can be found on our homepage at: [www.kawemed.de](http://www.kawemed.de).

Manufacturer **KaWe**

**Warning:**

**Limitation of reuse preparation:** Frequent preparation has a minimal effect on the instrument. The end of the product's life is normally due to wear and tear and damage by use. Ultrasonic cleaning is not allowed! Flash autoclave treatment and hot air sterilization will damage the blades! It is recommended, depending on its frequency of application, that the handle's head piece be replaced after intensive wear in order to ensure that the blade fits properly!

**General information:** The information and maintenance instructions of the manufacturer of the cleaning -/ disinfection unit (RDG) and of the cleaning disinfection agents are to be followed.

### **Instructions**

**Storage and transport:** It is recommended that the blades be immediately stored after their use in an instrument bowl filled with the appropriate cleaning/disinfection agents in order to prevent any residues from drying onto the material (protein fixation). It is recommended that the instruments be reprepared as soon as possible after their use. The instruments are to be transported to their place of treatment in a closed instrument bowl.

**Cleaning preparation:** It is recommended that the blades be immediately rinsed with running water or a mild alkaline solution immediately after their use in order to prevent any residues (such as blood) from drying onto the material. Satisfactory steam sterilization, as described below (see Cleaning, Sterilization) necessitates previous thorough cleaning

of the blades! Disassembling the instrument is not required. The laryngoscope blades are to be separated removed from the battery and/or charging handles.

**Manual cleaning:** In accordance with the Robert-Koch-Institute (RKI) it is recommended that the further preparation be preferentially carried out mechanically. Any manual cleaning will not be recommended.

**Automatic cleaning:** Equipment:

CM 310 (Maquet), cleaning agent (Neodisher FA forte 0,4 % / Neodisher Z 0.2%)

G7828 (Miele) cleaning agent (Mucapur XL 0.4 %/ Mucapur Z 0.5 %)

WD 390 (Belimed), cleaning agent (Mucapur AF 0.5 %/ Mucapur Z 0.1 %)

**1.** Prior to their mechanical preparation immediately rinse the instruments with running water in order to prevent any residues of the cleaning/disinfection agent from drying onto the machine. **2.** Place the instruments into an appropriate instrument stand. **3.** Place the instrument stand into the RDG in such a way that the full spray jet is directly centered on the instruments. **4.** Add the cleaning agents according to the instructions of the manufacturer. **5.** Start of the Vario TD program including thermal disinfection. The thermal disinfection is carried out in accordance with the A0-value and the national laws and provisions (EN/ISO 15883). **6.** At the end of program remove the instruments from the RDG and dry them (according to the RKI recommendations preferably with compressed air). The areas of the instrument stands that are difficult to reach are especially to be dried. **7.** Check for the intactness and cleanliness of the instrument with an appropriate magnifying glass (according to experience an eightfold magnification allows for an optical sight check). In the case that after the automatic preparation residual contaminations are to be found on the instrument, repeat the cleaning and the disinfection until no contaminations are still to be found. **8.** Caution! In the case of simple automatic cleaning (without evident disinfection) a terminal thermal disinfection in the steam-sterilizer – whereby the instruments are to be packaged in appropriate stands or sieve bowls – is required.

**Check and performance test:** All instruments: visual check for damage and wear. Keep all components clean to avoid influences on the electrical conductivity.

**Packaging:** Individual packaging: a standard polyethylene bag can be used. The size of the bag must correspond to the size of the instrument so that it closes properly.

**Sterilization:** A satisfactory sterilization of the blades is ensured with the following sterilization parameters: Fractionated prevacuum (3-fold), Sterilization temperature: min. 132 °C, max. 137 °C, Holding time: 3 min. minimum (complete cycle), Drying time: 10 min

**Storage:** Store the instruments in a dustproof, moisture-proof and recontamination-proof location.

**Additional information:** The actual legal regulations for the reprocessing of medical products must be heeded (see e.g. [www.rki.de](http://www.rki.de)). The effectivity of the above procedure for the processing of the instrument groups in question must be ensured by the manufacturer. The reprocessor must ensure that the realized reprocessing with equipment, materials and personnel will achieve the desired results. Therefore, routine checks of the validated automatical or standardized manual preparation processes will be required. Furthermore, all deviations from the realized processes and their possible adverse affects must be evaluated carefully by the reprocessor.

**Disposal:** The product is to be disposed of at a separate waste disposal location in accordance with local ordinances.

#### **Contacting the manufacturer:**

Address or phone number of the dealer or dial +49-7141-68188-0.

#### **User's Manual/ Maintenance Instructions for Charging Batteries NiMH/Li-Ion**

**Purpose:** The charging batteries are to be used as an energy source for the relevant KaWe products.

**Unsuited use/contraindication of the products:** Any use other than that described here is not in accordance with the intended use of the product. The manufacturer is not liable for any resulting damages. The user alone bears the risk.

**Possible risks:** **1.** Avoid short-circuiting the charging batteries. **2.** Never bring the charging batteries into contact with solutions.

**Further information, maintenance, storage:** The charging batteries can be cleaned with a damp cloth and then be wiped off with a dry fabric. Please ensure that no short-circuit is produced.

**Further maintenance instructions that may prolong the lifetime of the charging battery:** Please heed the following instructions: (The non-observance of these instructions may result in the accumulator losing its full capacity within a short period of time.)

NiMH technology can stimulate a lazy battery effect that is analogous to the classic memory effect, which is triggered by the energy that has not been completely spent prior to recharging the battery. This effect can be nearly eliminated by completely discharging the accumulator from time to time and then by recharging it two or three times in quick succession. Subsequently, the accumulator will regain its complete capability. In the case that the accumulator does no longer provide sufficient energy to operate the unit, immediately turn off the unit and charge the accumulator (in order to prevent overcharging and the resulting destruction of the accumulator). Please keep in mind that the charging batteries can spontaneously discharge. When the charging batteries are not used (charged) for a longer period of time, a self-discharge may occur. The capacity will be lost and the accumulators may become permanently defective. The manufacturer is not liable for any resulting damages. Therefore, when stored for more than one month, the accumulators should be partially charged to approx. 50% .

Please always respect the given charging time in order to prevent overcharging and the resulting destruction of the accumulator (the charging station KaWe MedCharge® 3000 and 4000 has an integrated automatic shut-off). NiMH accumulators with low temperatures show a loss of capacity (already at temperatures of 0°C).

**1.** Keep the accumulators away from heating installations and open flames. **2.** Accumulators must not be dipped into water or salt water. **3.** Never try to disassemble or modify the accumulators or expose them to heat. **4.** Avoid dropping the accumulators and protect them against being hit or crushed as this may damage the housing. **5.** Please only use the accumulators and accessory parts recommended by KaWe. **6.** When using the accumulators for the first time, please make sure to charge them completely.

The use of accumulators that are not explicitly approved for these units may cause explosions or leakages and the risk of fire, injuries and damages to surroundings. In the case that an accumulator obtains a leak and spilled substances come into contact with eyes, mouth, skin or clothing, immediately rinse the relevant spot with water and consult a doctor.

**1.** When charging is completed and the charging unit and the compact adapter are no longer in use, please separate them from the unit and the socket in order to avoid fires and other risks. **2.** Do not cover the charging unit with any kind of material (table clothes, carpets, bedclothes, pillows etc.) during operation. If you use the adapter for a longer period of time, it may become overheated, lose its original shape or catch fire. **3.** Please only charge the accumulator in the appropriate accumulator charging unit. The charging unit and the compact adapter have been developed especially for being used with your product – please do not use them with other products or accumulators. All other products may cause the risk of excessive heat and distortion and thus risk of fire and electric shock. **4.** Prior to disposal of the accumulator, ensure that the contacts are masked with adhesive strips or other insulating material in order to prevent their direct contact with other matters. Contact with other metallic material in trash cans may cause the risk of fire and explosions.

### Charging time:

REF	12.80110.712	12.80110.722	12.80120.712	12.80220.722
Capacity [mAh]	1600	3000	750	2200
Charging time [h]	approx. 8*	approx. 14*	approx. 4*	approx. 9*

\*Charging time after complete discharge of the accumulator in the charging station KaWe MedCharge® 3000 (not Lilon) or 4000.

**Disposal:** Dispose of the lithium-ion batteries at a separate collection site for electric and electronic devices.

**Disposal:** Please dispose of expired NiMH rechargeable batteries in a manner that is compliant with legal regulations.

### Explanation of the symbols:

	Manufacturer		Separate disposal of electric and electronic instruments
	Date of manufacture		Disposal of batteries
	Please heed the user's manual		Meets relevant EU standards
	Warning		

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**KIRCHNER & WILHELM** GmbH + Co. KG