
Platelia *Aspergillus* Ag

1 khay – ▾ 96

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1. **WHO fungal priority pathogens list to guide research, development and public health action.** Geneva: World Health Organization. Licence: CC BY-NC-SA 3.0 IGO; 2022.
2. Stynen D, Goris A, Sarfati J, Latgé JP: **A new sensitive sandwich enzyme-linked immunosorbent assay to detect galactofuran in patients with invasive aspergillosis.** *J Clin Microbiol* 1995, **33**:497–500.
3. Desai R, Ross LA, Hoffman JA: **The role of bronchoalveolar lavage galactomannan in the diagnosis of pediatric invasive aspergillosis.** *Pediatr Infect Dis J* 2009, **28**:283–286.
4. Agarwal R, Aggarwal AN, Sehgal IS, Dhoria S, Behera D, Chakrabarti A: **Performance of serum galactomannan in patients with allergic bronchopulmonary aspergillosis.** *Mycoses* 2015, **58**:408–412.
5. Çağlar İ, Özkerim D, Tahta N, Düzgöl M, Bayram N, Demirağ B, Hilkey Karapınar T, Sorguç Y, Gözmen S, Dursun V, et al.: **262. Assessment of Serum Galactomannan Test Results of Pediatric Patients with Hematologic Malignancies According to Different Threshold Levels and Consecutive Positivity in Terms of Invasive Aspergillosis Diagnosis: Cross-Sectional Research in a Tertiary Care Hospital.** *Open Forum Infect Dis* 2019, **6**:S145–S146.
6. Couchepin J, Brunel A, Jaton K, Meylan P, Bochud P, Lamothe F: **Role of bi-weekly serum galactomannan screening for the diagnosis of invasive aspergillosis in haematological cancer patients.** *Mycoses* 2018, **61**:350–354.
7. Okuturlar Y, Ozkalemkas F, Ener B, Serin SO, Kazak E, Ozcelik T, Ozkocaman V, Ozkan HA, Akalin H, Gunaldi M, et al.: **Serum galactomannan levels in the diagnosis of invasive aspergillosis.** *Korean J Intern Med* 2015, **30**:899–905.
8. Warris A, Lehrnbecher T, Roilides E, Castagnola E, Brüggemann RJM, Groll AH: **ESCMID-ECMM guideline: diagnosis and management of invasive aspergillosis in neonates and children.** *Clinical Microbiology and Infection* 2019, **25**:1096–1113.
9. Dichtl K, Seybold U, Ormanns S, Horns H, Wagener J: **Comparison of a novel galactomannoprotein ELISA and the Platelia galactomannan ELISA for the diagnosis of invasive aspergillosis.** *Mycoses* 2019, **62**:15.
10. Wang X, Guo G, Cai R, He P, Zhang M: **Utility of serum galactomannan antigen testing combined with chest computed tomography for early diagnosis of invasive pulmonary aspergillosis in patients with hematological malignancies with febrile neutropenia after antifungal drug treatment.** *J Int Med Res* 2019, **47**:783–790.

11. Mikulska M, Furfaro E, De Carolis E, Drago E, Pulzato I, Borghesi ML, Zappulo E, Raiola AM, Grazia CD, Del Bono V, et al.: **Use of *Aspergillus fumigatus* real-time PCR in bronchoalveolar lavage samples (BAL) for diagnosis of invasive aspergillosis, including azole-resistant cases, in high risk haematology patients: the need for a combined use with galactomannan.** *Medical Mycology* 2019, **57**:987–996.
12. Badiie P, Hashemizadeh Z, Ramzi M, Karimi M, Mohammadi R: **Non-Invasive Methods to Diagnose Fungal Infections in Pediatric Patients with Hematologic Disorders.** *Jundishapur J Microbiol* 2016, **9**.
13. Bellanger AP, Millon L, Berceanu A, Grenouillet F, Grenouillet FE, Larosa F, Deconinck E: **Combining *Aspergillus* mitochondrial and ribosomal QPCR, in addition to galactomannan assay, for early diagnosis of invasive aspergillosis in hematology patients.** *Med Myco* 2015, **53**:760–764.
14. Wheat LJ, Nguyen MH, Alexander BD, Denning D, Caliendo AM, Lyon GM, Baden LR, Marty FM, Clancy C, Kirsch E, et al.: **Long-Term Stability at –20°C of *Aspergillus* Galactomannan in Serum and Bronchoalveolar Lavage Specimens.** *J Clin Microbiol* 2014, **52**:2108–2111.
15. Da Silva TVM, Carneiro LC, Dos Santos Ramos F, Baethgen LF, Paz AA, Larentis DZ, Daudt LE, Tusset C, Lopes-Bezerra LM, Pasqualotto AC: **PCR as a Screening Test for Invasive Aspergillosis in Haematological Patients: A Pilot Study.** *Mycopathologia* 2014, **177**:111–114.
16. Taremi M, Kleinberg ME, Wang EW, Gilliam BL, Ryscavage PA: **Galactomannan antigen detection using bronchial wash and bronchoalveolar lavage in patients with hematologic malignancies.** *Ann Clin Microbiol Antimicrob* 2015, **14**:50.
17. Zheng J, Gui X, Cao Q, Yang R, Yan Y, Deng L, Lio J: **A Clinical Study of Acquired Immunodeficiency Syndrome Associated *Penicillium Marneffe* Infection from a Non-Endemic Area in China.** *PLoS ONE* 2015, **10**:e0130376.
18. Singh N, Paterson DL: ***Aspergillus* Infections in Transplant Recipients.** *Clin Microbiol Rev* 2005, **18**:44–69.
19. Barnes PD, Marr KA: **Risks, diagnosis and outcomes of invasive fungal infections in haematopoietic stem cell transplant recipients.** *Br J Haematol* 2007, **139**:519–531.
20. Weinbergerova B, Kocmanova I, Racil Z, Mayer J: **Serological Approaches.** In *Human Fungal Pathogen Identification*. Edited by Lion T. Springer New York; 2017:209–221.
21. Gefen A, Zaidman I, Shachor-Meyouhas Y, Avidor I, Hakim F, Weyl Ben-Arush M, Kassis I: **Serum Galactomannan Screening for Diagnosis of Invasive Pulmonary Aspergillosis in Children After Stem Cell Transplantation or with High-Risk Leukemia.** *Pediatric Hematology and Oncology* 2015, **32**:146–152.
22. Denning DW: **Invasive aspergillosis.** *Clin Infect Dis* 1998, **26**:781–803; quiz 804–805.
23. Ramanan P, Wengenack NL, Theel ES: **Laboratory Diagnostics for Fungal Infections.** *Clinics in Chest Medicine* 2017, **38**:535–554.
24. Badiie P, Amirghofran AA, Ghazi Nour M: **Evaluation of noninvasive methods for the diagnosis of fungal endocarditis.** *Medical Mycology* 2014, **52**:530–536.
25. de Azevedo PZ, Sylvestre TF, Cavalcante R de S, de Carvalho LR, Moris DV, de Oliveira MLCS, Mendes RP: **Evaluation of the Double Agar Gel Immunodiffusion Test and of the Enzyme Linked Immunosorbent Assay in the Diagnosis and Follow-Up of Patients with Chronic Pulmonary Aspergillosis.** *PLoS One* 2015, **10**:e0134841.
26. Leelahavanichkul A, Pongpirul K, Thongbor N, Worasilchai N, Petphuak K, Thongsawang B, Towannang P, Lorvinitnun P, Sukhontasing K, Katavetin P, et al.: **(1→3)-β-d-Glucan and Galactomannan for Differentiating Chemical “Black Particles” and Fungal Particles Inside Peritoneal Dialysis Tubing.** *Perit Dial Int* 2016, **36**:402–409.
27. Knoth H, Maywald D, Walter W: **In-vitro detection of mannan and galactomannan in components of total parenteral nutrition (TPN).** *Pharmazie* 2016, doi:10.1691/ph.2016.5797.

28. De Pauw B, Walsh TJ, Donnelly JP, Stevens DA, Edwards JE, Calandra T, Pappas PG, Maertens J, Lortholary O, Kauffman CA, et al.: **Revised definitions of invasive fungal disease from the European Organization for Research and Treatment of Cancer/Invasive Fungal Infections Cooperative Group and the National Institute of Allergy and Infectious Diseases Mycoses Study Group (EORTC/MSG) Consensus Group.** *Clin Infect Dis* 2008, **46**:1813–1821.
29. Maertens J, Verhaegen J, Lagrou K, Van Eldere J, Boogaerts M: **Screening for circulating galactomannan as a noninvasive diagnostic tool for invasive aspergillosis in prolonged neutropenic patients and stem cell transplantation recipients: a prospective validation.** *Blood* 2001, **97**:1604–1610.
30. Yoshimura K, Suzuki Y, Inoue Y, Nishimoto K, Mori K, Karayama M, Hozumi H, Furuhashi K, Enomoto N, Fujisawa T, et al.: **Utility of serum Aspergillus-galactomannan antigen to evaluate the risk of severe acute exacerbation in chronic obstructive pulmonary disease.** *PLoS ONE* 2018, **13**:e0198479.
31. Moura S, Cerqueira L, Almeida A: **Invasive pulmonary aspergillosis: current diagnostic methodologies and a new molecular approach.** *Eur J Clin Microbiol Infect Dis* 2018, **37**:1393–1403.
32. Cabana ÂL, Mendes JF, Klafke GB, Brandolt TM, Melo AM, Meireles MCA, Xavier MO: **Can Aspergillus fumigatus conidia cause false-positive results in the galactomannan enzyme immunoassay test?** *Rev Soc Bras Med Trop* 2018, **51**:387–389.
33. De Repentigny L, Kaufman L, Cole GT, Kruse D, Latgé JP, Matthews RC: **Immunodiagnosis of invasive fungal infections.** *J Med Vet Mycol* 1994, **32 Suppl 1**:239–252.
34. Erjavec Z, Verweij PE: **Recent progress in the diagnosis of fungal infections in the immunocompromised host.** *Drug Resistance Updates* 2002, **5**:3–10.
35. Latgé JP: **Tools and trends in the detection of Aspergillus fumigatus.** *Curr Top Med Mycol* 1995, **6**:245–281.
36. Verweij PE, Meis JFGM: **Microbiological diagnosis of invasive fungal infections in transplant recipients: Diagnosis of invasive fungal infections.** *Transplant Infectious Disease* 2000, **2**:80–87.
37. Yeo SF, Wong B: **Current Status of Nonculture Methods for Diagnosis of Invasive Fungal Infections.** *Clin Microbiol Rev* 2002, **15**:465–484.
38. Prattes J, Heldt S, Eigl S, Hoenigl M: **Point of Care Testing for the Diagnosis of Fungal Infections: Are We There Yet?** *Curr Fungal Infect Rep* 2016, **10**:43–50.
39. Maertens J, Maertens V, Theunissen K, Meersseman W, Meersseman P, Meers S, Verbeken E, Verhoef G, Van Eldere J, Lagrou K: **Bronchoalveolar Lavage Fluid Galactomannan for the Diagnosis of Invasive Pulmonary Aspergillosis in Patients with Hematologic Diseases.** *CLIN INFECT DIS* 2009, **49**:1688–1693.
40. Salonen J, Lehtonen OP, Teräsjärvi MR, Nikoskelainen J: **Aspergillus antigen in serum, urine and bronchoalveolar lavage specimens of neutropenic patients in relation to clinical outcome.** *Scand J Infect Dis* 2000, **32**:485–490.
41. Sanguinetti M, Posteraro B, Pagano L, Pagliari G, Fianchi L, Mele L, La Sorda M, Franco A, Fadda G: **Comparison of real-time PCR, conventional PCR, and galactomannan antigen detection by enzyme-linked immunosorbent assay using bronchoalveolar lavage fluid samples from hematology patients for diagnosis of invasive pulmonary aspergillosis.** *J Clin Microbiol* 2003, **41**:3922–3925.
42. Verweij PE, Latgé JP, Rijs AJ, Melchers WJ, De Pauw BE, Hoogkamp-Korstanje JA, Meis JF: **Comparison of antigen detection and PCR assay using bronchoalveolar lavage fluid for diagnosing invasive pulmonary aspergillosis in patients receiving treatment for hematological malignancies.** *J Clin Microbiol* 1995, **33**:3150–3153.
43. Meersseman W, Lagrou K, Maertens J, Wilmer A, Hermans G, Vanderschueren S, Spriet I, Verbeken E, Van Wijngaerden E: **Galactomannan in Bronchoalveolar Lavage Fluid: A Tool for Diagnosing Aspergillosis in Intensive Care Unit Patients.** *Am J Respir Crit Care Med* 2008, **177**:27–34.

44. Musher B, Fredricks D, Leisenring W, Balajee SA, Smith C, Marr KA: **Aspergillus galactomannan enzyme immunoassay and quantitative PCR for diagnosis of invasive aspergillosis with bronchoalveolar lavage fluid.** *J Clin Microbiol* 2004, **42**:5517–5522.
45. Sulahian A, Boutboul F, Ribaud P, Leblanc T, Lacroix C, Derouin F: **Value of antigen detection using an enzyme immunoassay in the diagnosis and prediction of invasive aspergillosis in two adult and pediatric hematology units during a 4-year prospective study.** *Cancer* 2001, **91**:311–318.
46. Becker MJ, Lugtenburg EJ, Cornelissen JJ, Van Der Schee C, Hoogsteden HC, De Marie S: **Galactomannan detection in computerized tomography-based broncho-alveolar lavage fluid and serum in haematological patients at risk for invasive pulmonary aspergillosis.** *Br J Haematol* 2003, **121**:448–457.
47. Kwak EJ, Nguyen MH: **Galactomannan detection in bronchoalveolar lavage fluid.** *Curr Fungal Infect Rep* 2008, **2**:206–213.
48. Clancy CJ, Jaber RA, Leather HL, Wingard JR, Staley B, Wheat LJ, Cline CL, Rand KH, Schain D, Baz M, et al.: **Bronchoalveolar Lavage Galactomannan in Diagnosis of Invasive Pulmonary Aspergillosis among Solid-Organ Transplant Recipients.** *J Clin Microbiol* 2007, **45**:1759–1765.
49. Lahmer T, Neuenhahn M, Held J, Rasch S, Schmid RM, Huber W: **Comparison of 1,3- β -d-glucan with galactomannan in serum and bronchoalveolar fluid for the detection of Aspergillus species in immunosuppressed mechanical ventilated critically ill patients.** *Journal of Critical Care* 2016, **36**:259–264.
50. Husain S, Clancy CJ, Nguyen MH, Swartzentruber S, Leather H, LeMonte AM, Durkin MM, Knox KS, Hage CA, Bentsen C, et al.: **Performance Characteristics of the Platelia Aspergillus Enzyme Immunoassay for Detection of Aspergillus Galactomannan Antigen in Bronchoalveolar Lavage Fluid.** *Clin Vaccine Immunol* 2008, **15**:1760–1763.
51. Paterson DL, Singh N: **Invasive aspergillosis in transplant recipients.** *Medicine (Baltimore)* 1999, **78**:123–138.
52. Latgé JP, Kobayashi H, Debeaupuis JP, Diaquin M, Sarfati J, Wieruszeski JM, Parra E, Bouchara JP, Fournet B: **Chemical and immunological characterization of the extracellular galactomannan of Aspergillus fumigatus.** *Infect Immun* 1994, **62**:5424–5433.
53. Stynen D, Sarfati J, Goris A, Prévost MC, Lesourd M, Kamphuis H, Darras V, Latgé JP: **Rat monoclonal antibodies against Aspergillus galactomannan.** *Infect Immun* 1992, **60**:2237–2245.
54. Wheat LJ, Walsh TJ: **Diagnosis of invasive aspergillosis by galactomannan antigenemia detection using an enzyme immunoassay.** *Eur J Clin Microbiol Infect Dis* 2008, **27**:245–251.
55. Latgé J-P: **Aspergillus fumigatus and Aspergillosis.** *Clin Microbiol Rev* 1999, **12**:310–350.
56. Maertens JA, Klont R, Masson C, Theunissen K, Meersseman W, Lagrou K, Heinen C, Crepin B, Eldere JV, Tabouret M, et al.: **Optimization of the Cutoff Value for the Aspergillus Double-Sandwich Enzyme Immunoassay.** *Clinical Infectious Diseases* 2007, **44**:1329–1336.
57. Verweij PE, Stynen D, Rijs AJ, De Pauw BE, Hoogkamp-Korstanje JA, Meis JF: **Sandwich enzyme-linked immunosorbent assay compared with Pastorex latex agglutination test for diagnosing invasive aspergillosis in immunocompromised patients.** *J Clin Microbiol* 1995, **33**:1912–1914.
58. Gangneux J-P, Lavarde D, Bretagne S, Guiguen C, Gandemer V: **Transient aspergillus antigenaemia: think of milk.** *Lancet* 2002, **359**:1251.
59. Herbrecht R, Letscher-Bru V, Oprea C, Lioure B, Waller J, Campos F, Villard O, Liu K-L, Natarajan-Amé S, Lutz P, et al.: **Aspergillus Galactomannan Detection in the Diagnosis of Invasive Aspergillosis in Cancer Patients.** *JCO* 2002, **20**:1898–1906.
60. Mennink-Kersten MASH, Ruegebrink D, Klont RR, Warris A, Gavini F, Op Den Camp HJM, Verweij PE: **Bifidobacterial Lipoglycan as a New Cause for False-Positive Platelia Aspergillus Enzyme-Linked Immunosorbent Assay Reactivity.** *J Clin Microbiol* 2005, **43**:3925–3931.

61. Siemann M, Koch-Dörfler M, Gaude M: **False-positive results in premature infants with the Platelia® Aspergillus sandwich enzyme-linked immunosorbent assay.** *Mycoses* 1998, **41**:373–377.
62. Verweij PE, Weemaes CM, Curfs JHAJ, Bretagne S, Meis JFGM: **Failure To Detect Circulating Aspergillus Markers in a Patient with Chronic Granulomatous Disease and Invasive Aspergillosis.** *J Clin Microbiol* 2000, **38**:3900–3901.
63. Schroeder I, Dichtl K, Liebchen U, Wagener J, Irlbeck M, Zoller M, Scharf C: **Digestive enzymes of fungal origin as a relevant cause of false positive Aspergillus antigen testing in intensive care unit patients.** *Infection* 2021, **49**:241–248.
64. De Heer K, Gerritsen MG, Visser CE, Leeflang MM: **Galactomannan detection in bronchoalveolar lavage fluid for invasive aspergillosis in immunocompromised patients.** *Cochrane Database of Systematic Reviews* 2019, **2020**.
65. Marr KA, Laverdiere M, Gugel A, Leisenring W: **Antifungal Therapy Decreases Sensitivity of the Aspergillus Galactomannan Enzyme Immunoassay.** *Clinical Infectious Diseases* 2005, **40**:1762–1769.
66. Knox KS, Rose AS, Hage CA: **Rapid Fungal Diagnosis: The Utility of Bronchoalveolar Lavage for Pneumocystis and Endemic Mycoses.** *Current Fungal Infection Reports* 2007,
67. Narreddy S, Chandrasekar PH: **False-positive Aspergillus galactomannan (GM) assay in histoplasmosis.** *Journal of Infection* 2008, **56**:80–81.
68. Swanink CM, Meis JF, Rijs AJ, Donnelly JP, Verweij PE: **Specificity of a sandwich enzyme-linked immunosorbent assay for detecting Aspergillus galactomannan.** *J Clin Microbiol* 1997, **35**:257–260.
69. Wheat LJ, Hackett E, Durkin M, Connolly P, Petraitiene R, Walsh TJ, Knox K, Hage C: **Histoplasmosis-Associated Cross-Reactivity in the BioRad Platelia Aspergillus Enzyme Immunoassay.** *Clin Vaccine Immunol* 2007, **14**:638–640.
70. Chambon-Pautas C, Costa J-M, Chaumette M-T, Cordonnier C, Bretagne S: **Galactomannan and Polymerase Chain Reaction for the Diagnosis of Primary Digestive Aspergillosis in a Patient with Acute Myeloid Leukaemia.** *Journal of Infection* 2001, **43**:213–214.
71. Ansorg R, van den Boom R, Rath PM: **Detection of Aspergillus galactomannan antigen in foods and antibiotics.** *Mycoses* 1997, **40**:353–357.
72. Letscher-Bru V, Cavalier A, Pernot-Marino E, Koenig H, Eyer D, Waller J, Candolfi E: **Recherche d'antigène galactomannane Aspergillaire circulant par Platelia® Aspergillus: Antigenémies positives persistantes en l'absence d'infection.** *J Mycol Med* 1998, **8**:112–113.
73. Maertens J, Verhaegen J, Demuyneck H, Brock P, Verhoef G, Vandenberghe P, Van Eldere J, Verbist L, Boogaerts M: **Autopsy-Controlled Prospective Evaluation of Serial Screening for Circulating Galactomannan by a Sandwich Enzyme-Linked Immunosorbent Assay for Hematological Patients at Risk for Invasive Aspergillosis.** *J Clin Microbiol* 1999, **37**:3223–3228.
74. Verdaguer V, Walsh TJ, Hope W, Cortez KJ: **Galactomannan antigen detection in the diagnosis of invasive aspergillosis.** *Expert Review of Molecular Diagnostics* 2007, **7**:21–32.
75. Blijlevens NMA, Donnelly JP, Meis JFGM, Verweij PE, de Pauw BE: **Aspergillus galactomannan antigen levels in allogeneic haematopoietic stem cell transplant recipients given total parenteral nutrition.** *Transpl Infect Dis* 2002, **4**:64–65.
76. Ortiz-Brizuela E, Rangel-Cordero A, Ponce-de-León A, Sierra-Madero J: **False-positive results in the galactomannan Platelia™ Aspergillus assay with generic piperacillin/tazobactam.** *Revista Iberoamericana de Micología* 2019, **36**:51–52.
77. Mikulska M, Furfaro E, Viscoli C: **Non-cultural methods for the diagnosis of invasive fungal disease.** *Expert Review of Anti-infective Therapy* 2015, **13**:103–117.

78. Sulahian A, Porcher R, Bergeron A, Touratier S, Raffoux E, Menotti J, Derouin F, Ribaud P: **Use and Limits of (1-3)- β -D-Glucan Assay (Fungitell), Compared to Galactomannan Determination (Platelia Aspergillus), for Diagnosis of Invasive Aspergillosis.** *J Clin Microbiol* 2014, **52**:2328–2333.
79. Aubry A, Porcher R, Bottero J, Touratier S, Leblanc T, Brethon B, Rousselot P, Raffoux E, Menotti J, Derouin F, et al.: **Occurrence and kinetics of false-positive Aspergillus galactomannan test results following treatment with beta-lactam antibiotics in patients with hematological disorders.** *J Clin Microbiol* 2006, **44**:389–394.
80. Vergidis P, Razonable RR, Wheat LJ, Estes L, Caliendo AM, Baden LR, Wingard JR, Baddley J, Assi M, Norris S, et al.: **Reduction in False-Positive Aspergillus Serum Galactomannan Enzyme Immunoassay Results Associated with Use of Piperacillin-Tazobactam in the United States.** *J Clin Microbiol* 2014, **52**:2199–2201.
81. Racil Z, Kocmanova I, Lengerova M, Winterova J, Mayer J: **Intravenous PLASMA-LYTE as a Major Cause of False-Positive Results of Platelia Aspergillus Test for Galactomannan Detection in Serum.** *Journal of Clinical Microbiology* 2007, **45**:3141–3142.
82. Surmont I, Stockman W: **Gluconate-Containing Intravenous Solutions: Another Cause of False-Positive Galactomannan Assay Reactivity.** *J Clin Microbiol* 2007, **45**:1373–1373.
83. Nguyen MH, Jaber R, Leather HL, Wingard JR, Staley B, Wheat LJ, Cline CL, Baz M, Rand KH, Clancy CJ: **Use of Bronchoalveolar Lavage To Detect Galactomannan for Diagnosis of Pulmonary Aspergillosis among Nonimmunocompromised Hosts.** *J Clin Microbiol* 2007, **45**:2787–2792.
84. Upton A, Gugel A, Leisenring W, Limaye A, Alexander B, Hayden R, Marr KA: **Reproducibility of Low Galactomannan Enzyme Immunoassay Index Values Tested in Multiple Laboratories.** *J Clin Microbiol* 2005, **43**:4796–4800.
85. Penack O, Rempf P, Graf B, Blau IW, Thiel E: **Aspergillus galactomannan testing in patients with long-term neutropenia: implications for clinical management.** *Annals of Oncology* 2008, **19**:984–989.
86. **NCCLS. Evaluation of Precision Performance of Quantitative Measurement Methods; Approved Guideline-Second Edition.** NCCLS document EP5-A2 [ISBN 1-56238-542-9]. NCCLS, 940 West Valley Road, Suite 1400, Wayne, Pennsylvania 19087-1898 USA; 2004.
87. Ascioğlu S, Rex JH, de Pauw B, Bennett JE, Bille J, Crokaert F, Denning DW, Donnelly JP, Edwards JE, Erjavec Z, et al.: **Defining opportunistic invasive fungal infections in immunocompromised patients with cancer and hematopoietic stem cell transplants: an international consensus.** *Clin Infect Dis* 2002, **34**:7–14.
88. Donnelly JP, Chen SC, Kauffman CA, Steinbach WJ, Baddley JW, Verweij PE, Clancy CJ, Wingard JR, Lockhart SR, Groll AH, et al.: **Revision and Update of the Consensus Definitions of Invasive Fungal Disease From the European Organization for Research and Treatment of Cancer and the Mycoses Study Group Education and Research Consortium.** *Clinical Infectious Diseases* 2020, **71**:1367–1376.

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(EN)

Danger

Causes severe skin burns and eye damage. May cause an allergic skin reaction. Harmful to aquatic life with long lasting effects. Wear protective gloves/protective clothing/eye protection/face protection. IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing. IF SWALLOWED: rinse mouth. Do NOT induce vomiting. IF ON SKIN (or hair): Remove/ Take off immediately all contaminated clothing. Rinse skin with water/shower. If skin irritation or rash occurs: Get medical advice/attention. Avoid release to the environment. Dispose of contents/ container in accordance with local/regional/national/international regulations.

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BIO-RAD là nhãn hiệu của Bio-Rad Laboratories, Inc.
EVOLIS, EVOLIS TWIN PLUS và PLATELIA là nhãn hiệu của Bio-Rad Europe, GmbH tại một số khu vực tài phán.