# HEMATOLOGY REAGENT FOR CD 1800,1700,1600,1400,1300 ANALYZERS

## 【INTENDED USE】

To be used for blood cell counting, sizing, and WBC differentiation, hemoglobin determination on ABBOTT CD 1800,1700,1600,1400,1300 hematology analyzers. This product is for In Vitro Diagnostic Use only.

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30-10-020	Diluent CD3	20L
33-20-020	Detergent CD3	20L
33-30-004	Lyse CD3	3.8L
33-40-E001	EZ-Cleaner	50ml×2

## 【KIT COMPOSITION】

DILUENT CD3	NaCl≤9.4%, Buffer≤1.6%, Anti-Microbial Agent≤1.0%
DETERGENT CD3	Sodium Sulfate<1.5%, Anti-microbial Agent<0.1%
LYSE CD3	Quaternary Ammonium Salt<10%, Sodium Lauryl Sulfate<1.5%
EZ-CLEANER	Subtilisin<5.0%

## 【APPLICABLE INSTRUMENT】

For use with CD1800,1700,1600,1400,1300 hematology analyzers.

## 【PREPARATION】

Reagent is ready to use.

## 【PROCEDURES】

- Using the function keys press the following menus:  
MAIN-SPECIAL                      PROTOCOLS-REAGENT  
RESERVOIRE-EMPTY DILUENT
- Remove the reagent inlet tube from the Diluent.
- Wipe the tube's surface with a lint free tissue, and then rinse with distilled water.
- Put the reagent tubes into the boxes and/or bottles of Labnovation reagents.
- Using the function keys press the following menu:  
FILL DILUENT  
This step will fill the reagent system with Labnovation reagent.
- Repeat steps 1-5 with Lyse and Sheath.
- Run at least 5 backgrounds. Check for the acceptable values.
- Run all three levels of a dedicated commercial control. Check for the acceptable values.
- Recalibrate the unit if necessary.  
(MAIN-CALIBRATION-ENTER-FACTOR)

## 【PRECAUTIONS】

- Re-calibration is recommended when changing the origin of the reagents.
- It is not recommended using different brands of solutions together on one unit of analyzer.
- Do not ingest. Avoid contact with skin and eyes.
- Skin contact: Flush with plenty of water.
- Eyes contact: Flush with plenty of water. Obtain immediate medical attention.
- This kit is for professionals and in vitro diagnostic use only.
- Never put distilled water on the analyser (except clean for shipping menu), otherwise the pneumatic unit can be damaged seriously.

## 【STORAGE】

30-10-020	Diluent CD3	24 months
33-20-020	Detergent CD3	24 months
33-30-004	Lyse CD3	24 months
33-40-E001	EZ-Cleaner	12months

Store the product at 5-35 °C in a dust free place with good

ventilation. Temperature outside the range may damage the product.

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# HEMATOLOGY REAGENT FOR NIHON KHODEN 3 PART DIFF ANALYZERS

## 【INTENDED USE】

To be used for blood cell counting, sizing, and WBC differentiation, hemoglobin determination on Nihon Khoden 3 part diff hematology analyzers. This product is for In Vitro Diagnostic Use only.

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41-10-020	Diluent NK	20L
41-30-A05	Lyse NK	500ml×3
41-40-005R	Rinse NK	5L
L-0401005	Cleaner-C	500ml

## 【KIT COMPOSITION】

DILUENT NK	NaCl≤9.4%, Buffer≤1.6%, Anti-Microbial Agent≤1.0%
CLEANER-C	Sodium Hypochlorite 5%
LYSE NK	Quaternary Ammonium Salt<10%, Sodium Lauryl Sulfate<1.5%
RINSE NK	Buffer<0.9%, Sufurctant<1.0%

## 【APPLICABLE INSTUMENT】

For use with Nihon Kohden MEK-6318, MEK-6400, MEK-6700, MEK-6410, MEK-6420 hematology analyzers.

## 【PREPARATION】

Reagent is ready to use.

## 【PROCEDURES】

1. Remove all the reagent aspiration tubes from the currently used reagent boxes.
2. Wipe the tube's surface with a lint free tissue, and then rinse with distilled water. Put the tubes into Labnovation reagents boxes or bottles.
3. Run at least 50 backgrounds. Check for the acceptable background values. At first background counts can be seen, that should lower to normal after a few runs.
4. Run all three levels of a dedicated commercial control. Check for the acceptable values.
5. Recalibrate the unit if necessary.

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Store the product at 5-35 °C in a dust free place with good ventilation. Temperature outside the range may damage the product.

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# HEMATOLOGY REAGENT FOR SYSMEX 3 PART DIFF ANALYZERS

## 【INTENDED USE】

To be used for blood cell counting, sizing, and WBC differentiation, hemoglobin determination on SYSMEX KX-21 hematology analyzers. This product is for In Vitro Diagnostic Use only.

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20-10-020	Diluent-ST	20L
21-30-SWH	Lyser-WH	500ml×3
21-CL-50	Cell-Cleaner	50ml

## 【KIT COMPOSITION】

DILUENT-ST	NaCl≤9.4%, Buffer≤1.6%, Anti-Microbial Agent≤1.0%
LYSER-WH	Quaternary Ammonium Salt<10%, Sodium Lauryl Sulfate<1.5%
CELL-CLEANER	Sodium Hypochlorite 5%

## 【APPLICABLE INSTRUMENT】

For use with Sysmex KX-21, KX-21N, POCH-100i, XP-300, XP-100 hematology analyzers.

## 【PREPARATION】

Reagent is ready to use.

## 【PROCEDURES】

1. Remove the diluents tubing from the currently used Cellpack.
2. Wipe the tube's surface with a lint free tissue, and then rinse with distilled water.
3. Press the "Start" switch which is behind the aspiration probe.
4. The following message will appear: "Replenish Diluent".
5. Put the reagent aspiration tubes into the appropriate Labnovation box.
6. Press the "Help" key on the keypad followed by "1"-Reaspiring Diluent.
7. Remove the Stromatolyser WH tubing from the currently used Stromatolyser WH reagent box.
8. Wipe the tube's surface with a lint free tissue and then rinse with distilled water.
9. Put the Lyse tube into the appropriate Labnovation box.
10. Press the "Select" key on the keypad followed by "4" – Replace Lyse to prime the lyse reagent system.
11. Repeat this step (step no. 10) 4 times.
12. Run at least 5 backgrounds by pressing the "Start" switch. Check for the acceptable values.
13. Run all three levels of a dedicated commercial control. Check for the acceptable values.
14. Recalibrate the unit if necessary.

## 【PRECAUTIONS】

1. Re-calibration is recommended when changing the origin of the reagents.
2. It is not recommended using different brands of solutions together on one unit of analyzer.
3. Do not ingest. Avoid contact with skin and eyes.
4. Skin contact: Flush with plenty of water.
5. Eyes contact: Flush with plenty of water. Obtain immediate medical attention.
6. This kit is for professionals and in vitro diagnostic use only.

## 【STORAGE】

20-10-020	Diluent-ST	24 months
21-30-SWH	Lyser-WH	18months
21-CL-50	Cell-Cleaner	12months

Store the product at 5-35°C in a dust free place with good ventilation. Temperature outside the range may damage the product.

The open vial life of the product is 60 days. Unused reagents should be discarded after 60 days. Do not mix residual with a new carton.

## 【GLOSSARY OF SYMBOLS】

 In Vitro Diagnostics ;  Lot Number ;

 Temperature Limit ;  Reference Number ;

 Avoid Sunlight ;  Expiry Date ;

 Take Operation Manual for reference ;

**Wellkang Ltd**  
Suite B, 29 Harley Street  
LONDON, W1G 9QR, U.K.  
www.CE-Marking.EU

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20-10-020	DILUENT-ST	20L
21-30-SWH	LYSER-WH	500mL x 3
21-CL-50	CELL-CLEANER	50mL

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DILUENT-ST	NaCl≤9.4%, Buffers≤1.6%, Anti-Microbial Agent≤1.0%
LYSER-WH	Quaternary Ammonium Salt<10%, Sodium Lauryl Sulfate<1.5%
CELL-CLEANER	Sodium Hypochlorite 5%

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