Automation NIHON KOHDEN™ hematology CELL COUNTER REAGENTS IVD

NIHON KOHDEN Cell counters

REAGENT	RIF
ISOTONIC CT DIFF	4590.3
LYSE CT DIFF	4593.9 - 4593.9/5
ENZYDET	4592.1

1.0 PRODUCT IDENTIFICATION

Name : ISOTONIC CT DIFF , LYSE CT DIFF , ENZYDET

2.0 INTENDED USE

In vitro diagnostic use only. Reagents to be used on NIHON KOHDEN® mod. MEK-6400, MEK-6318, WBC 5-part differential.

3.0 STORAGE AND SHELF LIFE AFTER FIRST OPENING

Store at +15 / +30 ° C. Install each reagent after checking package integrity and validity of the device. Once opened, the reagent can be used until exhausted within 60 days. In case of temporary suspension of use, use the cap to close the undamaged original packaging.

4.0 COMPOSITION

Active ingredients : See table 1

5.0 WARNING AND PRECAUTION

Labelling information: See table 2

6.0 DISPOSAL INFORMATION

Unused product and emptied container should be disposed in accordance with local environmental regulations. Instrument wastes, as potential bio hazardous, must be handled with care and disposed in accordance with local environmental regulations.

7.0 SPECIMEN

The devices have been designed to be used with whole blood samples. Refer to Operator's reference manual for sample collection and storage conditions.

8.0 PROCEDURE

ISOTONIC CT DIFF is an isotonic solution (diluent) ready to use.

LYSE CT DIFF is a HGB cyanide solution ready for use.

ENZYDET is an enzymatic detergent ready to use.

9.0 INSTALLATION PROCEDURE

Refer to the Operator's reference manual for reagents installation and priming and for instrument operation. When installing a new batch of dilu-

ent or lysing reagent, check the background according to Operator' Reference manual. Refer to Operator's Reference manual for working temperature. Use each reagents in its own packaging (do not transfer). Follow GLP recommendations. If the instrument connections were not compatible, work according to Good Laboratory Practice (GLP), in order to avoid any kind of contamination, either from external environment to the device, or from the device to the external environment.

When installing LYSE for the first time, prime the instrument with the reagent, after completely draining the line, in order to avoid any kind of contamination which could affect system performance. After priming the instrument, run some background test until the background for WBC and HgB is within Operator's Reference manual limits. It is also recommended to verify instrument calibration.

10.0 METHODOLOGY AND CALCULATION OF ANALYTICAL RESULTS

Refer to Operator' Reference manual for information on methodology and calculation of analytical results.

11.0 INTERNAL QUALITY CONTROL

11.1 - Background verification

Before running testing samples check background count according to Operator's Reference Manual, and verify that all the values for each parameter are within Operator's reference manual limits. If not follow Operator's Reference manual prescriptions. In case of problems attributed to one of the reagents, install a new unit of reagent and verify the background again. Do not run testing samples until the background is within the specified limits.

11.2 - Calibration and calibration verification procedures

Instrument calibration should be verified in the following conditions:

1. after first installation
2. at least every month (or according to manufacturer's requirements for calibration frequency, if more stringent)
3. after a major repair or a major change (replacement of critical parts that may influence test performance, software upgrade, complete change of reagents)
4. when quality control data reflect a significant change of baseline accuracy (drift) or an out-of-control condition.

Calibration procedures should be preceded by a preventive maintenance (verification of baseline accuracy and precision: ref. NCCLS document H38-P).

Calibration and calibration verification must be performed using appropriate calibration material (fresh whole blood samples which have been analysed by reference methods, or appropriate calibrator), FOLLOWING Operator' Reference manual prescriptions. When calibration is required, verify instrument calibration immediately after calibrating, according to Operator's Reference Manual instructions (refer also to NCCLS document H38-P).

11.3 - Daily quality control procedure

In order to assure the accuracy and reliability of the test results, it is recommended to perform quality control procedures, according to Operator' Reference Manual instructions, using appropriate control material (retained patient specimens or stabilised control blood), immediately before, during

or immediately after each analysis run (intended as an analysis sequence). It is recommended to run at least two replicates of at least two levels of stabilised control blood (the Normal and the Low level: refer to NCCLS document H38-P) to the beginning of each run and to repeat the Normal level also during instrument operation (at least every 8 hours), at the end of each run and when installing a new batch of reagent.

Table 1

Reagent Name	REF	Active ingredients	Content
ISOTONIC CT DIFF	4590.3	Sodium chloride	< 5 g/l
LYSE CT DIFF	4593.9	Quaternary ammonium salt	< 4 g/l
LYSE CT DIFF	4593.9/5	Quaternary ammonium salt Potassium cyanide	< 4 g/l < 0.4 g/l
ENZYDET	4592.1	Protease	< 1 %

Table 2

Reagent name	REF	Information on labelling
ISOTONIC CT DIFF	4590.3	No hazard label required according to the directive 1999/45/EC
LYSE CT DIFF	4593.9	No hazard label required according to the directive 1999/45/EC
LYSE CT DIFF	4593.9/5	No hazard label required according to the directive 1999/45/EC
ENZYDET	4592.1	No hazard label required according to the directive 1999/45/EC

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