Document Detail

Type: 0290-MKT-DFU
Document No.: 0290-MKT-DFU-0007-ML[2]
Title: Retainer Retriever IFU
Comment:
Status: CURRENT
Effective Date: 30-Jul-2021

Approval

<u>Owner Role</u>	<u>Sign-off By</u>	<u>Sign-off Date</u>	
0281-Document Control	Tessa Lee	30-Jul-2021 6:34 pm	GMT
0281-Document Control			
0280-Document Owner Document Owner	Margaret Jacob	30-Jul-2021 6:38 pm	GMT
0281-Document Control	Tessa Lee	30-Jul-2021 8:10 pm	GMT
0281-Document Control			

Reference

<u>Document No.</u>	<u>Title</u>	<u>Relation</u>
0290-DMR-041 [9]	Device Master Record for Retainer Retriever	Related

Retainer Retriever/Appliance Remover Tool

ENG

DIRECTIONS FOR USE

INTENDED USE:

Retainer Retriever and/or Appliance Removal Tool helps/aids the patient/user in taking removable orthodontic appliances out with minimum effort while protecting the appliance from damage and/or scratches. The device is intended to be used with other, traditional removable orthodontic devices.

INDICATIONS FOR USE:

Indicated for patients wearing thermoformed aligners or retainers.

PERFORMANCE CHARACTERISTICS AND CLINICAL BENEFITS:

This device has a small hook on one end and a handle on the other to aid the patient in the removal of their aligner or retainer.

Clinical benefits of orthodontic treatment may include overall improvement of tooth health and correction of malocclusion.

CONTRAINDICATIONS:

None known.

WARNINGS:

Single patient use only.

Use of excessive force to remove retainer or aligner may cause damage to it.

PRECAUTIONS:

Do not reuse the product on another person. Risks of reuse include infection and contamination (i.e. bacterial, viral).

Report serious incidents that have occurred in relation to the device to your healthcare professional or the device manufacturer and/or the Competent Authority in Member State of residence, or other applicable regulatory agency.

ADVERSE REACTIONS:

None known.

STEP BY STEP INSTRUCTIONS:

1. Grasp the Retainer Retriever by the end with the hole in it.
2. Insert the hooked end of the Retainer Retriever into the mouth.
3. Using the hook, gently grasp an edge of the retainer or aligner near the front of the mouth and pull away from the teeth.
4. If resistance is felt, it may be necessary to repeat step 3 in several locations along the teeth. Do not use excessive force to try and remove the retainer or aligner.

SAFE DISPOSAL OF DEVICE:

Dispose of in the regular household trash, as you would your toothbrush or used dental floss.

WARRANTY AND LIMITED REMEDY

Manufacturer warrants that this product is free from defects in materials and manufacture.

Manufacturer's sole obligation and customer's sole remedy in the event of a claimed defect shall be limited to replacement merchandise or refund of the purchase price.

All implied warranties, including implied warranties of merchantability or fitness for any particular purpose or use are disclaimed.

LIMITATION OF LIABILITY

Except where prohibited by law, Manufacturer will not be liable for any loss or damages arising from this product, whether direct, indirect, special, incidental, or consequential, regardless of the theory asserted, including warranty, contract, negligence or strict liability. This limitation does not apply to third party personal injury claims.

 Manufacturer: **Dentsply Sirona Orthodontics Inc.**
7290 26th Court East Sarasota, FL 34243 USA

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