

NON MAGNETIC STOOL ADJUSTABLE HEIGHT 3T COMPATIBLE



Référence du produit : IRM1062

Specifications

Non magnetic stool with adjustable height and mobile step for MR rooms up to 3 Tesla.

TECHNICAL DATA:

- Reinforced nylon five stars mobile step
- Five twins castors without brake
- Inox steel screw for height adjustment
- Seat min. Height 520mm
- Seat max. Height 700mm
- Seat diameter 350 mm
- Seat thickness 70 mm
- Upholstery blue color
- Upholstery fire resistant M1
- Upholstery with antibacterial treatment
- Weight : 5,0 Kg