


HỆ THỐNG ĐÈN ĐỘI ĐẦU PHẪU THUẬT
HÃNG: SUNOPTIC TECHNOLOGIES LLC, MỸ
XUẤT SỨ: MỸ

MẪU NHÃN HÀNG HÓA ĐƯỢC LƯU HÀNH TẠI VIỆT NAM

SUNOPTIC TECHNOLOGIES®
 6018 BOWDENDALE AVENUE -C -D
 JACKSONVILLE, FL 32216, USA
REF P/N: TITANX450ST -A -K -J
DESC: 400W XENON LIGHTSOURCE
S/N: 450-XXXXY 
 +M678TX0H
 +\$\$+7450-XXXXY/16D20190520H1

EC REP AJW Technology Consulting GmbH **CE**
 Königsallee 106
 40215 Düsseldorf, Germany
 +49 (0) 211 3013 2232
IF NOT INCLUDED, INSTRUCTIONS FOR USE (VUL) DOCUMENTS
 FOR THIS PRODUCT ARE AVAILABLE AT WWW.SUNOPTIC.COM UNDER THE BROCHURE/FU
 HEADERS. OTHER COPIES MAY BE OBTAINED BY CONTACTING CUSTOMER SUPPORT AT
 877677.2822
FORM: CCF-304

	Shipping	Storage
Temp F° / C°	88° - 104° F (20° C - 40° C)	-4° - 140° F (-20° C - 60° C)
Relative Humidity	30% - 85%	0% - 95%
Altitude	700 to 1060 hPa	700 to 1060 hPa

SUNOPTIC TECHNOLOGIES®
 6018 BOWDENDALE AVENUE -A -D
 JACKSONVILLE, FL 32216 -C -K
REF P/N: TITANX350ST -C -K
DESC: 350W XENON LIGHTSOURCE
S/N: 350-XXXXD 
 +M678TX0H
 +\$\$+7350-XXXXD/16D20190520HM

EC REP AJW Technology Consulting GmbH **CE**
 Königsallee 106
 40215 Düsseldorf, Germany
 +49 (0) 211 3013 2232
IF NOT INCLUDED, INSTRUCTIONS FOR USE (VUL) DOCUMENTS
 FOR THIS PRODUCT ARE AVAILABLE AT WWW.SUNOPTIC.COM UNDER THE BROCHURE/FU
 HEADERS. OTHER COPIES MAY BE OBTAINED BY CONTACTING CUSTOMER SUPPORT AT
 877677.2822
FORM: CCF-304

	Shipping	Storage
Temp F° / C°	68° - 104° F (20° C - 40° C)	-4° - 140° F (-20° C - 60° C)
Relative Humidity	30% - 85%	0% - 95%
Altitude	700 to 1060 hPa	700 to 1060 hPa

SUNOPTIC TECHNOLOGIES™
 6018 BOWDENDALE AVENUE -A -D -J
 JACKSONVILLE, FL 32216, USA -C -K
REF SG35-110-VHB
DESC: HEADLIGHT
S/N XXXXX 
 +M678HLO
 +\$\$+7XXXXX-O

EC REP AJW Technology Consulting GmbH **CE**
 Königsallee 106
 40215 Düsseldorf, Germany
 +49 (0) 211 3013 2232
IF NOT INCLUDED, INSTRUCTIONS FOR USE (VUL) DOCUMENTS
 FOR THIS PRODUCT ARE AVAILABLE AT WWW.SUNOPTIC.COM UNDER THE BROCHURE/FU
 HEADERS. OTHER COPIES MAY BE OBTAINED BY CONTACTING CUSTOMER SUPPORT AT
 877677.2822
FORM: CCF-304

FORM: CCF-304

SUNOPTIC TECHNOLOGIES®
 6018 BOWDENDALE AVENUE -A -D -J
 JACKSONVILLE, FL 32216, USA -C -K
REF P/N: SSL-4000T-SC
DESC: LED 4000 LIGHTSOURCE
S/N: XXXXXXXXY 
 +M678LED400007
 +\$\$+7XXXXXXXXY/16D2019052072

EC REP AJW Technology Consulting GmbH **CE**
 Königsallee 106
 40215 Düsseldorf, Germany
 +49 (0) 211 3013 2232
IF NOT INCLUDED, INSTRUCTIONS FOR USE (VUL) DOCUMENTS
 FOR THIS PRODUCT ARE AVAILABLE AT WWW.SUNOPTIC.COM UNDER THE BROCHURE/FU
 HEADERS. OTHER COPIES MAY BE OBTAINED BY CONTACTING CUSTOMER SUPPORT AT
 877677.2822
FORM: CCF-304

	Shipping	Storage
Temp F° / C°	32°F - 104°F (0°C - 40°C)	-4°F - 140°F (-20°C - 60°C)
Relative Humidity	0% to 85% non-condensing	0% to 95% non-condensing
Altitude	700 to 1060 hPa	700 to 1060 hPa


TÂN TEO LAY
GIÁM ĐỐC