



Beckman Coulter Ireland Inc.,
Lismeehan,
O'Callaghan's Mills,
Co. Clare,
Ireland.
Tel. +353-65-6831100.
Fax. +353-65-6831122.

Technical Specification

iQ Lamina

Part Number: 800-3236

- Intended use:** iQ Lamina is intended for in vitro diagnostic use with iQ200 Series Analyzers.
- Packaging**
 - 2 x 7L iQ Lamina
 - 1 x 5 μ m Filter
- Storage:** 20 – 28°C
- Physical and Chemical Specification**

Attribute	Specification
Physical State	Liquid
Color	Colorless
Odor	Odorless
pH	6.5 – 7.5
Specific Gravity at 20°C	1.007 \pm 0.001
Clarity	Clear
Counts	\leq 50/0.5 mL

Authorized Signatory & Date:

Marguerita Sweeney 27th June 2019
Marguerita Sweeney
Regulatory Affairs Manager
Beckman Coulter Ireland Inc.
Phone: 353 65 6831495

Registered Office: Lismeehan, O'Callaghan's Mills, Co. Clare. Directors: Piers Devereux, Paul Henry, Orlaith A. Lawler

Incorporated in Republic of Ireland with limited liability. Co. Reg. No. 501471.