

mẫu hộp Băng cá nhân
QC: 9,2 x 10,0 x 4,5 cm

102 plaster
1.9 x 7.2 cm

HOUGOSTICK
Durable Elastic Fabric

HOUGOSTICK
Durable
Elastic Fabric



102 plaster
1.9 x 7.2 cm

HOUGOSTICK
Durable Elastic Fabric

Imported by: Tan Cuong
trade and phar co.,ltd

HOUGOSTICK
Durable Elastic Fabric
1.9 x 7.2 cm
102 MIENG



HOUGOSTICK
Durable Elastic Fabric

- DÍNH CHẮC KHI SỬ DỤNG
- THÔNG THOÁNG
- MỠ BĂNG KHÔNG ĐAU

Composition:

- Ribbon: Viscose and Polyamide elastic fabric
- Dress: Viscose with topmost layer made by polyethylene prevent sticky
- Sticker: Zinc oxide without solvent
- One piece put in small pack

Indications:

individual dressing clear, high agglutinated, elastic fabric to protect small injury,graze, chapped skin

Directions for use:

- Hygiene and dry the skin before application the dressing
- Change the dressing daily

Storage: Cool and dry place

Thành phần:

- Băng: Vải Viscose và Polyamide co giãn
- Gạc: Viscose, có lớp màng trên cứng bằng polyethylene không gây dính
- Chất dính: Ôxít kẽm không dung môi
- Mỗi miếng đựng trong từng bao riêng

Cử dụng:

Băng dính cá nhân thông thoáng, độ dính cao, độ co giãn tốt bảo vệ các vết thương nhỏ, vết trầy xước, rách da

Chăm sóc:

- Làm vệ sinh da, lau khô trước khi dán băng
- Thay băng hàng ngày

Lưu ý: Nơi thoáng mát và khô ráo

Manufactured by:
Zhejiang Bangli Medical Products Co., Ltd.P.R.China

