

Technical Specification

SFRI Cleaner 3.1

Part Number: HSG305

1. Intended use: Cleaning reagent for the cleaning of blood cells counters.
2. Packaging Specification
Bottle: PEBD Colorless Cubitainer
Cap: Colorless plastic cap
Packaging: Carton white box
3. Physical and Chemical Specification
Storage: Store at room temperature

Parameter	Test method	Specification
Physical Status	NA	Liquid
Colour	NA	Colorless
Odor	NA	Odorless
Volume	Weight	5L
pH	pH-meter	7,2 ± 0,3 (20°C)
Resistivity	Resistivimeter	70 ± 6
Osmotic pressure	Osmometer	320 ± 20

Signature