

Technical Specification

SFRI Clair 3.1

Part Number: HSC102

1. Intended use: Hypochloric solution for the cleaning of blood cells counters.
2. Packaging Specification
Bottle: PEHD White Bottle
Cap: Plastic white cap
Packaging: Kraft box
3. Physical and Chemical Specification
Storage: Store at room temperature

Parameter	Test method	Specification
Physical Status	NA	Liquid
Colour	NA	Pale Yellow
Odor	NA	Odorless
Volume	Weight	60mL
pH	pH-meter	12,5 ± 0,5 (20°C)
Resistivity	Resistivimeter	21± 6
Osmotic pressure	Osmometer	NA

Signature