

Technical Specification

Na Conditioner

Part Number: R04S028

1. Intended use: Cleaning solution for sodium electrode of SFRI ISE.

2. Packaging Specification

- Bottle:

Dimension: 62 x 30

Colour: Colorless

Material: Polypropylene

- Cap:

Dimension: High 25 mm

low 5 mm High: 25, Low: 0.5

Colour: White

Material: Plastic

- Packaging: 30mL

3. Physical and Chemical Specification

Storage: Store at room temperature

Parameter	Test method	Specification
Appearance	NA	Liquid
Colour	NA	Colorless
Odor	NA	Odorless
Volume	Weight	30mL
pH	pH meter	4 ± 1 (20°C)

Signature