

Technical Specification

Cuvette Ring V2 with Mixer 1.0×4.0

REF NUMBER: 80911020

INTENDED USE: Intended to be used with Dialab Coalyser

PRODUCT DESCRIPTION: Plastic ring containing 32 single use cuvettes for automated testing of Hemostasis parameters with Dialab's fully automated coagulation analyser Coalyser

TECHNICAL INFORMATION:

- Material: Polystyrol Styron 678E
- Packaging Size: 6 boxes containing 10 rings each together with a cuvette card
- Dimension: 13.7cm diameter, 2.1 cm high
- Path length: 7mm
- Wavelength Range (nm): 450nm – 740nm
- Fill Volume: 500µl
- Storage: dry, dustfree



Product Manager

6.6.2019

Date