

Eluent A

1320ML

For use with HA-1500 Automated HbA1c Analyzer.
Pour utilisation avec l'analyseur HA-1500.
Verwendung mit HA-1500 Analyzer.
Para uso en los analizadores HA-1500.
Con analizzatore HA-1500



Active Ingredients:

Phosphate: 40mmol/L;
Sodium chloride: 20mmol/L;
Perserver: ≤0.5mL/L.



EC REP MedNet GmbH
Borkstrasse 10 · 48163
Muenster · Germany



Mti-diagnostics GmbH
Limburger Straße 45 D-65510 Idstein, GERMANY
Tel: +49(0)6126-9595 262 Fax: +49(0)6126-9596 264



Ver.1709

104.86 mm

109.96 mm

Eluent B 920ML

For use with HA-1500 Automated HbA1c Analyzer
Pour utilisation avec l'analyseur HA-1500.
Verwendung mit HA-1500 Analyzer.
Para uso en los analizadores HA-1500.
Con analizzatore HA-1500

Active Ingredients:

Phosphate: 40mmol/L;
Sodium chloride: 170mmol/L;
Perserver: ≤0.5mL/L.



Mti-diagnostics GmbH

Limburger Straße 45 D-65510 Idstein, GERMANY
Tel: +49(0)6126-9595 262 Fax: +49(0)6126-9596 264

Ver.1709

98.09 mm

70.19 mm