

# TIANCHANG DONGAN PROTECTIVE EQUIPMENT CO., LTD.

HONGMIAO VILLAGE, TIANCHA ROAD, TIANCHANG CITY, CHINA

Date: July 6th, 2019

## TECHNICAL DOCUMENTS SPECIFICATION

DISPOSABLE MEDICAL FACE MASK EAR LOOP

DISPOSABLE MEDICAL FACE MASK TIE ON

1, The Nonwoven Medical Face Mask with Earloop & Nonwoven Medical Face Mask with Tie On are made of non-woven fabric, melt-blown non-woven fabric, plastic nose wire, nylon earloop and tie on;

2, It consists of outer ply, middle ply, inner ply, nose wire, ear loop, tie on;

3, The medical mask is hydrophobic and breathable, and has a remarkable filtering effect on micro-virus aerosol or harmful dust. It has anti-wet, anti-bacterial and anti-leakage functions. Commonly used in hospitals, clinics, ICUs, operating rooms and other medical facilities.

4, Applications: The medical mask is light, comfortable, tasteless, glass-free, non-irritating to the skin, and waterproof. The nose clip of the medical mask is soft, comfortable and adjustable. Breathing resistance is small and wears hygienic. Medical masks can absorb particulate matter, aerosols, dust, smoke, mist, poisonous gas and toxic vapor in the air through the filter material, and block it without being inhaled by humans.

5, Two Types: Ear Loop & Tie On.

6, Packing:

(1) 50PCS/INNER BOX, 2000PCS/CTN;

(2) 1PC/BAG, 50PCS/INNER BOX, 2000PCS/CTN;

陆宇佳



TIANCHANG DONGAN PROTECTIVE EQUIPMENT CO., LTD.