



Technical Specification
MicroScan MH BROTH with 3% LHB

Part Number: B1015-25

- 1. **Intended use:** For use with MicroScan Panels
- 2. **Packing:** 10 x 25 mL
- 3. **Storage:** 2 – 8°C
- 4. **Physical and Chemical Specification**

Attribute	Specification
Physical State	Liquid
Color	Tan
Odor	Odorless
Transparency	Clear
pH	≈ 7.2

- 5. **Concentration of main ingredients:**

Mixture	Specification
Mueller Hinton Broth	Per CLSI concentrations
Lysed Horse Blood	≈ 3%
HEPES	12 gms

Christina Thomas
Microbiology Regulatory Affairs