

Tài liệu mô tả chức năng thông số kỹ thuật TTBYT

Tài liệu được xác nhận bằng chữ ký số và có hiệu lực kể từ ngày ký.

Hà Nội, ngày 16 tháng 08 năm 2019

Người đại diện hợp pháp của cơ sở

Giám đốc
Trịnh Diệu Hương

Technical Specification

Product name: Magazine Pack Peripheral Blood Orange 10 pcs

Part number: XU-10133-01

Intended use: For use to load slides on automated cell-locating devices which automatically locates and presents images of blood cells on peripheral blood smears. The operator identifies and verifies the suggested classification of each cell according to type.

Storage condition: Room temperature

Parameter	Specification
General Inspection	The magazine barcode must be correct Orange magazines with barcodes consisting of prefix "PB" followed by at least three digits - a peripheral blood analysis will be performed.
Appearance	The color of magazines is orange; The magazines are supplied with barcode labels with magazine ID
Material	Plastic
Capacity	12 slides



Ann-Christin Johansson
QA & Regulatory Specialist
CellaVision AB